

FIA-MRMS BASED METABOLOMICS REVEALS ANTI-INFLAMMATORY EFFECTS OF HOP BITTER ACIDS IN DENDRITIC CELLS



FP 551, ASMS 2021

Emanuela Salviati¹, Eduardo Sommella¹, Giulio Verna¹, Aiko Barsch², Matthias Witt², Jochen Friedrich², Pietro Campiglia¹

¹ University of Salerno, Department of Pharmacy, Fisciano, Italy; ² Bruker Daltonics GmbH & Co. KG, Bremen, Germany

Introduction


Inflammation is a complex, evolutionarily conserved process, that involves immune and non-immune cells in the host for protection from harmful stimuli. Metabolic changes occur in macrophages and dendritic cells (DCs) leading to pro or anti-inflammatory phenotypes. In this regard, the growing field of immunometabolism aims to target specific metabolic pathways to modulate inflammation, and natural compounds are more and more used for preventive action [1]. To understand the metabolic shift toward anti-inflammatory phenotypes, metabolomics has emerged as leading approach. The metabolome profiling is usually carried out by liquid chromatography coupled to tandem mass spectrometry (LC-MS/MS). As an alternative to LC-MS/MS, direct infusion magnetic resonance mass spectrometry (DI-MRMS) or flow injection (FIA) MRMS delivers short analysis time for sample screening and features unmatched mass accuracy, resolution, sensitivity, and isotopic fine structure. Hence, the FIA- or DI-MRMS workflows provide a higher throughput in metabolic profiling compared to LC-MS methods [2]. In this study, MRMS was used to highlight the metabolic changes in dendritic cells stimulated with lipopolysaccharide (LPS) and co-treated with a Hop derived fraction rich in beta acids and prenylflavonoids and compared their effect with the flavonol quercetin.

Methods

Sample preparation and metabolite extraction

DCs from murine bone marrow (BMDCs) from six- to eight-week-old mice.

24 h later. LPS stimulation (1 µL/mL)



Pelleted cells thawed on ice


Sonicate 6 min, 25 °C + vortex for 30 sec

14680 rpm, 10 min, 4°C.

100 µL of ice-cold MeOH/H₂O (80:20 v/v).

The supernatants were dried under nitrogen and reconstituted in 100 µL of ACN/H₂O (70:30) + 0.1% HCOOH (v/v) before MRMS aXelerate analysis.


ELISA assay of cell supernatants



RNA extraction and qPCR analysis

GAPDH, *Slpi*, *Ap-1*, *Hmox1*, *Ferroportin*, *Nrf2*, *Nqo1*

MRMS aXelerate workflow using nESI



Automated multisample chip-based nESI sample ionization platformed using a TriVersa NanoMate (Advion BioSciences Ltd, Ithaca NY, U.S.A), coupled to a solarix XR 7T (Bruker Daltonics GmbH & Co. KG, Bremen, Germany) MRMS system.

Tuning and external calibration: NaTFA (0.1 mg/mL in 50% acetonitrile); **Mass range:** 90–800 *m/z*; **Ion accumulation time:** 10 ms; **Number of single scan added:** 32; **Data size:** 2 million data points (2M). Nebulizing (N₂) and drying gases (N₂): 1 and 4 L/min, respectively; Drying gas temperature: 200 °C.

Ionization mode: positive and negative ionization. Five measurement replicates of each sample were performed.

Data processing

Peak alignment, tentative annotation of compounds based on accurate MS measurements: **MetaboScape 2021** (Bruker Daltonics GmbH & Co. KG, Bremen, Germany). Adducts were combined and deisotoped during the feature generation.

MS database for compound identification: HMDB and LipidMAPS. **Statistical significance:** One-way ANOVA with the Tukey post hoc test analysis for aligned, log-transformed and scaled data. **Classification models:** Principal component analysis (PCA) and partial least square discriminant analysis (PLS-DA). Metaboanalyst 5.0 was used for statistical analysis.

Results

DCs stimulated with LPS pre-exposed to 25 µg/mL of HOP B and C fractions showed reduced production of several inflammatory cytokines (IL-6, IL-1α, IL-1β and TNF). The fraction C was significantly more effective than the fraction B. Based on the results of inflammatory cytokine and qPCR measurements, the DCs samples were divided into four groups: control (CTRL), LPS stimulated DCs (LPS), LPS-DCs plus HOP C (LPS+H) and plus quercetin (LPS+Q). As can be observed in the PLS-DA score plots in Figure 1A, metabolites in the different DCs groups were separated into distinct clusters. The metabolite abundance variations induced by the LPS stimulation and after the two treatments were visualized by a heatmap shown in Figure 1B.




Figure 1: PLS-DA score plot (A) and Heat map (B) of statistically relevant (ANOVA followed by a Tukey's post hoc, *p* < 0.05) DCs metabolites modulated by LPS, HOP C, and quercetin. Color changes reflect normalized intensity

This untargeted metabolomics study showed that quercetin and HOP C induced changes are segregated in two independent pathways. LPS-activated DCs are characterized by anaerobic glycolysis and a pro-inflammatory phenotype [3]. As shown in the Whisker box plots in Figure 2, the preventive administration of HOP C acts through a regulation of the oxidative and non-oxidative branch of pentose phosphate pathway, regulates the arginine-citrulline and arginine-succinate shunt, and normalizes the nucleotide metabolism, differently from quercetin. These results suggest that HOP C fraction is effective to impair the inflammatory response by regulating the metabolic reprogramming of DCs toward anti-inflammatory and resting conditions.




Figure 2: Normalized intensity of altered metabolites involved in metabolic reprogramming of LPS-DCs and modulated by HOP C and quercetin treatments: (A) Pentose phosphate pathway intermediates, (B) TCA cycle and citrulline, (C) Purine/Pyrimidines and Nucleotides

Conclusions

- MRMS aXelerate is a fast and reliable method for metabolome profiling.
- MRMS aXelerate allows to take the snapshot of metabolic changes in BMDCs stimulated by LPS and co-treated with a Hop bitter acids and prenylflavonoids fraction.
- This data further highlights the potential role of natural compounds in targeting immunometabolism for the adjuvant and preventive therapies for the treatment of inflammatory diseases.

MRMS Metabolomics

References

- Immunometabolism: Cellular Metabolism Turns Immune Regulator. R. M. Loftus and D. K. Finlay, *J. Biol. Chem.*, **2016**, 291, 1–10.
- Analysis of the metabolic switch induced by the spirulina peptide SP6 in high fat diet ApoE^{-/-} mice model: A direct infusion FT-ICR-MS based approach. E. Sommella, A. Carrizzo, F. Merciai, V. Di Sarno, D. Carbone, M. De Lucia, S. Musella, C. Vecchione and P. Campiglia, *J. Pharm. Biomed. Anal.*, **2021**, 195, 113865.
- Toll-like receptor-induced changes in glycolytic metabolism regulate dendritic cell activation. C. M. Krawczyk, T. Holowka, J. Sun, J. Blagih, E. Amiel, R. J. DeBerardinis, J. R. Cross, E. Jung, C. B. Thompson, R. G. Jones and E. J. Pearce, *Blood*, **2010**, 115, 4742–4749.