

Trapped ion mobility separation (TIMS) of Glucose-6-Phosphate and Fructose-6-Phosphate

Surendar Tadi¹, Xuejun Peng², Beixi Wang², Erica Forsberg², Matthew Lewis³. ¹Bruker Scientific LLC, 40 Manning Road, Billerica, MA 01821, USA, ²Bruker Scientific LLC, 101 Daggett Drive, San Jose, CA 95134, USA, ³Bruker Daltonics GmbH & Co. KG, Fahrenheitstraße 4, 28359 Bremen, Germany

Introduction



Evaluating central carbon metabolism in living systems is challenging due to the dynamic range of metabolite concentrations. Among sugars and organic acids, phosphorylated sugars glucose-6-phosphate (G6P) and fructose-6-phosphate (F6P) are biologically relevant isomeric compounds. Although there are effective chromatographic separation methods [1], it is often difficult to separate in high throughput studies. In this study, G6P and F6P are rapidly separated based on ion mobility on a Bruker timsTOF Pro 2. Trapped ion mobility spectrometry (TIMS) enables rapid separation of isomeric species based on size and shape in the gas phase for high-throughput applications

Here we applied ion mobility separation using TIMS as a 4th dimension of separation and annotation tool using collisional cross section (CCS). Each metabolite was isolated with baseline mobility peak separation of G6P and F6P with definitive CCS values. This enables rapid identification and quantitation of G6P and F6P in biological matrices using rapid flow injection analysis and short chromatographic gradients

Methods

G6P and F6P were purchased from Sigma-Aldrich (St Louis, MO, USA), LC-MS grade water was purchased from Fisher Scientific. A 1 mg/mL stock solution was prepared in LC-MS grade water further diluted to 10 µg/mL working solution. Data was acquired with a Bruker Elute UHPLC on Phenomenex Kinetex F5 (150 x 2.1 mm, 2.6 µm, 100Å) with column temperature at 40°C. Mobile phase A was 0.1% formic acid in H₂O and B was 0.1% formic acid in acetonitrile. Gradient conditions are 0% B at 1 to 2 min, up to 95% at 8 min with a 2 min hold and equilibrated for 5 min on a VIP-HESI timsTOF Pro 2 with TIMS enabled in ESI negative mode. A previously established 4D-Metabolomics method with a mass range of 20 – 1300 m/z and mobility range of 1/K₀ 0.55 – 0.85 V-s/cm² and ramp time of 300 ms was used. Other acquisition parameters were further optimized by adjusting ion mobility resolution mode, accumulation time, duty cycle. Mass and CCS calibration were performed prior to data analysis. Data processing was conducted in Bruker Data Analysis 5.3 software.


Results



The sugar phosphate reference standard mixture using LC-VIP-HESI-timsTOF Pro 2 in negative mode co-elute chromatographically yet had baseline mobility peak separation and therefore can be distinguished. The average CCS deviation is 1.6% compared to reported reference values in the CCS Compendium [2].

Table 1. CCS value comparison to published reference data

Name	[M-H] ⁻	Mobility 1/K ₀	CCS (Å) ² measured	CCS (Å) ² reported	CCS (Å) ² Error %
Glucose-6-phosphate	259.0244	0.697	148.7	146.7	1.3
Fructose-6-phosphate	259.0244	0.707	146.9	144.3	1.8



References

[1] Wang, Y., et al., "Baseline Separation of Six Hexose Phosphate Isomers by Liquid Chromatography-Mass Spectrometry from Tissues" Bio-protocol Preprint, 2022, bio-protocol.org/prep1541.

[2] Picache, J.A., et al., "Collision cross section compendium to annotate and predict multi-omic compound identities", Chem. Sci., 2019, 10, 983-993.

[3] Song J., et al., "The Multiple Roles of Glucose-6-Phosphate Dehydrogenase in Tumorigenesis and Cancer Chemoresistance", Life, 2022, 12, 271.

Conclusion

- TIMS-MS based method for separation of Glucose-6-phosphate and Fructose-6-phosphate
- Mobilgrams can provide relative quantitation of isoforms in biological matrices .
- This method is compatible with direct infusion or short LC gradients
- Base peak separation of isoforms will distinguish between two analytes by comparing to reference CCS values reported in database or can be validated with reference standards

LC-TIMS-MS/MS Isomer Separation