



REUNION DES UTILISATEURS RMN BRUKER

From **biologics** characterisation to **bioproduction** analytics

Martial Piotto, Ph.D.

Outlook

01 Peptides
Therapeutic peptides: where are we today?

02 Oligonucleotides
Therapeutic oligonucleotides – where can we help?

03 Antibodies
NMR analysis of biotherapeutics

04 Bioproduction
Understanding and Monitoring

05 PIPAc project
Intelligent production of Active Pharmaceutical Ingredient (API)

01 Peptides

Therapeutic peptides, the 'sweet spot' of NMR

High dispersion spectral data

NMR applications

- Secondary and tertiary structure [1, 2, 3]
- Dynamics
- Binding [1,2,3]
- Similarity assessment [4]
- Quantification (absolute and relative) [5]

1. Chiva C, Barthe P, Codina A ... Sakakibara S, Albericio F, Giralt E, *JACS*, 2003;125:1508-1517
2. Codina A, Love JD, Li Y, Lazar MA, Neuhaus D, Schwabe JW, *Proc Natl Acad Sci U S A*. 2005;102(17):6009-6014.
3. Codina A, Benoit G, Gooch J T, Neuhaus D, Perlmann, Schwabe JWR, *JBC*, 2004:279, 53338
4. Haxhom GW, Bent O, Malmstrom J, *J Pharm Sci*, 2019;108: 3029 (2019)
5. Bradley SA, Jackson WC Jr, Mahoney PP, *Anal Chem*. 2019; 5,91(3):1962-1967



FDA encourages the use of state-of-the art technology

Development of Therapeutic Protein Biosimilars: Comparative Analytical Assessment and Other Quality-Related Considerations

Guidance for Industry

DRAFT GUIDANCE

218 Despite improvements in analytical techniques, current analytical methodology may not be able
219 to detect or characterize all relevant structural and functional differences between the two protein
220 products. A thorough understanding of each analytical method's limitations will be critical to a
221 sponsor's successful identification of residual uncertainties and, in turn, to the design of
222 subsequent testing. In addition, there may be incomplete understanding of the relationship
223 between a product's structural attributes and its clinical performance. FDA encourages the use of
224 available **state-of-the-art** technology. Sponsors should use appropriate analytical methodologies
225 that have adequate sensitivity and specificity to detect and characterize differences between the
226 proposed product and the reference product.

<https://www.fda.gov/media/125484/download>

https://www.casss.org/docs/default-source/hos/2021-hos-speaker-presentations/speaker-presentation-chen-kang-cder-fda-2021-.pdf?sfvrsn=8f7418ed_3

- What we don't measure or control may lead to unintended consequences such as loss of potency or efficacy. Or even pathological action.
- High resolution/high sensitivity data => improved product risk assurance (product knowledge)
- MS and NMR are considered high sensitivity and high-resolution state-of-the-art (modern) techniques


<https://www.fda.gov/drugs/news-events-human-drugs/cder-scientists-use-modern-measurement-tools-quality-assurance-and-comparability-complex-drugs>

HOS Characterization of Pharmaceutical Proteins: Insulin and Glucagon like peptide 1 (GLP-1)


Higher-Order Structure Characterization of Pharmaceutical
Proteins by 2D Nuclear Magnetic Resonance Methyl Fingerprinting

Gitte W. Haxholm, Bent O. Petersen, Joan Malmstrøm*

Novo Nordisk A/S, R&D, Novo Nordisk Park, Maaloev, Denmark




Reference insulin analog (teal) and 2 related insulin variants (red, blue) varying by exchange of only one amino acid in the backbone.



Reference GLP-1 analog (blue) with 40% α -helix and 2 synthetic stereoisomers with 25% (red) and 20% (green) α -helix

Glucagon like peptide (GLP) proton NMR spectra @ 600 MHz



- Insulin hyperglycaemia, stable in solution
- Glucagon hypoglycaemia, GLP and analogues treat type II diabetes, obesity
- Glucagon fibrilizes rapidly at the acidic pH required for solubility => formulated as a lyophilized powder that is reconstituted in an acidic solution immediately before use

02 Oligonucleotides

Therapeutic Oligonucleotides

Eighteen nucleic acid therapeutics approved by the FDA as of 2023

■ The challenge:

- Characterization of modified therapeutic oligonucleotide fragments (15 to 30 nucleotides)
- Quality control
- Biosimilars
- ...

5	SPINRAZA® (nusinersen), 2016, spinal muscular atrophy (SMA)	SMN2, CNS	18-mer SSO, PS, 2'-O-MOE, splicing modification (exon inclusion)	<p>C: 5-Methyl-C</p>	12.5 mg once every 4 months, IT
6	ONPATTRO® (patisiran), 2018, hereditary ATTR (hATTR)	TTR, liver	21-mer/21-mer siRNA, 2'-OH/2'-OMe/2'-H, RISC Ago2	<p>Passenger 5'-GUAAGCAAGAGUAUUGCAU(TT)-3' Guide 3'-TTCAUUGG(U)UCUCAUAAGGUA-5'</p>	0.3 mg/kg once every 3 weeks, IV, LNP (DLin-MC3-DMA, PEG-DMG lipid, cholesterol, DSPC)
7	TEGSEDI® (inotersen), 2018, hereditary ATTR (hATTR)	TTR, liver	20-mer gapmer ASO, PS, 2'-O-MOE/2'-H, RNase H1	<p>C: 5-Methyl-C</p>	300 mg once weekly, SC
8	GIVLAARI® (givosiran), 2019, acute hepatic porphyria	ALAS1, liver	21-mer/23-mer siRNA, enhanced stabilization chemistry (ESC), PS, 2'-OMe/2'-F, RISC Ago2	<p>Passenger 5'-CAGAAAGAGUGUCUCAUCUUA-3' Guide 3'-UGG(U)GUUCUCACAGAGUAAGAAU-5'</p>	2.5 mg/kg once a month, SC, GalNAc conjugate

Table 1. Summary of approved oligonucleotide-based therapeutics



Therapeutic oligonucleotides



NMR has a role to play:

- Characterize batches of therapeutic oligonucleotides containing mixtures of phosphorothioate diastereoisomers
- Composition is believed to affect the pharmaceutical activity




FIG. 1. Schematic representation of phosphorothioate-modified oligonucleotides with Sp (red, left) and Rp (blue, right) configuration at the phosphorous atom.

NUCLEIC ACID THERAPEUTICS
Volume 32, Number 4, 2022
Mary Ann Liebert, Inc.
DOI: 10.1089/nat.2021.0089


Open camera or QR reader and scan code to access this article and other resources online.



Inhomogeneous Diastereomeric Composition of Mongersen Antisense Phosphorothioate Oligonucleotide Preparations and Related Pharmacological Activity Impairment

Lorenzo Arrico,^{1*} Carmine Stolfi,^{2*} Irene Marafini,² Giovanni Monteleone,² Salvatore Demartis,³ Salvatore Bellinvia,³ Francesca Viti,³ Marie McNulty,⁴ Irene Cabani,⁵ Anita Falezza,⁵ and Lorenzo Di Bari¹

FIG. 3. The ³¹P-NMR spectra of a subset of mongersen batches, revealing differences in the fine structures of the ³¹P resonances.



FDA guidance on Givosiran Sodium

Contains Nonbinding Recommendations

Draft – Not for Implementation



Draft Guidance on Givosiran Sodium

May 2023

This draft guidance, when finalized, will represent the current thinking of the Food and Drug Administration (FDA, or the Agency) on this topic. It does not establish any rights for any person and is not binding on FDA or the public. You can use an alternative approach if it satisfies the requirements of the applicable statutes and regulations. To discuss an alternative approach, contact the Office of Generic Drugs.

The test API sequence, chemical structure and composition including strand composition, duplex vs residual single strands, diastereomeric composition, and P=S to P=O ratios should be compared to those of the API from the RLD using a broad range of orthogonal analytical methods with sufficient sensitivity, discriminating, and resolving power, that could include but are not limited to the following:

- Mass spectrometry (MS), including tandem mass spectrometry (MS/MS)
- Nuclear magnetic resonance (NMR) spectroscopy
- Liquid chromatography (LC)
- Flame atomic absorption spectroscopy (FAAS)
- Duplex melting temperature (T_m)



O⁻ denotes phosphodiester linkage
S⁻ denotes phosphorothioate linkage
Dashed lines denote Watson-Crick base pairing

Abbreviations: Af = adenine 2'-F ribonucleoside; Cf = cytosine 2'-F ribonucleoside; Uf = uracil 2'-F ribonucleoside; Am = adenine 2'-OMe ribonucleoside; Cm = cytosine 2'-OMe ribonucleoside; Gf = guanine 2'-F ribonucleoside; Gm = guanine 2'-OMe ribonucleoside; Um = uracil 2'-OMe ribonucleoside; L96 = triantennary GalNAc (N-acetylgalactosamine)

https://www.accessdata.fda.gov/drugsatfda_docs/psg/PSG_212194.pdf



<https://www.bachem.com/news/galnac-delivering-promise-of-oligonucleotides/>

30 nt DNA aptamer in interaction with testosterone

Ascend Evo 1.0 GHz



Data courtesy
Jose Martins (Ghent University)
Miquel Pons & Margarida Gairí
(University of Barcelona)



© 2022 Bruker

2D NOESY 100 ms mixing time 200 nmole of DNA (0.4 mM, 1.9 mg)




03 Antibodies

Biologics: mAbs (monoclonal antibodies)

Goal is to obtain an accurate fingerprint of the mAb (HOS) and analyze both drug substance and drug product

The challenge

- MW ca. 150 kDa for a full antibody
- Fast T2 relaxation
- Natural abundance only !
- Acquire mainly 1D ^1H spectra and 2D $^1\text{H}/^{13}\text{C}$ spectra of CH_3 groups



<https://en.wikipedia.org/wiki/Antibody>


mAb fingerprinting by NMR at very high field 1GHz data acquired with a TCI cryoprobe

Sheep IgG


- IgG1 sample MW 150 KDa
- 286 μ M in NIST Buffer
- Temperature 323 K
- Natural abundance

Data courtesy

Miquel Pons & Margarida Gairí
(University of Barcelona)



1D Spectrum showing all the ^1H of the mAb
Experiment time 5 min



$^1\text{H}/^{13}\text{C}$ correlation: Each peak represents a different methyl group of the mAb (Amino acids observed: Met, Ile, Thr, Ala, Val, Leu)
Experiment time 7 hours

Higher Order Structure (HOS) analysis for mAbs studies Guidelines & acquisition parameters



Parameters guidelines documents: sample preparation, data acquisition and data analysis

1D & 2D (HOS*) Topspin acquisition parameter sets (from TS 4.1.4)

Filter by file names: Exclude: Clear

Class = Dim = Show Recommended

Type = SubType = SubTypeB =

HOS_1DNOESY	HOS_1DSTE	HOS_AFHMQC	HOS_CPMGES1D	HOS_HSQC
HOS_T1RES1D	HOS_T1RESPE1D	HOS_XLAFHMQC	HOS_ZGESGP	HOS_ZGESGPPE

Analysis by BioHOS package (Mestrelab Research)

1D & 2D spectra analysis, link between NMR spectra and statistical tools (including bucketing)



1D Analysis (PCA, PLS, SIMCA & PROFILE)

2D Analysis (ECHOS, CCSD, PCA, PLS & SIMCA)




NMR spectroscopy as a characterization tool enabling biologics formulation development

Junhe Ma^{a,*}, Charles Pathirana^a, David Q. Liu^b, Scott A. Miller^a

^a Chemical Process Development, Bristol Myers Squibb Company, 1 Squibb Drive, New Brunswick, NJ 08903, United States

^b Vertex Pharmaceuticals, 50 Northern Avenue, Boston, MA 02210, United States



04 Bio-Production

From bioprocess understanding to control with the **Fourier 80**

Brings the detector to the (bio)reactor. Monitors 'where it happens as it happens'

Online monitoring by NMR

- Quantitative by definition
- Calibrates other techniques
- On-the-fly quantitative build-up curves
- Increased selectivity, resolution and sensitivity compared to vibrational spectroscopy
- Reduced risk - direct knowledge transfer from mid/high-field NMR to low-field benchtop

InsightMR for **Fourier 80**



Control




InsightMR for mid/high field NMR



Understanding


Integrating NMR into Process Analytical Technology with **synTQ**

- Enabling simultaneous acquisition of data from different instruments with on-the-fly data visualisation
- Data driven results for increase understanding, increase quality, decrease risk and cost
- Integration with system controller (e.g., Temp. and pH in bioreactors) facilitating feedback control
- Complete data journey integrating with MVA models, historian, MES, EKM, LIMS



PEAXACT Automating data analysis with hard modelling

- Implementation of PEAXACT module in synTQ software for data analysis
- 7 days bioproduction of the IgG-like bispecific antibody
- Transfected mammalian cell using Expi CHO Expression system
- 3 different sample line
 - Control (no transfection); Transfected; Transfected with supplemented nutrients (from day 2)



05 PIPAc:
Production Intelligente de Principes Actifs


Fourier enabling API Production on Demand (PoD)

PIPAC Production Intelligente de Principes Actifs

- Intelligent production of API via creation of **“compact, mobile, autonomous and frugal production units”**
- flow chemistry** with **artificial intelligence** enabled process control using NMR and mIR as detectors

Artificial Intelligence


Continuous flow synthesis routes




Inline analytics & PAT

Additive manufacturing (3D printed reactors)

To bring flow chemistry to the edge, we designed an AI agent that learns to control the reaction



PIPAC






Acknowledgements

- BPIfrance: PIPAc Project
- UCB Pharma: Emmanuel Cornet and Cédric Schaefer
- Bruker BioSpin colleagues: Pharma Team, Software, Applications and R&D

Webinar “Rational Design of a New Generation of Therapeutic Oligonucleotide Tools” – 5 December, 17:00 CET



Presenting new chemically modified oligonucleotide tools with the potential to tackle multifactorial disorders (among them, cancer and related drug resistance issues), while overcoming the limitations of oligonucleotide therapeutics.



Dr. Montserrat Terrazas

Associate Professor Department of Inorganic and Organic Chemistry (Organic Chemistry Section), Faculty of Chemistry, University of Barcelona, Barcelona, Spain





Innovation with Integrity