

● InsightMR

The Solution for online Process Monitoring
-NMR在线监控化学反应过程的解决方案



吕娟

Bruker Webinar, August 2016



Content



- Why reaction monitoring?
为什么要进行反应监测?
- Reaction monitoring by NMR
用NMR来进行反应监测
- InsightMR
 - Software part 软件功能介绍
 - Hardware part 硬件配置介绍
- Conclusions 总结



Lab



The plant

Why Reaction Monitoring



- Purpose:

- Reaction completion determination 反应完成
- In-situ yield determination 原位产率
- **Reaction understanding**: mechanism 反应机制
- Extracting kinetic parameters from the time course data: activation energy, rate constant
提取动力学参数：活化能，速度常数（物化）
- **Kinetic modelling** 动力学模拟
- Process optimisation 工艺优化



- Applications:

- Chemistry (e.g. catalysis) 化学
- Stability, degradation studies, crystallisation 稳定性研究，降解研究，晶体化研究
- Dissolution 分解
- Bio transformations, fermentation 生物转化，发酵

- Classic method – gained popularity recently.
e.g. PAT – process control







Reaction Monitoring by NMR

Advantages



- Quantitative by default: no need to calculate response factors 定量
- Information rich: possible to elucidate the structures of starting materials, products, and intermediates 定性



Reaction Monitoring by NMR

Type of measurements



- **At-line** sampling then measuring in a tube (or transfer)
- **In-line**
 - monitoring reaction in an NMR tube
 - problem*: dead-time between mixing reagents(outside the tube)and acquisition. Not good for fast reactions
 - not always real conditions(e.g. temperature)
 - reagents flow and mix just before entering the NMR probe or in the probe
- **On-line** reaction vessel next to the magnet. Sample flows out and then back

Special issue research article

Received: 19 August 2014

Revised: 11 March 2015

Accepted: 15 April 2015

Published online in Wiley Online Library

(wileyonlinelibrary.com) DOI 10.1002/mrc.4259

Reaction monitoring using online vs tube NMR spectroscopy: seriously different results

David A. Foley,^{a*} Anna L. Dunn^{a,b} and Mark T. Zell^a

We report findings from the qualitative evaluation of nuclear magnetic resonance (NMR) reaction monitoring techniques of how each relates to the kinetic profile of a reaction process. The study highlights key reaction rate differences observed between the various NMR reaction monitoring methods investigated: online NMR, static NMR tubes, and periodic inversion of NMR tubes. The analysis of three reaction processes reveals that rates derived from NMR analysis are highly dependent on monitoring method. These findings indicate that users must be aware of the effect of their monitoring method upon the kinetic rate data derived from NMR analysis. Copyright © 2015 John Wiley & Sons, Ltd.

Keywords: NMR reaction monitoring; online NMR; reaction mechanism




Introduction

Nuclear magnetic resonance (NMR) is an extremely powerful tool for the analysis of reaction mixtures. Not only can NMR be used as a quantitative method of monitoring reaction processes, it also has the advantage of providing detailed structural information

vessel allows replication of reaction conditions while allowing for analysis in an essentially unperturbed state. There are examples of online NMR reaction monitoring being applied at both high and low fields.^[6]

The third method, stopped-flow NMR (iii), is typically used for the detection of rapid kinetics (within 25–100 ms of reaction time) in


NMR



Flow vs Tube:

- ✓ Better captures the start of the reaction
- ✓ Mimics real reaction conditions: stirring, Te, P
- ✓ Enables the simultaneous acquisition IR, pH, MS

IR





Reaction Monitoring by NMR

Unparalleled Information for Reaction Understanding

JOC The Journal of Organic Chemistry

Article

pubs.ac.org/joc

ReactNMR and ReactIR as Reaction Monitoring and Mechanistic Elucidation Tools: The NCS Mediated Cascade Reaction of α -Thioamides to α -Thio- β -chloroacrylamides

David A. Foley,¹ Christopher W. Doecke,¹ Jonas Y. Baser,¹ Jeremy M. Merritt,¹ Linda Murphy,² Marie Kissane,² Stuart G. Collins,² and Stuart G. Collins,² and Stuart G. Collins,²

MAGNETIC
Magn. Reson. C
Published online
(www.interscience.com)

Optim
indus

Michael A
AstraZeneca P

TOP TE
Analytical Chem

analytical
chemistry

Letter

pubs.ac.org/ac

Online NMR and HPLC as a Reaction Monitoring Platform for Pharmaceutical Process Development

David A. Foley,* Jian Wang, Brent Maranzano, Mark T. Zell, Brian L. Marquez, Yanqiao Xiang, and George L. Reid

Organic Process
Research &
Development

COMMUNICATION

pubs.acs.org/OPRD

A Safe and Efficient Synthetic Route to a 2,5-Dimethyl-1-aryl-1H-imidazole Intermediate

Ana C. Barrios Sosa,* R. Thomas Williamson, Ryan Conway, Ashish Shankar, Rosline Sumpter, and

What's missing?

1. Integration and feed-back: on-the-fly acquisition control based on real-time data processing and analysis
积分和反馈信息: 基于实时数据处理和分析的动态采样控制
2. Hardware for online monitoring in real-time under real condition
真实条件下进行在线监测的硬件

Received: 16 February 2010

Revised: 19 March 2010

Accepted: 25 March 2010

Published online in Wiley InterScience: 29 April 2010

(www.interscience.com) DOI 10.1002/mrc.2610

A simple flowcell for reaction monitoring by NMR

M. Khajeh,^a M. A. Bernstein^b and G. A. Morris^{a*}

Reaction Monitoring




InsightMR

The Solution for Process Monitoring




*InsightMR flow tube**



*InsightMR software:
TopSpin+IconNMR+Dynamic center
(TS 3.5 Windows)*



InsightMR Software



InsightMR software part

Key features




- A single interface for automated acquisition control, interactive processing and analysis, resulting in real-time kinetic profiles
自动采样控制、交互式处理和分析的单一界面，可完成实时动力学轮廓
- Acquisition and real-time analysis of a series of 1D NMR spectra, using different nuclei and interleaved experiments
支持对一系列采用不同核和隔行扫描的1D NMR谱图进行采样和实时分析
- On-the-fly acquisition parameters adjustment based on real-time kinetic data
基于实时数据处理和动力学轮廓计算，进行即时采样参数调整
- Simultaneous monitoring of multiple reactions at the same time using parallel acquisition and analysis capabilities
运用平行采样和分析的功能在同一时间对多个采样进行同步监控
- Default kinetics parameters provided for both deuterated and non-deuterated solvents
提供默认的动力学参数设置，能够容易地设置实验来观测氘代溶剂或非氘代溶剂中的反应过程

InsightMR software part




- A single platform for acquisition control, interactive processing and data analysis


Click Analysis



A Quick Demonstration Simple Diels Alder



Key Message: Fully integrated acquisition, data processing and analysis





● **InsightMR Software 1.0**

Functionality Overview





Answers to Chemical Questions

Concentration Profiles



Concentration Profiles




Multiple Samples Run in Parallel

The screenshot displays the InsightMR software interface, version 1.0.b14 (2015 Mar/03). The interface is divided into several sections:

- Method Center:** A sidebar on the left containing 'InsightMR', 'Window Control', and 'Kinetics'.
- Sample Configuration:** A panel with 'Sample Position' set to 1, 'Name' (empty), and 'Solvent' options including Acetone, D2O, CD3CN, DMSO, EtOD, Tol, MeOD, C6D6, and CDCB.
- Experiment Configuration:** A panel with radio buttons for 'K_PROTON', 'K_CARON', 'K_NOESY1D', and 'K_WATER'.
- Reaction Monitoring:** A panel with controls for 'Modify Acquisition Parameters', 'Number of Experiments' (set to 10), 'Delay between Experiments (H:M:S)' (set to 00:00:10), 'Use Start Time' (checked, set to 22:12:43 Thu 15), 'Start Reaction Stopwatch', and 'Stopwatch time (seconds)' (set to 0).
- Queue:** A panel with an 'Update from selected' checkbox and a 'Queue Experiments' button.
- Experiment Table:** A table with columns: Hol..., Type, Status, Disk, Name, N..., Solvent, Experiment, Pri, Par, Title/Orig, Time, User, Start Time. It shows five rows, all with 'Available' status.
- Buttons:** A row of buttons including 'Submit', 'Cancel', 'Edit', 'Delete', 'Add', 'Copy', and 'Change User'.
- Preceding Experiments:** A table with columns: #, Date, Holder, Name, No., Experiment, Load, ATM, Rotation, Lock, Shim, Acq, Proc, User, Disk, Title/Orig. It is currently empty.
- Status Bar:** At the bottom, it shows 'Search Preceding', a search icon, and system information: 'BACS/SampleCase/Pro (RS232) | Busy until: No Jobs! | Day Experiments: 00:00 | Night Experiments: 00:00 | User: ReactionMonitoring'.

Interleaved Experiments



InsightMR, Version: 1.0.b15 (2015 Mar/10)

Help

Method Center
 InsightMR
 Window Control
 Kinetics

Start

Sample
 Sample Position: 1
 Name:
 Solvent:
 Acetone D2O
 CD3CN DMSO
 EtOD Tol MeOD
 C6D6 CDCB3

Experiment
 K_PROTON K_NOESY1D K_WATER
 K_C13IG

Title

Reaction Monitoring
 Modify Acquisition Parameters
 Number of Experiments: 10
 Delay between Experiments (H:M:S): 00:00:10
 Use Start Time: 13:17:42 Thu Mar 12 2016
 Start Reaction Stopwatch
 Stopwatch time (seconds): 0

Queue
 Update from selected
 Queue Experiments

Experiment Table

Hol...	Type	Status	Disk	Name	No.	Solvent	Experiment	Pri	Par	Title/Orig	Time	User	Start Time
▶ 1		Available											
▶ 2		Available											
▶ 3		Available											
▶ 4		Available											
▶ 5		Available											
▶ 6		Available											
▶ 7		Available											
▶ 8		Available											
▶ 9		Available											
▶ 10		Available											
▶ 11		Available											
▶ 12		Available											
▶ 13		Available											
▶ 14		Available											
▶ 15		Available											
▶ 16		Available											

Submit Cancel Edit Delete Add 1 Cppy 1 Change User



Preceding Experiments

#	Date	Holder	Name	No.	Experiment	Load	ATM	Rotation	Lock	Shim	Acq	Proc	User	Disk	Title/Orig	Remarks

Search Preceding

BACS/SampleCase/Pro (RS232) Busy until: No Jobs! Day Experiments: 00:00 Night Experiments: 00:00 User: ReactionMonitoring

Interleaved Experiments






InsightMR Flow Tube



InsightMR hardware part

laboratory setup





InsightMR hardware part

Flow tube features

- Real-time monitoring (实时监测)
- Compatible with Bruker 5 mm probes (与Bruker 5mm 探头兼容)
- **Temperature control** transfer lines (控温传输线)
- Withstand **pressures over 10 bar** (承受压力超过10bar)
- **Interchangeable glass tube** (可更换玻璃管)
- Proven robustness - industry tested (工业测试证明稳定)
- Simultaneous acquisition of : NMR, IR, pH, MS... (可与其它分析仪器同时进行检测)
- Samples closer to the start of the reaction (样品接近反应的开始)



InsightMR hardware part

Demo lab UK



InsightMR hardware part

Demo lab UK




Flow interface

Enabling Parallel Acquisition of Data from Different Techniques




NMR Signal Stability Testing






Line Shape Test – 10 mM sucrose





Diels Alder Chemistry, NON-Deuterated Solvent



Diels Alder Chemistry, NON-Deuterated Solvent

Process Understanding – Structural Information of Intermediates



Conclusions



- Answers to key chemical questions: reaction yield, reaction kinetics
- Reaction understanding, identification of reaction intermediates and mechanistic information
- Rapid and straightforward generation of data to build kinetic models
- Intuitive – straightforward acquisition and processing workflow - makes NMR an accessible PAT tool for all audiences (no NMR experience required)
- Enables the facile use of NMR data to make strategic process chemistry decisions, ultimately leading to **cost savings**






Conclusions


1. InsightMR Software


- Download:
<https://www.bruker.com/service/support-upgrades/software-downloads/nmr.html>

- 3 months demo:
<https://www.bruker.com/service/support-upgrades/license-requests/insightmr-license-request-form.html>



2. InsightMR (Flow Tube + Software)





www.bruker.com

InsightMR@bruker.com

Q & A



- 是否有问题？
- 请您在Q&A面板中提交您的问题
- 我们做的如何？
- 当您退出webinar的时候, 请填写您对本次webinar的评价, 我们非常感谢您的反馈。

- 非常感谢！

Thank you!

Juan.lv@bruker.com

Lu.shan@bruker.com

Would you like to learn more? Contact a customer service representative.

