

What's new in MetaboScape 2024


One integrated solution for processing and interpreting of MS based Phenomics, non-targeted Metabolomics, Lipidomics and MALDI Imaging data



Outline

- 01 Combine Confidence and Speed using the Sample Characterization Workflow
- 02 Improved Accuracy with CCS-Predict Pro 2024 & Prediction for Single Structures & CCS-Predict via REST API for custom use
- 03 Use the REST API to create MetaboScape projects from your favourite programming language
- 04 T-ReX Improvements for even higher quality


Combine Confidence and Throughput




The Sample Characterization Table-based Annotation workflow allows separation of highly confident annotation and highly confident quantitation/high throughput.

- Create a **Feature Table A** from **pooled samples** acquired with a high-dimensional mode like LC-TIMS-PASEF, for the sake of sample characterization. Use the CCS and MSMS information for highly confident annotation.
- Create a second **Feature Table B** from **individual samples**, acquired with e.g. LC-MS.
- Use the highly confident annotations from **Feature Table A** to annotate **Feature Table B**. This results in annotations in **Feature Table B**, referencing the CCS and MSMS annotation scores from **Feature Table A**.

Combine Confidence and Throughput



Review annotation criteria from Feature Table A without leaving Feature Table B:




MS/MS spectra, EICs or EIMs that are references to Feature Table A are painted in turquoise and decorated with an SC icon, that allows to quickly navigate to the referenced Feature Table and Feature.

If MS/MS spectra, EICs or EIMs are present in Feature Table B, they will not be overwritten, but displayed as found in Feature Table B.


Annotate with CCS-Predict Pro 2024 confidence!

The CCS-Predict Pro models for $[M+H]^+$, $[M-H]^-$, and $[M+Na]^+$ have been trained with larger reference data sets, which have grown by 90%, 72% and 84% respectively, compared to the previous models.

Class-specific CCS prediction errors for the old and the new CCS-Predict Pro model for $[M+H]^+$. Predictions were matched against references from both, training and test data.



 **CCS-Predict Pro**



 **CCS-Predict Pro 2024**

*Only six reference structures for compound class OrgHa

- Improvement in predictive accuracy across all main compound classes
- The three models for $[M+H]^+$, $[M+Na]^+$ and $[M-H]^-$ benefit from increased training data.

CCS-Predict Pro 2024 fully integrated into your workflow

MetaboScape CCS-enables annotation with

- TL Target List
- SL Spectral Library
- BT Biotransformer

if structural information is available

MS/MS libraries without CCS values

HMDB Metabolite Library
BRUKER Daltonics


NIST Library
BRUKER Daltonics

MetaboBASE Personal Library
BRUKER Daltonics


MS-DIAL*

* MS-DIAL is not a Bruker product

MS/MS reference



compound structure



Met
SL + CCS-Predict Pro

Name	Annotations	AQ	AQ Details	MS/MS score	ΔCCS [%]
Caffeine	SL TL			952.5	-0.1


CCS-enabled MS/MS library annotation based on CCS prediction

Try it out!

Easy CCS-Predict with drag-and-drop

In addition to the integration into automated workflows, MetaboScape's CCS-Predict in the Compound ID perspective allows to predict CCS values for structures via drag-and-drop. Choose from the CCS-Predict Pro 2024 and the CCS-Predict 4D-Lipidomics™ model.

InChI=1S/C6H12O6/c7-2-6(11)5(10)4(9)3(8)1-12-6/h3-5,7-11H,1-2H2/t3-,4+,5-,6+/m0/s1
 InChI=1S/C6H12O6/c7-1-2-3(8)4(9)5(10)6(11)12-2/h2-11H,1H2/t2-,3+,4+,5-,6+/m1/s1
 InChI=1S/C6H12O6/c7-1-2-3(8)4(9)5(10)6(11)12-2/h2-11H,1H2/t2-,3-,4-,5-,6-/m0/s1
 InChI=1S/C6H12O6/c7-1-3-4(9)5(10)6(11,2-8)12-3/h3-5,7-11H,1-2H2/t3-,4-,5+,6-/m1/s1
 InChI=1S/C6H12O6/c7-1-2(8)4(10)6(12)5(11)3(1)9/h1-12H/t1-,2-,3-,4+,5-,6+
 InChI=1S/C6H12O6/c7-1-3-4(9)5(10)6(11,2-8)12-3/h3-5,7-11H,1-2H2/t3-,4-,5+,6-/m1/s1
 InChI=1S/C6H12O6/c7-1-2-3(8)4(9)5(10)6(11)12-2/h2-11H,1H2/t2-,3-,4+,5-,6+/m1/s1
 InChI=1S/C6H12O6/c7-1-3(9)5(11)6(12)4(10)2-8/h3,5-9,11-12H,1-2H2/t3-,5+,6-/m1/s1
 InChI=1S/C6H12O6/c7-1-2(8)4(10)6(12)5(11)3(1)9/h1-12H/t1-,2-,3-,4+,5-,6+
 InChI=1S/C6H12O6/c7-1-2-3(8)4(9)5(10)6(11)12-2/h2-11H,1H2/t2-,3-,4+,5-,6-/m1/s1



Name	Molecular Form...	[M+Na]+ Prediction	Δ CCS [%]	Database	Ion Notation	Measured	Prediction	Δ CCS [%]
C6H12O6	C ₆ H ₁₂ O ₆	147.8	-1.2	manual (CCS-P	[M+H] ⁺	-	138.0	-
C6H12O6	C ₆ H ₁₂ O ₆	146.0	-0.0	manual (CCS-P	[M+Na] ⁺	146.0	146.0	-0.0
C6H12O6	C ₆ H ₁₂ O ₆	145.3	0.5	manual (CCS-P	[M-H] ⁻	-	119.8	-
C6H12O6	C ₆ H ₁₂ O ₆	144.7	0.9	manual (CCS-P				
C6H12O6	C ₆ H ₁₂ O ₆	144.7	0.9	manual (CCS-P				
C6H12O6	C ₆ H ₁₂ O ₆	144.7	0.9	manual (CCS-P				
C6H12O6	C ₆ H ₁₂ O ₆	144.7	0.9	manual (CCS-P				
C6H12O6	C ₆ H ₁₂ O ₆	142.0	2.8	manual (CCS-P				

Drag and drop into CCS-Predict:



InChI
text



SMILES
text



.mol
file



.sdf
file


Use CCS-Predict Pro in your own code

MetaboScape's CCS-Predict models now are accessible via the REST API via new endpoints.

GET /ccspredict/models

POST /ccspredict/models/{modelId}

Your code (for example with python):



```
models = ccspredict.list_ccs_predict_models()
...

structure = metaboscape_client.Structure("SMILES", "CC(C)C(CCCN(C)CCC1=CC(=C(C=C1)OC)OC)(C#N)C2=CC(=C(C=C2)OC)OC")
ion_predictions = ccspredict.predict_ccs("CCS-Predict_Pro_2024", body=structure)
for ion in ion_predictions:
    print(ion.ion_notation, ': ', ion.predicted_ccs)

[M+H]+ : 206.58969106016252
[M-H]- : 222.53750677292976
[M+Na]+ : 226.36295962368797

structure = metaboscape_client.Structure("INCHI", "InChI=1S/C27H38N2O4/c1-20(2)27(19-28,22-10-12-24(31-5)26(18-22)33)
ion_predictions = ccspredict.predict_ccs("CCS-Predict_Pro_2024", body=structure)
for ion in ion_predictions:
    print(ion.ion_notation, ': ', ion.predicted_ccs)

[M+H]+ : 206.58969106022255
[M-H]- : 222.5375067728079
[M+Na]+ : 226.36295962364431
```


Connect your Lab to MetaboScope

Connect your data repository to MetaboScope with custom code, using the REST API.

POST /projects/create

GET /projects/taskstatus

POST /experiments/create


POST /featuretable/create



Improved Feature Finding Workflow



- Improved retention time alignment
- Improved merging of MSMS spectral peaks across scans
- Improved determination of peak start and end in T-ReX 3D
- If required, use calibration results created with Data Analysis



Innovation with Integrity