




UPPSALA
UNIVERSITET

Chemical Metabolomics - Advanced chemoselective probe-based mass spectrometric investigation of the human carbonyl-metabolome




Alejandro Torregrosa Chinillach^[1]; Weifeng Lin^[1]; Maria Theodorou^[1]; Wawrzyniec Haberek^[1]; Nikolas Kessler^[2]; Sofie Weinkouff^[2]; Aiko Barsch^[2]; Cristian De Gobba^[2]; Matthew R. Lewis^[2]; and Daniel Globisch^[1]


^[1] Department of Chemistry, Biomedical Center, Uppsala University, Sweden; ^[2] Bruker Daltonics GmbH & Co. KG, Bremen, Germany
Daniel.globisch@kemi.uu.se, Homepage: <https://www.globischlab.org> Follow us on Bluesky: @GlobischLab

New approaches are required for the discovery of unknown biomarkers. Carbonyl-containing metabolites are a highly reactive compound class that can form conjugates with DNA and proteins. They are also involved in a wide range of biochemical pathways and their dysregulation has been linked to pathological conditions as well as linked to microbiome-metabolism. Six carbonyl-containing metabolites have been classified as carcinogens by the International Agency for Research of Cancer. We have recently developed a new methodology based on chemoselective probes to investigate this compound class after derivatization at enhanced mass spectrometric sensitivity due to removal of the captured metabolites from the matrix.¹⁻⁴ To improve the current bottleneck of metabolomics analyses, we have enhanced our analysis with the timsTOF technology for structure elucidation of metabolites through in-silico derivatization, CCS prediction and annotation. We have utilized our Quantitative Sensitive CHEmoselective MetAbolomics (quant-SCHEMA) methodology for the investigation of biological samples.¹⁻³ This method allows for the simultaneous and semi-quantitative analysis at the femtomole level as well as qualitative analysis at attomole quantities that allows for detection of more than 200 metabolites in diverse human samples. To improve our analysis and metabolite structure elucidation workflow, we treated NIST[®] SRM[®] 1950 plasma samples with the quant-SCHEMA method and analyzed the derivatized carbonyl metabolites with a UHPLC coupled to DDA-TIMS-MS/MS. The data was acquired on a timsTOF instrument (Bruker) in positive ionization mode. Data was processed with a development version of MetaboScape 2026, using its in-silico derivatization workflow to automatically annotate and identify the captured compounds.



Chemoselective Probe Design, Workflow and Results^[1-7]



Dietary intervention study^[3]




In-silico derivatization with CCS-Predict Pro



Name	CCS	Neutral Formula	CAS	HMDB
2-Ethylbutyraldehyde	205.1	C10H20O	97-96-1	HMDB0012320
Isomonic acid	217.0	C20H38O4	694-38-8	HMDB0012797
trans-Cinnamaldehyde	217.0	C9H8O	104-90-2	HMDB0000441
Isotonic acid	197.4	C18H34O4	133-90-2	HMDB0000730
Naringenin	221.2	C15H12O5	480-41-1	HMDB0002670
2-Octanone	214.4	C8H16O	111-51-7	HMDB0001234
3-Methyl-2-oxovaleric acid	206.1	C7H12O4	14809-86-3	HMDB0000461
Sulfasal	218.0	C12H10N2O4	136-26-7	HMDB0000261
Pyruvic acid	186.7	C3H4O3	137-13-3	HMDB0000343
trans-2-Heptanal	208.4	C7H14O	6728-26-3	HMDB0001496
4-Hydroxybenzaldehyde	209.1	C7H6O3	123-02-0	HMDB0012179
Ethylpyruvate	193.0	C6H10O4	96-26-4	HMDB0001892
2,4-Dihydroxyacetophenone	205.5	C8H8O4	99-84-9	HMDB0002639
2-Deoxy-D-ribose	205.0	C5H10O4	533-50-5	HMDB0001099
Benzaldehyde	204.3	C7H6O	100-52-7	HMDB0001115
Raspberry ketone	209.1	C10H18O	5471-51-2	HMDB0003729
o-Tolualdehyde	208.5	C8H8O	100-24-4	HMDB0002626
Acetone	188.8	C3H6O	67-64-1	HMDB0001039
p-Tolualdehyde	209.9	C8H8O	104-67-0	HMDB0002638
	209.4	C8H8O	107-88-0	HMDB0002634

CCS reference library for >90 metabolites created for targeted annotation at a high confidence level

Carbonyl-Chemprobe Standards timsMetabo vs. tims QTOF



References

- W. Lin, L.P. Conway, M. Vujanovic, J.-M. Löhr, D. Globisch, *Angew. Chem. Int. Ed.* **2021**, *60*, 23232-23240.
- Lin, Gerullat, Braadland, Fournier, Hov, Globisch *Angew. Chem. Int. Ed.* **2024**, *63*, e202318579.
- Lin, Mellinghaus, Rodriguez-Mateos, Globisch *Food Chem.* **2023**, *425*, 136481.
- Lin, Conway, Block, Sommi, Vujanovic, Corra, Globisch *Analyst* **2020**, *145*, 3822.
- N. Garg, L.P. Conway, C. Ballet, M.S.P. Correia, F. Olsson, M. Vujanovic, J.-M. Löhr, D. Globisch, *Angew. Chem. Int. Ed.* **2018**, *57*, 13805-13809.
- A. Kaur, W. Lin, V. Dohvalyuk, L. Driutti, M.L. Di Martino, M. Vujanovic, J.-M. Löhr, M. Sellin, D. Globisch, *Chem. Sci.* **2023**, *14*, 5291-5301.
- W. Lin, F. Romero Garcia, E.L. Norin, D. Kart, L. Engstrand, J. Du, D. Globisch, *Chem. Comm.* **2023**, *59*, 5843-5846.

Acknowledgements

