


Field based multiplex detection of biotreat agents



15th Medical Biodefense Conference

Christopher Pöhlmann and Thomas Elssner
 Bruker Daltonik GmbH, Permoserstrasse
 15, 04318 Leipzig, Germany
 C.Poehlmann@bruker.com

Introduction

Due to the potential risk of terrorist attacks using highly toxic or pathogenic biotreat agents (bacteria, viruses and toxins) there is a need for rapid, sensitive and specific detection of all kinds of biotreat agents in the field. In general, nucleic acid based detection platforms are applied for pathogen detection. However, PCR is unable to detect proteinaceous toxins and requires a clean sample demanding elaborate sample preparation. In contrast, immunoassay based technologies offer the possibility for universal detection of bacteria, viruses and toxins.

Methods

pBDi platform:

- Biochip based detection platform
- Integrated in a robust suitcase for use in the field by first responders (Fig. 1 A)
- Based on an automated electrochemical sandwich-ELISA (Fig. 1 B)
- Capture antibodies against specific biotreat agents are immobilized on electrodes of a biochip (16 interdigitated gold electrodes)
- Multiplex detection of up to 6 biotreat agents within 20 minutes
- Inherent positive and negative control
- Very sensitive due to strongly amplified electrochemical detection system including a redox cycling process of enzyme product between interdigitated electrode structure (Fig. 1 B)
- Fully automated fluidics
- Software controlled data analysis applying a ruggedized tablet PC via USB or bluetooth connection

Biochips: Three different biochip panels detecting each up to six biotreat agents in parallel including inherent positive and negative control electrodes are available (Tab. 1). Furthermore, biochips for site-specific immobilization of customer antibodies are available.




Fig. 1 The pBDi platform.

Biotreat agent	LOD ¹	Chip
SEA	0.25 ng mL ⁻¹	Toxin chip II
SEB	0.1 ng mL ⁻¹	Toxin chip I
BoNT/A	0.75 ng mL ⁻¹	Toxin chip I
BoNT/B	0.5 ng mL ⁻¹	Toxin chip I
BoNT/C	10.0 ng mL ⁻¹	Toxin chip II
BoNT/D	2.5 ng mL ⁻¹	Toxin chip II
BoNT/E	17.5 ng mL ⁻¹	Toxin chip II
BoNT/F	1.0 ng mL ⁻¹	Toxin chip I
Ricin	0.5 ng mL ⁻¹	Toxin chip I
Abrin	1.0 ng mL ⁻¹	Toxin chip II
<i>B. anthracis</i> (spores)	1x10 ⁵ cfu mL ⁻¹	Biothreat chip
<i>Y. pestis</i>	5x10 ⁴ cfu mL ⁻¹	Biothreat chip
<i>F. tularensis</i>	1x10 ³ cfu mL ⁻¹	Biothreat chip
<i>Brucella spp.</i>	1x10 ⁴ cfu mL ⁻¹	Biothreat chip
<i>Burkholderia mallei</i>	1x10 ⁴ cfu mL ⁻¹	Biothreat chip
Vaccinia virus	9x10 ⁴ pfu mL ⁻¹	Biothreat chip

Tab. 1 Detection limits obtained with pBDi platform. ¹LOD = Blank + 3 x SD

Results

- pBDi is usable in the field by first responders using personal protective equipment (Fig. 2).
- pBDi was evaluated by German and French fire fighters using 20 unknown samples containing varying amounts of toxin (C-terminal part of the heavy chain of BoNT/C) or bacteria (inactivated *F. tularensis*) in buffer or various sample matrices, respectively. pBDi identified all 20 unknown samples correctly.
- pBDi shows excellent limit of detections (LODs) for rapid multiplex detection of bacteria, *orthopoxvirus* and toxins (Tab. 1).
- Specificity of biochip detection is demonstrated applying different amounts of bacteria or vaccinia virus, respectively (Fig. 3).
- Demonstration of assay robustness detecting ricin or *F. tularensis*, respectively, in sample matrices with high recovery rate (Fig. 4).



Fig. 2 End user test in course of GEFREASE project. pBDi operated with personal protective equipment by French and German fire fighters.




Fig. 3 Specificity of microorganism detection.




Fig. 4 (A) Detection of 10 ng mL⁻¹ ricin in different sample matrices. (B) Detection of 5x10⁵ cfu mL⁻¹ *F. tularensis* in various sample matrices.

Summary

The results demonstrate the potential of electrochemical biochips for sensitive multiplex on-site detection of biotreat agents applying a fully automated procedure. The assays are simple to perform, show no cross-reactivity with tested agents, are rapid, and require only minimal sample preparation.

Acknowledgments

Parts of the research is funded by the Federal Ministry of Education and Research (BMBF) joint research project GEFREASE, support code 13N12222.

Conclusions

- Biotreat agents analysis in approx. 20 minutes with an automated procedure
- Parallel detection of up to six biotreat agents including positive and negative control
- Portable and vehicle integrable due to specific suitcase design
- Excellent sensitivity
- Extendable to new agents

pBDi