

Application Note FTMS-48

Microbial Imaging Mass Spectrometry with Fourier Transform Ionization Mass Spectrometry

Abstract

In this application note we describe a straightforward approach for the preparation and molecular analysis of a bacterial-fungal interaction and the visualization of metabolic exchange factors using an FT-ICR high resolution mass spectrometer.

Introduction

Since the discovery of penicillin it has been widely recognized that microbes produce biologically active molecules, some of which have been found beneficial for treating infections (e.g. penicillin, vancomycin and daptomycin), enabling transplantations (e.g. rapamycin) and lowering cholesterol (e.g. lovastatin). Despite the importance of these molecules in agriculture, diagnostic and therapeutic areas, there are few methods that capture molecular information directly from microbial cultures.

In 2009, microbial imaging mass spectrometry was developed as a way to capture the molecular information from cultured microbes directly from agar surfaces¹. Microbial imaging mass spectrometry is an extension of the field of histological imaging mass spectrometry² and has resulted in the capture of otherwise hidden molecular microbial phenotypes and metabolic exchange of compounds that promise new therapeutic and biochemical

Authors

Vanessa Phelan¹, Pieter C. Dorrestein¹,
Shannon Cornett²
University of California San Diego¹, Bruker Daltonics²

Keywords	Instrumentation and Software
MALDI Imaging	solariX
SmartFormula	
Bacteria	
Microbiology	

activities^{1, 2-13}. Often however, these bioactive compounds can be small-molecules that are isobaric with intense matrix ions or other endogenous compounds. When imaging with MALDI-TOF systems it is often useful to supplement MALDI imaging analyses with other techniques such as TOF/TOF imaging or LC-MS in order to fully characterize molecular structures and/or molecular distributions. However these approaches are not without disadvantages. MALDI-TOF/TOF imaging provides high molecular specificity which lends itself quite well to verification of a suspected molecule but is limited for discovery experiments. On the other hand, LC-MS strategies provide access to a wide dynamic range of compounds but at the sacrifice of enriching specific molecules via extraction and the LC procedure itself, as many molecules are difficult to extract and/or chromatograph poorly, as well as missing the spatial information.

Although MALDI-TOF based imaging is incredibly useful for microbiology, Bruker's solariX high-performance mass spectrometer offers an enhanced strategy for MALDI imaging that is well suited to studying microbial metabolites and other small molecules¹². It's high mass-resolution capabilities in MS-mode facilitates discovery studies by separating metabolite images from the interferences presented by ions that are very close in mass. At the same time, solariX's high measurement accuracy can often provide unambiguous molecular formula assignment for each of these resolved images. The same resolution and accuracy capabilities are also present in MS/MS mode if more detailed structural confirmation is needed. In this note we show the power of solariX for metabolic imaging of a *Pseudomonas aeruginosa* and *Aspergillus fumigatus* interaction, two organisms found in cystic fibrosis lungs.³

Methods


Detailed descriptions of the microbial sample preparation can be found in the references and several minor variations exist. Typically, 10-11 ml of 1-2% agar with necessary nutrients was poured into a Petri-dish to form a layer of agar ~1mm thick. *Pseudomonas aeruginosa* and *Aspergillus fumigatus* were subsequently added to the agar and incubated until colonies began to interact. During the incubation period it is important to seal the Petri-dish in order to minimize evaporation. When the microbes reached the desired phenotype the agar was cut and peeled from the Petri-dish and placed on top of an unpolished MALDI target plate. A digital photo was taken of the sample before matrix was applied. While the sample was wet, it was 'dusted' with crushed crystals of 50:50 DHB:HCCA (Sigma-Aldrich #50149-1G-F) from a manually-agitated 53 μm sieve. After the sample became visibly saturated with matrix, it was dried at 370C for 1-3 hours. As it dried the agar shrank from 1 mm to ~120 μm thick. Residual matrix

was then removed using a gentle stream of compressed air. It is important to note that the plate can be analyzed immediately or stored in a desiccator. Although molecular decomposition during prolonged storage is possible, in our hands, samples stored in a desiccator for two weeks still produce strong ion signal with minimal degradation. A video of the sample preparation is available on youtube (<http://youtu.be/RIQ4BtjSCvI>) and the experimentalist wanting to perform these experiments should read ref 11.

MALDI imaging was performed using a 9.4T solariX equipped with dual ESI/MALDI ion source and Smartbeam II laser set to minimum focus. The m/z range acquired was 200 – 360 with a transient length of 700ms for a broad band resolving power of ~150,000. Single-scan spectra consisting of 500 accumulated laser shots were acquired at a rate of <1s/spectrum and calibrated using the DHB cluster ion at m/z 273.03936 as lock mass before saving. Data were post-processed using a proprietary data reduction algorithm which eliminates baseline data points but keeps intact full peak information. Fleximaging software was used to register the digital photo of the sample plate to the X-Y stage of the instrument so that measurement regions could be drawn around the colonies directly on the imported digital image of the sample. This enabled the direct correlation of the chemistry of microbes with the phenotypes of the microbes and surrounding agar. For the present study, high spatial resolution wasn't a requirement and so MALDI images were acquired at a spatial resolution of 300 μm .

Results and discussions

When two microbes interact, metabolic exchange factors are secreted as a means of chemical communication resulting in synergy, commensalism and antagonism. Key to identifying and characterizing compounds with high activity is to map relationships between compound distribution and any underlying microbial response. Figure 1 illustrates the performance benefits of microbial imaging with solariX. The example shows a narrow segment at m/z 274 from the microbial interaction in which at least 18 distinct ion images are resolved. Five of these ions present clear differentiation to either *P. aeruginosa* or *A. fumigatus* as seen from the images shown in the figure. When investigating the activity of known target ions or discovering new compounds considerable molecular specificity is gained by separating isobaric ion images. The observation that some of the individual images are clearly specific to an individual microbial colony would otherwise go undetected in lower resolution systems. It is important to note here that ion desorption in solariX is decoupled from the mass measurement, unlike traditional MALDI-TOF, and we therefore observe minimal degradation of resolution or measurement accuracy from the non-conductive agar substrate.



Higher molecular specificity is accompanied by higher measurement accuracy which helps identify the molecular entity being imaged by allowing assignment of a molecular formula to each detected peak. For many metabolites, a measurement accuracy of 1 ppm or less is more than sufficient to determine a single molecular formula for an ion. An example of this is illustrated in Figure 1 using the ion images denoted as 1-5. Each of these five peaks was analyzed by SmartFormula to identify combinations of C, H, N and O would sum to within 5 ppm of the measured m/z and the results are presented in Table 1. With exception of ion 5, m/z 274.2075, a search window as large as 5 ppm returns a single molecular formula for each ion image. Enlarging the search window to 10 ppm returned a single formula, listed in Table 1 as ion 5a, but on examination, the difference of 10 ppm observed between the calculated and measured m/z for 5a is out of line with differences for peaks 1-4. One must consider the possibility that peak 5 does not represent a monoisotopic ion but rather one that contains 1 or more heavy isotopes. Indeed, the peak labelled M was measured at m/z 274.04249, a difference of 0.80 ppm from the calculated mass of the A+1 peak of [DHB-H₂O]₂+H⁺. When we included heavy isotopes in the potential formulas peak 5 it matches very well to the A+2 peak of nitrohydroquinone (Row 5b, Table 1) a known molecular product of *Pseudomonas*. The difference between the calculated and measured m/z is consistent

with peaks 1-4. Further, assignment to this compound is also supported by its image, with its distribution being localized to this colony.

peak	Meas. m/z	only formula in search win.	Calc m/z -meas m/z	Search window
1	274.1051	C ₁₁ H ₁₂ N ₇ O ₂	-1.57 ppm	5 ppm
2	274.1436	C ₁₆ H ₂₀ NO ₃	0.62 ppm	5 ppm
3	274.1800	C ₁₇ H ₂₄ NO ₂	0.53 ppm	5 ppm
4	274.1875	C ₁₁ H ₂₄ N ₅ O ₃	-0.34 ppm	5 ppm
5a	274.2075	C ₉ H ₂₄ N ₉ O	8.40 ppm	10 ppm
5b	274.2075	¹² C ₁₆ ¹³ C ₂ H ₂₆ NO	-0.36 ppm	5 ppm

Conclusion

Microbial imaging is a strategy that enables the chemical analysis of microbes grown on agar surfaces under controlled conditions. Using solarix to carry out MALDI-MS imaging, one can utilize its high resolving power to differentiate these species, especially in the lower mass range where the matrix signals and many of the relevant microbial metabolic exchange factors, virulence factors, and quorum sensors are found. In addition improving the specificity of images solarix can provide molecular formulas to ion images of interest. These capabilities can be invaluable tools for understanding the molecular discussions associated with microbial interactions.



References

- ¹Yang, Y.L., Xu, Y., Straight, P., Dorrestein, P.C., Translating metabolic exchange with imaging mass spectrometry. *Nature Chemical Biology*, 2009 Dec;5(12):885-7.
- ²Schwamborn K., Caprioli, R.M., Molecular imaging by mass spectrometry — looking beyond classical histology, *Nature Rev Cancer* 10, 639-646 (September 2010).
- ³Liu WT, Yang YL, Xu Y, Lamsa A, Haste NM, Yang JY, Ng J, Gonzalez D, Ellermeier CD, Straight PD, Pevzner PA, Pogliano J, Nizet V, Pogliano K, Dorrestein PC., Imaging mass spectrometry of intraspecies metabolic exchange revealed the cannibalistic factors of *Bacillus subtilis*., *Proc Natl Acad Sci U S A*. 2010 Sep 14;107(37):16286-90.
- ⁴Yang YL, Xu Y, Kersten RD, Liu WT, Meehan MJ, Moore BS, Bandeira N, Dorrestein PC., Connecting chemotypes and phenotypes of cultured marine microbial assemblages by imaging mass spectrometry. *Angew Chem Int Ed Engl*. 2011 Jun 20;50(26):5839-42.
- ⁵Gonzalez DJ, Haste NM, Hollands A, Fleming TC, Hamby M, Pogliano K, Nizet V, Dorrestein PC., Microbial competition between *Bacillus subtilis* and *Staphylococcus aureus* monitored by imaging mass spectrometry. *Microbiology*. 2011 Sep;157(Pt 9):2485-92.
- ⁶Watrous JD, Dorrestein PC., Imaging mass spectrometry in microbiology., *Nat Rev Microbiol*. 2011 Aug 8;9(9):683-94.
- ⁷Liu WT, Kersten RD, Yang YL, Moore BS, Dorrestein PC., Imaging Mass Spectrometry and Genome Mining via Short Sequence Tagging Identified the Anti-Infective Agent Arylomycin in *Streptomyces roseosporus*. *J Am Chem Soc*. 2011 Nov 16;133(45):18010-3.
- ⁸Gonzalez DJ, Okumura CY, Hollands A, Kersten R, Akong-Moore K, Pence MA, Malone CL, Derieux J, Moore BS, Horswill AR, Dixon JE, Dorrestein PC, Nizet V., Novel phenol-soluble modulin derivatives in community-associated methicillin-resistant *Staphylococcus aureus* identified through imaging mass spectrometry. *J Biol Chem*. 2012 Apr 20;287(17):13889-98.
- ⁹Xu Y, Kersten RD, Nam SJ, Lu L, Al-Suwailem AM, Zheng H, Fenical W, Dorrestein PC, Moore BS, Qian PY., Bacterial biosynthesis and maturation of the didemnin anti-cancer agents, *J Am Chem Soc*. 2012 May 23;134(20):8625-32.
- ¹⁰Gonzalez DJ, Xu Y, Yang YL, Esquenazi E, Liu WT, Edlund A, Duong T, Du L, Molnár I, Gerwick WH, Jensen PR, Fischbach M, Liaw CC, Straight P, Nizet V, Dorrestein PC., Observing the invisible through imaging mass spectrometry, a window into the metabolic exchange patterns of microbes, *J Proteomics*. 2012 Aug 30;75(16):5069-76.
- ¹¹Yang JY, Phelan VV, Simkovsky R, Watrous JD, Trial RM, Fleming TC, Wenter R, Moore BS, Golden SS, Pogliano K, Dorrestein PC., Primer on agar-based microbial imaging mass spectrometry, *J. Bacteriol*. 2012 Nov;194(22):6023-8.
- ¹²Moree WJ, Phelan VV, Wu CH, Bandeira N, Cornett DS, Duggan BM, Dorrestein PC., Interkingdom metabolic transformations captured by microbial imaging mass spectrometry. *Proc Natl Acad Sci U S A*. 2012 Aug 21;109(34):13811-6.
- ¹³Rath CM, Alexandrov T, Higginbottom SK, Song J, Milla ME, Fischbach MA, Sonnenburg JL, Dorrestein PC., Molecular analysis of model gut microbiotas by imaging mass spectrometry and nanodesorption electrospray ionization reveals dietary metabolite transformations. *Anal Chem*. 2012 Nov 6;84(21):9259-67.
- ¹⁴Watrous JD, Phelan VV, Hsu CC, Moree WJ, Duggan BM, Alexandrov T, Dorrestein PC., Microbial metabolic exchange in 3D, *ISME J*. 2013 Apr;7(4):770-80.
- ¹⁵Hoefler BC, Gorzelnik KV, Yang JY, Hendricks N, Dorrestein PC, Straight PD., Enzymatic resistance to the lipopeptide surfactin as identified through imaging mass spectrometry of bacterial competition. *Proc Natl Acad Sci U S A*. 2012 Aug 7;109(32):13082-7.

For research use only. Not for use in diagnostic procedures.

● Bruker Daltonik GmbH

Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

www.bruker.com/ms

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

Fremont, CA · USA
Phone +1 (510) 683-4300
Fax +1 (510) 687-1217
cam-sales@bdal.com