




America

CERTIFICATE

No. U8 094770 0004 Rev. 01

Holder of Certificate: **Tornado Spectral Systems Inc.**
5155 Spectrum Way, Unit 6
Mississauga ON L4W 5A1
CANADA

Certification Mark:



Product: **Laboratory equipment**
Process Guardian Raman

This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL for USA and a Standards Council of Canada ISO/IEC 17065 accredited Certification body for Canada.

Test report no.: 7169010267-001

Date, 2022-04-13

(Peter Keith)



America

CERTIFICATE

No. U8 094770 0004 Rev. 01

Model(s):

300100

Tested

CAN/CSA-C22.2 No. 61010-1:2012/R:2019-07

according to:

UL 61010-1:2012/R:2019-07

Parameters:

Rated Input Voltage: 120/240 V~

Rated Frequency: 50/60 Hz

Rated Input Current: 1.1 A

Protection Class: I

Laser Class: 3B