



# Laser-induced post-ionization for the enhanced MALDI-2-MS analysis of N-glycans

Bram Heijs<sup>1,2</sup>; Alexander Potthoff<sup>1</sup>; Hans Dalebout<sup>2</sup>; Jens Soltwisch<sup>1,3</sup>; Klaus Dreisewerd<sup>1,3</sup>

<sup>1</sup> Institute for Hygiene, University of Münster, Münster, Germany. <sup>2</sup> Center for Proteomics & Metabolomics, Leiden University Medical Center, Leiden, The Netherlands.

<sup>3</sup> Interdisciplinary Center for Clinical Research (IZKF), University of Münster, Münster, Germany.



b.p.a.m.heijs@lumc.nl

AIM: Assess the feasibility and optimize conditions for the analysis of N-glycans by MALDI-2 mass spectrometry imaging


## Introduction

The analysis of N-linked glycosylation has garnered significant interest in the biomedical and clinical research communities. N-glycans, a class of common post-translational protein modifications, are involved in numerous cellular processes, including cell-cell interactions, and signaling. Aberrant glycosylation patterns have been associated with disease, including autoimmune diseases, and multiple cancer types. One of the most recent analysis platforms for N-glycans is matrix-assisted laser desorption/ionization mass spectrometry imaging (MALDI-MSI). Laser-induced post-ionization coupled to MALDI-MS, MALDI-2-MS, has recently been shown to drastically enhance ion yields and measurement sensitivity for the MALDI-MS-based analysis of (glyco-)lipids [1]. Here we aim to investigate the possible benefits of MALDI-2-MS for the analysis of polysaccharides, to ultimately improve the MALDI-MSI analysis of N-glycans.

[1] Soltwisch *et al.* (2015) *Science*


## Experimentals

### MALDI-2 parameter optimization



Optimization of laser pulse energies, inter-laser delay, and cooling gas pressure

### N-glycan MALDI-2-MSI



On-tissue digestions with PNGase F

Overnight at 37 °C in a saturated humid environment

### Bruker Daltonics MALDI-2 timsTOF flex

QTOF-MS with trapped ion mobility separation module, modified for MALDI-2

Ablation laser: 10 kHz Bruker smartbeam 3D (355 nm)

PI laser: 1 kHz Ekspla N204-1k-FH (266 nm, 7 ns pulse length)

Pressure: 200 - 350 Pa



## Results & Discussion

### MALDI-2 parameter optimization on maltoheptaose standards




Figure 1: (A) In positive ion-mode sodium adducts did not benefit from MALDI-2. (B) In negative ion-mode deprotonated ions showed a beneficial intensity boost by MALDI-2. Optimization of the inter-laser delay (ILD) at low (C) and high (D) ablation laser intensity (ALI) revealed the presence of two ion populations. One fast moving population of single ions formed through thermal desorption (low ALI: < 100 Pa, ILD 10-20 μs) and one slower moving population of ablation clusters (low ALI: > 1 mbar, ILD 20-40 μs; high ALI: > 150 Pa, ILD 20-40 μs). Post-ionizing the ablation clusters resulted in the highest ion yields.

### MALDI-2 measurement sensitivity assessment




Figure 2: A dilution series of maltoheptaose was sprayed on glass slides and analyzed by (A) MALDI-2-MS in negative ion-mode, (B) MALDI-MS in negative ion-mode and (C) MALDI-MS positive ion-mode. The lower limit of detection (LLOD: S/N ≥ 3) for MALDI-2 (A) was three orders of magnitude lower in MALDI-2 compared to MALDI in negative ion-mode (B), and one order of magnitude compared to the current positive ion-mode gold standard (C).




Figure 3: The smartbeam 3D “M5 small” ablation pattern. The ablation area is 1280 μm<sup>2</sup>. Amounts of maltoheptaose per pixel in Figure 2 are based on this ablation pattern.

### Negative ion-mode MALDI-2-MSI and in-situ MALDI-2-MS/MS of N-glycans in formalin-fixed, paraffin-embedded human brain (cerebellum)




Figure 4: Comparison of positive and negative ion-mode N-glycan MSI average spectra (A) and ion distributions normalized to internal standard (B&C). (D&E) On-tissue MALDI-2-MS/MS spectra of N-glycans. Negative ion-mode provides many fragments that provide valuable structure information. More so than positive ion-mode.

## Conclusions

- MALDI-2-MS of deprotonated oligosaccharides provides substantial boosts in ion yields
- Optimal MALDI-2 conditions are:
  - Cooling gas pressure: 250 – 300 Pa
  - Inter-laser delay: 30 μs
- MALDI-2 provides high quality MALDI-MSI and on-tissue MS/MS data for N-glycans directly from tissue

## Acknowledgements

The authors would like to thank Bruker Daltonics for funding and assistance during the project