

Identification of Synthetic Heparan Sulfate Glycosaminoglycan Hexasaccharides Epimers by Capillary Zone Electrophoresis Negative Electron Transfer Dissociation Tandem Mass Spectrometry

UNIVERSITY OF
GEORGIA

Marshall Liss¹; Lauren Pepi¹; Pradeep Chopra²; Geert-Jan Boons²; I. Jonathan Amster¹

¹University of Georgia, Chemistry Department, Athens, GA; ²University of Georgia, Complex Carbohydrate Research Center, Athens, GA

Introduction

- Electron detachment dissociation (EDD) and negative electron transfer dissociation (NETD) have previously been demonstrated for analysis of highly sulfated glycosaminoglycans (GAGs)
- NETD rivals EDD fragmentation patterns but can easily be paired with separation techniques.
- Past work has addressed the ability of EDD to distinguish synthetic heparan sulfate (HS) hexamers differing in glucuronic acid (GlcA) vs iduronic acid (IdoA).
- MS/MS of HS hexamers in EDD and NETD produced abundant glycosidic products and cross-ring fragmentation.




Figure 1. Backbone structure for synthetic HS epimers. Red highlights regions which sulfation and stereochemistry vary. Eighteen unique HS hexamers were synthesized using a modular approach and investigated for this work.

Experimental

- NETD and EDD experiments were performed on a 12T Bruker solarix XR FT-ICR MS. (Figure 2).
- Each sample was injected at a rate of 10 μ L/h in 50:50 MeOH: H₂O and ionized by nanoESI at a concentration of 0.2mg/mL.
- Capillary-Zone Electrophoresis (CZE) -MS was performed on an Agilent G1600AX 3D CE coupled to a Thermo Scientific Velos Orbitrap Elite by an EMAS-II interface.




Figure 2. Diagram of the 12T Bruker solarix XR FT-ICR MS

Results and Discussion




Figure 4. GlcA-GlcNS6S-IdoA-GlcNS6S-GlcA/IdoA-GlcNS6S-O(CH₂)₅NH₂ (A) scores plot and (B) loadings plot from a PCA analysis. Left cluster contains GlcA at reducing end, right cluster contains IdoA at reducing end. B₅ ions were found to be diagnostic for GlcA, Y₄, Z₅, Y₁ and ^{2,4}X₀ ions to be diagnostic for IdoA




Figure 5. (A) GlcA sample fragment map using EDD fragmentation of precursor [M-6H+Na]⁵⁻. B₅ fragment ion determined to be diagnostic for GlcA as the reducing end uronic. (B) IdoA sample fragment map using EDD fragmentation of precursor [M-6H+Na]⁵⁻. The Z₅ fragment ion is determined by PCA to be diagnostic for IdoA.

Capillary-Zone Electrophoresis




Figure 3. Electropherogram of a HS hexamer. Peak at 17.64 depicts the expected structure, with the 17.94 peak depicting one less sulfate from the expected structure.

- CZE separates analytes based on their differences in effective electrophoretic mobilities. Since glycosaminoglycans can be highly charged substances, CZE becomes an effective separation technique.
- HS Hexamer samples can contain other unwanted GAGs or organic compounds.
- CZE sample introduction requires minimal sample for tandem MS, (~1 μ L).




Figure 6. NETD of GlcA-GlcNS6S-GlcA-GlcNS6S-IdoA-GlcNS6S. Fragments found in the PCA of EDD spectra were also produced in this NETD spectrum such as Y₄, and B₅ ions.

EDD and NETD Analysis

- PCA analysis of EDD spectra shows the effectiveness of electron-based ion activation methods. NETD techniques produce similar fragmentation patterns with respect to cross-ring and glycosidic bond cleavages.
- Diversity of ions (shown as in Figure 4) are diagnostic for GlcA/IdoA.
- Due to scarce sample quantities and impure samples, further analysis is not effective without CZE separations.
- Radial repulsion between a negatively charged ion cloud and electron beam make EDD ineffective on short time scales

Conclusions

- EDD displays restricted capabilities when coupled with electrophoretic separation techniques such as CZE.
- NETD is an effective technique when used with shorter analyte time scales.
- CZE elucidates a method to purify and conserve samples.
- Future work will pair CZE sample introduction with NETD methods in FT-ICR MS utilizing PCA analysis.

References

- Leach, F. E., Wolff, J. J., Xiao, Z., Ly, M., Laremore, T. N., Arungundram, S., ... Amster, I. J. (2011). Negative Electron Transfer Dissociation Fourier Transform Mass Spectrometry of Glycosaminoglycan Carbohydrates. *European Journal of Mass Spectrometry*, 17(2), 167–176. doi: 10.1255/ejms.1120
- Smith, R. D., Olivares, J. A., Nguyen, N. T., & Udseth, H. R. (1988). Capillary zone electrophoresis-mass spectrometry using an electrospray ionization interface. *Analytical Chemistry*, 60(5), 436–441. doi: 10.1021/ac00156a013

Acknowledgements

This project would not be possible without the help from my supervisor and professor Dr. Jonathan Amster, and my graduate student mentor Lauren Pepi. We gratefully acknowledge the financial support of the National Institutes of Health, grant U01 CA231074