



**PERFORMANCE OF NEW 400 MHZ HTS  
POWER-DRIVEN MAGNET NMR  
TECHNOLOGY ON CINACALCET HCL  
(TYPICAL API IN PHARMA)**

**MARIA VICTORIA SILVA ELIPE\*, NEIL DONOVAN, ROBERT  
KRULL, DONALD POOKE, KIM COLSON**

PANIC NMR 2018

**AMGEN**<sup>®</sup>


Pioneering science delivers vital medicines™

# OUTLINE


- **Historical highlights of magnet development**
- **Symbiotic scientific collaboration**
- **Scheme of the instrument**
- **Cinacalcet HCl**
- **NMR data (1D 1H, 2D COSY, 2D NOESY, 2D 1H-13C HSQC)**
- **HTS NMR system installed at Amgen chemistry laboratory**
- **Pros and cons of HTS magnet system**
- **Acknowledgements**
- **Thank you**

# HISTORICAL HIGHLIGHTS OF MAGNET DEVELOPMENT



Increase on MHz over the years



Electromagnets



# SYMBIOTIC SCIENTIFIC COLLABORATION




Maria is happy!



Provided March 6, 2018, as part of an oral presentation and is qualified by such, contains forward-looking statements, actual results may vary materially; Amgen disclaims any duty to update.

# SCHEMATICS OF THE INSTRUMENT




Silva Elipe, MV, Donovan, N, Krull, R, Pooke, D, Colson, K, 2018, submitted

Provided March 6, 2018, as part of an oral presentation and is qualified by such, contains forward-looking statements, actual results may vary materially; Amgen disclaims any duty to update.



# CINACALCET HCL

- API of commercial oral formulated tablets of Sensipar®/Mimpara®
- Approved in 2004 by FDA as Sensipar® and by EMEA as Mimpara®
- Therapeutic indications:
  - Secondary hyperparathyroidism (sHPT) for patients with chronic kidney disease (CKD) on dialysis
  - Primary hyperparathyroidism and parathyroid carcinoma



cinacalcet HCl



# COMPARING 1D 1H NMR ALIPHATIC REGION OF HTS NMR SYSTEM WITH CONVENTIONAL NMR INSTRUMENT



Silva Elipe, MV, Donovan, N, Krull, R, Pooke, D, Colson, K, 2018, submitted

Provided March 6, 2018, as part of an oral presentation and is qualified by such, contains forward-looking statements, actual results may vary materially; Amgen disclaims any duty to update.

# COMPARING 1D 1H NMR AROMATIC REGION OF HTS NMR SYSTEM WITH CONVENTIONAL NMR INSTRUMENT




Silva Elipe, MV, Donovan, N, Krull, R, Pooke, D, Colson, K, 2018, submitted

Provided March 6, 2018, as part of an oral presentation and is qualified by such, contains forward-looking statements, actual results may vary materially; Amgen disclaims any duty to update.




# COMPARING 2D COSY OF HTS NMR SYSTEM WITH CONVENTIONAL NMR INSTRUMENT



Silva Elipe, MV, Donovan, N, Krull, R, Pooke, D, Colson, K, 2018, submitted

Provided March 6, 2018, as part of an oral presentation and is qualified by such, contains forward-looking statements, actual results may vary materially; Amgen disclaims any duty to update.


# COMPARING 2D NOESY OF HTS NMR SYSTEM WITH CONVENTIONAL NMR INSTRUMENT



Silva Elipe, MV, Donovan, N, Krull, R, Pooke, D, Colson, K, 2018, submitted

Provided March 6, 2018, as part of an oral presentation and is qualified by such, contains forward-looking statements, actual results may vary materially; Amgen disclaims any duty to update.


# COMPARING 2D 1H-13C HSQC OF HTS NMR SYSTEM WITH CONVENTIONAL NMR INSTRUMENT



Silva Elipe, MV, Donovan, N, Krull, R, Pooke, D, Colson, K, 2018, submitted

Provided March 6, 2018, as part of an oral presentation and is qualified by such, contains forward-looking statements, actual results may vary materially; Amgen disclaims any duty to update.



# COMPARING 1D 1H NMR OF HTS NMR SYSTEM (BBFO/BOSS3 AND QNP/BOSS1) WITH CONVENTIONAL NMR INSTRUMENT



Silva Elipe, MV, Donovan, N, Krull, R, Pooke, D, Colson, K, 2018, submitted

Provided March 6, 2018, as part of an oral presentation and is qualified by such, contains forward-looking statements, actual results may vary materially; Amgen disclaims any duty to update.

# COMPARING 1D 1H NMR ALIPHATIC REGION HTS NMR SYSTEM (QNP/BOSS1) WITH CONVENTIONAL NMR INSTRUMENT



Silva Elipe, MV, Donovan, N, Krull, R, Pooke, D, Colson, K, 2018, submitted

Provided March 6, 2018, as part of an oral presentation and is qualified by such, contains forward-looking statements, actual results may vary materially; Amgen disclaims any duty to update.



# HTS NMR SYSTEM INSTALLED AT AMGEN CHEMISTRY LAB



Provided March 6, 2018, as part of an oral presentation and is qualified by such, contains forward-looking statements, actual results may vary materially; Amgen disclaims any duty to update.

# PROS AND CONS OF HTS MAGNET SYSTEM

## Pros:

- **No need to fill with liquid cryogenics (no liquid N2 and He needed) and no O2 sensors**
- **Easy access to chemist (in the same lab with no restricted access)**
- **Small footprint magnet (can site in hood, easy sample insertion, small stray field)**
- **Quick installation (couple of days to field and acquiring data)**
- **Ability to shutdown magnet during times of non-use (easy to do by the user, no need for engineer, magnet can be easily moved to other location if needed, no need for extra money)**

## Cons:

- **Sensitive to power outages (UPS/generator needed)**
- **Compressor & power supply must be nearby (real-estate)**
- **Water cooling or air-cooling needed**
- **Noisy equipment (He compressor)**

# ACKNOWLEDGEMENTS

- **Amgen Sponsors: Janet Cheetham & David Semin**
- **Cinacalcet HCl sample from Amgen: Minhui Ma & Robert Munger**
- **Bruker BioSpin: Clemens Anklin**
- **Others: many people helped through the process**



# THANK YOU



**Elvis is alive and has been working with a cryogenic magnet for decades**

But now ...



No cryogens????



Oops!!!