


Upgrade Your MultiMode!

You'll Be Surprised at the New Things Your Old AFM Can Do



● **Recent Bruker innovations now available with the MultiMode 8 include:**

ScanAsyst® Imaging Mode

- Makes operation easier for both expert and novice users with automatic scan optimization
- Delivers tremendously easier operation in liquid, with no need to “tune” the cantilever and no setpoint drift

PeakForce QNM® Mode

- Enables direct quantitative mapping of nanoscale sample modulus and adhesion
- Provides unambiguous and non-destructive property mapping on a wide range of sample types

Fast Scanning for the MultiMode®

- Offers 20X faster survey scans and 6X faster scans with no loss of resolution in air with ScanAsyst-HR

Quantitative Nanoscale Characterization

- Renders the highest resolution quantitative surface potential imaging with PeakForce KPFM™
- Provides high-resolution, non-destructive, quantitative conductivity mapping with PeakForce TUNA™

Advanced EC with Turnkey Glovebox Solution

- Supplies sub-ppm water and oxygen control for organic photovoltaics, battery research, and other applications

New Service Contract Programs

- Replace older NanoScope® controllers at or near the end of service life
- Keep your MultiMode running long into the future

Every MultiMode can be upgraded.

For more information on improving the performance of your MultiMode, please contact your Bruker representative at productinfo@bruker-nano.com