


Valentin Göldner<sup>1,2,3</sup>, Linda Quach<sup>2,4</sup>, Egy Adhitama<sup>2,5</sup>, Arne Behrens<sup>1,6</sup>, Luisa Junk<sup>1</sup>, Martin Winter<sup>2,5,7</sup>, Tobias Placke<sup>2,5</sup>, Frank Glorius<sup>2,4</sup>, and Uwe Karst<sup>1,2</sup>

<sup>1</sup>Institute of Inorganic and Analytical Chemistry, University of Münster, <sup>2</sup>International Graduate School for Battery Chemistry, Characterization, Analysis, Recycling and Application, University of Münster,

<sup>3</sup>Department of Environmental Geosciences, Centre for Microbiology and Environmental Systems Science, University of Vienna, <sup>4</sup>Institute of Organic Chemistry, University of Münster;


<sup>5</sup>MEET Battery Research Center, University of Münster, <sup>6</sup>Bruker Daltonik GmbH; <sup>7</sup>Helmholtz Institute Münster

## Lithium ion batteries (LIBs)




- LIB electrolyte is degraded during charging
- Reduction products form the solid electrolyte interphase (SEI) on the negative electrode
- SEI prevents further electrolyte degradation
- A stable SEI is important for battery lifetime

## Silicon as active material for negative electrodes




- Silicon is a promising material for negative electrodes
- Specific capacity is nearly ten times higher than with graphite<sup>1</sup>
- Severe volume change upon (de-)lithiation (up to 300%)<sup>1</sup>
- Continuous (re-)formation of SEI and consumption of active Li<sup>1</sup>
- More flexible SEI is required to withstand mechanic stress

## Film-forming electrolyte additives




- Electrolyte additives are used to tailor SEI properties
- N*-Carboxyanhydrides (*N*-CAs) enhance battery lifetime<sup>2</sup>
- Polymerization after electrochemical reduction proposed<sup>2</sup>
- Polymers incorporated in the SEI may increase flexibility<sup>2</sup>
- No suitable detection method for these polymers available

- A timsTOF fleX mass spectrometer from Bruker Daltonics (Bremen, Germany) was used for all analyses
- Laser desorption/ionization was used in combination with trapped ion mobility spectrometry and time-of-flight mass spectrometry in positive ion mode (LDI(+)-TIMS-TOF-MS)
- The samples were introduced into the mass spectrometer under Ar-atmosphere to avoid contact to ambient air




## LDI-MS analysis of electrodes cycled with different *N*-CAs:




- Detection of oligomers derived from *N*-CAs after cycling
- Substituents influence the size of the formed oligomers
- LDI-MS can support a more targeted additive design

## Fast structural elucidation by simultaneous fragmentation of oligomers using prm-PASEF:



## Mass spectrometry imaging (MSI) visualizes lateral distribution on the electrode surface:



- Spot size: 50  $\mu\text{m}$
- LDI-MSI visualizes the distribution of SEI components on electrode surfaces
- Substituents of *N*-CA additives influence the lateral distribution patterns of the resulting oligomers on electrodes
- PhGly-*N*-CA leads to more homogenous SEI formation than Ala-*N*-CA
- Homogeneity of the SEI layer influences battery performance and lifetime

- LDI-MS is suitable for the detection and characterization of SEI components in LIBs derived from electrolyte additives
- Obtained information is complementary to established methods like XPS, SIMS, FTIR and Raman Spectroscopy
- Upon reduction, *N*-CA-based electrolyte additives form oligomers which increase the flexibility of the SEI
- TIMS separation and prm-PASEF fragmentation enhance structural elucidation and speed up the data acquisition
- LDI-MSI analysis reveals the lateral distribution of SEI components on electrode surfaces
- The acquired data supports a more targeted approach on the design of electrolyte additives

**References:** <sup>1</sup>A. Franco Gonzalez, N. H. Yang, R. S. Liu, *J. Phys. Chem. C*, **2017**, *121*, 27775–27787. <sup>2</sup>J.-P. Schmiegel, R. Nölle, J. Henschel, L. Quach, S. Nowak, M. Winter, F. Glorius, T. Placke, *Cell Reports Phys. Sci.* **2021**, *2*, 100327.

**Acknowledgement:** This work was supported by the International Graduate School for Battery Chemistry, Characterization, Analysis, Recycling and Application (BACCARA), which is funded by the Ministry for Culture and Science of North Rhine-Westphalia, Germany.

**Publications:** V. Göldner, L. Quach, E. Adhitama, A. Behrens, L. Junk, M. Winter, T. Placke, F. Glorius, U. Karst, *iScience* **2023**, *26*, 107517. L. Quach, E. Adhitama, V. Göldner, A. Das, F. Demelash, M. Winter, U. Karst, T. Placke, F. Glorius, *ACS Appl. Energy Mater.* **2023**, *6*, 9837–9850.