



MALDI Biotyper®

MBT Compass HT CA Software

Empowering unequaled speed and optimal performance

Maximize the power of your MALDI Biotyper® CA System with MBT Compass HT CA Software

What's new?

- Accelerate your system's performance for unmatched speed
- Prioritize urgent samples with streamlined processing
- Ensure consistent results with the IDealTune™ auto-optimization feature
- Simplify compliance with built-in 21 CFR Part 11 support
- Stay protected with proactive cybersecurity patching
- Enhance quality control and inventory workflows with greater efficiency
- Unlock deeper epidemiological insights with advanced analytics
- Connect multiple workstations to optimize your entire workflow

Boost system speed for maximum throughput

Achieve analysis of up to 600 samples per hour - or just 5 minutes for a fully loaded MBT Biotarget 96 US IVD plate, including target exchange and report generation.

Prioritize urgent samples with confidence

Use the sample prioritization feature to deliver critical results in under 3 minutes. The interrupted run can be resumed seamlessly - no need to prepare samples again.

Maintain peak performance with IDEalTune™

This intelligent auto-tuning feature continuously monitors data quality and adjusts acquisition parameters automatically - no hands-on time or extra sample preparation required.

Simplify compliance with 21 CFR Part 11 support

Take advantage of enhanced user management, comprehensive audit trails, and integrated electronic signatures to meet regulatory requirements.

Stay protected with proactive cybersecurity updates

Safeguard your system and data integrity, benefit from regular threat scans and timely patch releases.

Streamline QC and inventory management

Improve traceability and efficiency with software-assisted reagent lot tracking and inventory control.

Gain actionable epidemiological insights

Leverage built-in organism tracking and trending tools to support local and regional surveillance efforts.

What is MBT Compass HT CA software?

MBT Compass HT CA is the next-generation software for the MALDI Biotyper® CA System, offering a modern graphical user interface (GUI) with built-in compliance functionality and significantly enhanced processing performance.

The software continues to support the two established workflows - Standard Identification of bacteria and yeast, and the valuable Sepsityper® workflow - with no changes to their functionality, except for a welcome shorter time-to-result.

Its modernized architecture enables the seamless addition of new features - without affecting existing workflows. This ensures long-term flexibility and scalability for your laboratory needs.

Only for sale in the USA and Puerto Rico.

IVD

Rx ONLY

MALDI Biotyper® is a registered trademark of the Bruker group of companies.

Online information
bruker.com/microbiology



Bruker Daltonics GmbH and Co. KG

Bremen · Germany
Phone +49 (0) 421-2205-0

info.md@bruker.com

Bruker Scientific LLC

Billerica, MA · USA
Phone +1 (978) 663-3660

