


High Power Filters

- RF Specifications
Technical Data
Version 09



Copyright © by Bruker Corporation

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form, or by any means without the prior consent of the publisher. Product names used are trademarks or registered trademarks of their respective holders.

© December 08, 2020 Bruker Corporation

Document Number: 1656205067

P/N: Z31481

Contents

1	NMR Filters	7
1.1	Quick Search - Typical Double Resonance Experiments	8
1.2	Typical Triple Resonance Experiments	12
1.2.1	Quick Search - Triple Resonance Experiments	12
1.2.1.1	300MHz Typical Experiments	14
1.2.1.2	400MHz Typical Experiments	15
1.2.1.3	500MHz Typical Experiments	16
1.2.1.4	600MHz Typical Experiments	17
1.2.1.5	700MHz Typical Experiments	18
1.2.1.6	750MHz Typical Experiments	19
1.2.1.7	800MHz Typical Experiments	20
1.2.1.8	850MHz Typical Experiments	21
1.2.1.9	900MHz Typical Experiments	22
1.2.1.10	950MHz Typical Experiments	23
1.2.1.11	1000MHz Typical Experiments	24
1.2.1.12	1100MHz Typical Experiments	25
1.2.1.13	1200MHz Typical Experiments	26
1.2.2	Triple Resonance 1H/X/Y	27
1.2.3	Triple resonance 1H/19F/X	28
1.3	AVANCE 300MHz	30
1.3.1	Lowpass Filter	30
1.3.2	Bandstop Filter	32
1.3.3	Bandpass Filter	34
1.3.4	Diplexer Filter	42
1.4	AVANCE 400MHz	43
1.4.1	Lowpass Filter	43
1.4.2	Bandstop Filter	45
1.4.3	Bandpass Filter	47
1.4.4	Diplexer Filter	56
1.5	AVANCE 500MHz	57
1.5.1	Lowpass Filter	57
1.5.2	Bandstop Filter	59
1.5.3	Bandpass Filter	64
1.5.4	Diplexer Filter	75
1.6	AVANCE 600MHz	77
1.6.1	Lowpass Filter	77
1.6.2	Bandstop Filter	79
1.6.3	Bandpass Filter	82
1.6.4	Diplexer Filter	92
1.7	AVANCE 700MHz	93
1.7.1	Lowpass Filter	93
1.7.2	Bandstop Filter	95
1.7.3	Bandpass Filter	97

1.7.4	Diplexer Filter	106
1.8	AVANCE 750MHz	107
1.8.1	Lowpass Filter	107
1.8.2	Bandstop Filter	109
1.8.3	Bandpass Filter	111
1.8.4	Diplexer Filter	118
1.9	AVANCE 800MHz	119
1.9.1	Lowpass Filter	119
1.9.2	Bandstop Filter	121
1.9.3	Bandpass Filter	123
1.10	AVANCE 850MHz	129
1.10.1	Lowpass Filter	129
1.10.2	Bandpass Filter	133
1.10.3	Diplexer Filter	140
1.11	AVANCE 900MHz	141
1.11.1	Lowpass Filter	141
1.11.2	Bandstop Filter	143
1.11.3	Bandpass Filter	144
1.11.4	Diplexer Filter	151
1.12	AVANCE 950MHz	152
1.12.1	Lowpass Filter	152
1.12.2	Bandstop Filter	154
1.12.3	Bandpass Filter	155
1.12.4	Diplexer Filter	162
1.13	AVANCE 1000MHz	163
1.13.1	Lowpass Filter	163
1.13.2	Bandstop Filter	165
1.13.3	Bandpass Filter	167
1.13.4	Diplexer Filter	173
1.14	AVANCE 1100MHz	174
1.14.1	Lowpass Filter	174
1.14.2	Bandpass Filter	176
1.14.3	Diplexer Filter	182
1.15	AVANCE 1200MHz	183
1.15.1	Lowpass Filter	183
1.15.2	Bandpass Filter	185
1.15.3	Diplexer Filter	192
1.16	Filters Category	193
1.16.1	Lowpass Filter	194
1.16.2	Bandstop Filter	196
1.16.3	Bandpass Filters	200
1.16.4	Diplexer Filter	214
2	Cases and Accessories	215
2.1	Notch Filters	215
2.1.1	Case 3S063	215
2.1.2	Case 4S063D	215
2.2	Bandpass Filters	216

2.2.1	Case 5R032	216
2.2.2	Case 5R045	216
2.2.3	Case 5R069	217
2.2.4	Case 5R100	217
2.2.5	Case 6R125	218
2.2.6	Case 7R155	218
2.2.7	Case 8R055P	219
2.2.8	Case 8R105P	219
2.2.9	Case 8R041P	220
2.3	Helical Filters	221
2.3.1	Case 4H073	221
2.3.2	Case 5H133	221
2.3.3	Case 5H069	222
2.3.4	Case 6H073	222
2.3.5	Case 8H073	223
2.3.6	Case 8H097	223
2.3.7	Case 8H139	224
2.3.8	Case 8H169	225
2.4	Low Pass Filters	226
2.4.1	Case 2D050	226
2.4.2	Case 2D095	226
2.5	Accessories	227
3	Typical Curves	229
3.1	Bandpass Filters Curves	229
4	Chemicool	237
4.1	Periodic Table	237
4.2	NMR Applications Frequencies	238
4.2.1	NMR Frequency Table	238
4.2.1.1	NMR Frequency Table (1)	239
4.2.1.2	NMR Frequency Table (2)	240
4.2.2	With increasing NMR Frequency	241
4.2.2.1	With increasing NMR Frequency (1)	242
4.2.2.2	With increasing NMR Frequency (2)	243
5	Contact	245
	List of Figures	247
	List of Tables	251

1 NMR Filters

1.1 Quick Search - Typical Double Resonance Experiments

Pass \ Stop	¹ H	¹⁹ F	³¹ P	¹³ C	² H	¹⁵ N
¹ H	*	W1346235	W1346239	W1346239	W1346239	W1346239
¹⁹ F	W1346236	*	W1346249	W1346249	W1346249	W1346249
³¹ P	W1346642	W1346642	*	W1346372	W1346372	W1346372
¹³ C	W1346642	W1346642	W1346271	*	W1346271	W1346275
² H	W1346642	W1346642	W1346273	W1346273	*	W1346273
¹⁵ N	W1346642	W1346642	W1346278	W1346278	W135032	*

Figure 1.1: AVANCE 300MHz [[▶ 30](#)]

Pass \ Stop	¹ H	¹⁹ F	³¹ P	¹³ C	² H	¹⁵ N
¹ H	*	W1346204	W1346240	W1346240	W1346240	W1346240
¹⁹ F	W1346205	*	W1346250	W1346250	W1346250	W1346250
³¹ P	W1346636	W1346636	*	W1346349	W1346349	W1346349
¹³ C	W1346636	W1346636	W1346488	*	W1346488	W1346226
² H	W1346636	W1346636	W1346260	W1346260	*	W1346260
¹⁵ N	W1346636	W1346636	W1346261	W1346628	W1346261	*

Figure 1.2: AVANCE 400MHz [[▶ 43](#)]

Pass \ Stop	¹ H	¹⁹ F	³¹ P	¹³ C	² H	¹⁵ N
¹ H	*	W1346247	W1346223	W1346223	W1346223	W1346223
¹⁹ F	W1346248	*	W1346251	W1346251	W1346251	W1346251
³¹ P	W1346637	W1346637	*	W1346351	W1346351	W1346351
¹³ C	W1346637	W1346637	W1346489	*	W1346489	W1346224
² H	W1346637	W1346637	W1346266	W1346266	*	W1346266
¹⁵ N	W1346637	W1346637	W1346267	W1346629	W1346267	*

Figure 1.3: AVANCE 500MHz [[▶ 57](#)]

Pass \ Stop	¹ H	¹⁹ F	³¹ P	¹³ C	² H	¹⁵ N
¹ H	*	W1346228	W1346241	W1346241	W1346241	W1346241
¹⁹ F	W1346229	*	W1346252	W1346252	W1346252	W1346252
³¹ P	W1346638	W1346638	*	W1346367	W1346367	W1346367
¹³ C	W1346638	W1346638	W1346490	*	W1346490	W1346230
² H	W1346638	W1346638	W1346282	W1346282	*	W1346282
¹⁵ N	W1346638	W1346638	W1346277	W1346630	W1346277	*

Figure 1.4: AVANCE 600MHz [[▶ 77](#)]

Pass \ Stop	¹ H	¹⁹ F	³¹ P	¹³ C	² H	¹⁵ N
¹ H	*	W1346288	W1346296	W1346296	W1346296	W1346296
¹⁹ F	W1346292	*	W1346300	W1346300	W1346300	W1346300
³¹ P	W1346639	W1346639	*	W1346304	W1346304	W1346304
¹³ C	W1346639	W1346639	W1346491	*	W1346491	W1346316
² H	W1346639	W1346639	W1346324	W1346324	*	W1346324
¹⁵ N	W1346639	W1346639	W1346328	W1346659	W1346328	*

Figure 1.5: AVANCE 700MHz [▶ 93]

Pass \ Stop	¹ H	¹⁹ F	³¹ P	¹³ C	² H	¹⁵ N
¹ H	*	W1346289	W1346297	W1346297	W1346297	W1346297
¹⁹ F	W1346293	*	W1346301	W1346301	W1346301	W1346301
³¹ P	W1346644	W1346644	*	W1346305	W1346305	W1346305
¹³ C	W1346644	W1346644	W1346313	*	W1346313	W1346317
² H	W1346644	W1346644	W1346493	W1346493	*	W1346493
¹⁵ N	W1346644	W1346644	W1346329	W1346660	W1346329	*

Figure 1.6: AVANCE 750MHz [▶ 107]

Pass \ Stop	¹ H	¹⁹ F	³¹ P	¹³ C	² H	¹⁵ N
¹ H	*	W1346290	W1346298	W1346298	W1346298	W1346298
¹⁹ F	W1346294	*	W1346302	W1346302	W1346302	W1346302
³¹ P	W1346640	W1346640	*	W1346306	W1346306	W1346306
¹³ C	W1346640	W1346640	W1346314	*	W1346314	W1346318
² H	W1346640	W1346640	W1346326	W1346326	*	W1346326
¹⁵ N	W1346640	W1346640	W1346330	W1346661	W1346330	*

Figure 1.7: AVANCE 800MHz [▶ 119]

Pass \ Stop	¹ H	¹⁹ F	³¹ P	¹³ C	² H	¹⁵ N
¹ H	*	W1346674	W1346655	W1346655	W1346655	W1346655
¹⁹ F	W1346676	*	W1346675	W1346675	W1346675	W1346675
³¹ P	W1346677	W1346677	*	W1346715	W1346715	W1346715
¹³ C	W1346677	W1346677	W1346710	*	W1346710	W1346657
² H	W1346677	W1346677	W1346697	W1346697	*	W1346697
¹⁵ N	W1346677	W1346677	W1346716	W1346656	W1346716	*

Figure 1.8: AVANCE 850MHz [▶ 129]

NMR Filters

Pass \ Stop	¹ H	¹⁹ F	³¹ P	¹³ C	² H	¹⁵ N
¹ H	*	W1346291	W1346299	W1346299	W1346299	W1346299
¹⁹ F	W1346295	*	W1346303	W1346303	W1346303	W1346303
³¹ P	W1346641	W1346641	*	W1346307	W1346307	W1346307
¹³ C	W1346641	W1346641	W1346315	*	W1346315	W1346319
² H	W1346641	W1346641	W1346327	W1346327	*	W1346327
¹⁵ N	W1346641	W1346641	W1346331	W1346662	W1346331	*

Figure 1.9: AVANCE 900MHz [[▶ 141](#)]

Pass \ Stop	¹ H	¹⁹ F	³¹ P	¹³ C	² H	¹⁵ N
¹ H	*	W1346869	W1346707	W1346707	W1346707	W1346707
¹⁹ F	W143774	*	W143775	W143775	W143775	W143775
³¹ P	W143786	W143786	*	W143776	W143776	W143776
¹³ C	W143786	W143786	W1346706	*	W1346706	W1346706
² H	W143786	W143786	W139785	W139785	*	W139785
¹⁵ N	W143786	W143786	W1346705	W143787	W1346705	*

Figure 1.10: AVANCE 950MHz [[▶ 152](#)]

Pass \ Stop	¹ H	¹⁹ F	³¹ P	¹³ C	² H	¹⁵ N
¹ H	*	W1346870	W1346665	W1346665	W1346665	W1346665
¹⁹ F	W1346871	*	W1346666	W1346666	W1346666	W1346666
³¹ P	W1346678	W1346678	*	W1346667	W1346667	W1346667
¹³ C	W1346678	W1346678	W1346669	*	W1346669	W1346669
² H	W1346678	W1346678	W1346671	W1346671	*	W1346671
¹⁵ N	W1346678	W1346678	W1346717	W1346672	W1346717	*

Figure 1.11: AVANCE 1000MHz [[▶ 163](#)]

Pass \ Stop	¹ H	¹⁹ F	³¹ P	¹³ C	² H	¹⁵ N
¹ H	*	W1346872	W1342132	W1342132	W1342132	W1342132
¹⁹ F	W1346873	*	W1346473	W1346473	W1346473	W1346473
³¹ P	W1342373	W1342373	*	W1346859	W1346859	W1346859
¹³ C	W1342373	W1342373	W1348165	*	W1348165	W1348165
² H	W1342373	W1342373	W1346865	W1346865	*	W1346865
¹⁵ N	W1342373	W1342373	W1346867	W1342624	W1346867	*

Figure 1.12: AVANCE 1100MHz [[▶ 174](#)]














Pass \ Stop	¹ H	¹⁹ F	³¹ P	¹³ C	² H	¹⁵ N
¹ H	*	W1346874	W1348174	W1348174	W1348174	W1348174
¹⁹ F	W1346879	*	W135027	W135027	W135027	W135027
³¹ P	W134925	W134925	*	W1346860	W1346860	W1346860
¹³ C	W134925	W134925	W135028	*	W135028	W135028
² H	W134925	W134925	W1346866	W1346866	*	W1346866
¹⁵ N	W134925	W134925	W1346868	W134926	W1346868	*

Figure 1.13: AVANCE 1200MHz [▶ 183]

1.2 Typical Triple Resonance Experiments

1.2.1 Quick Search - Triple Resonance Experiments

See also

-  [300MHz Typical Experiments \[▶ 14\]](#)
-  [400MHz Typical Experiments \[▶ 15\]](#)
-  [500MHz Typical Experiments \[▶ 16\]](#)
-  [600MHz Typical Experiments \[▶ 17\]](#)
-  [700MHz Typical Experiments \[▶ 18\]](#)
-  [750MHz Typical Experiments \[▶ 19\]](#)
-  [800MHz Typical Experiments \[▶ 20\]](#)
-  [850MHz Typical Experiments \[▶ 21\]](#)
-  [900MHz Typical Experiments \[▶ 22\]](#)
-  [950MHz Typical Experiments \[▶ 23\]](#)
-  [1000MHz Typical Experiments \[▶ 24\]](#)
-  [1100MHz Typical Experiments \[▶ 25\]](#)
-  [1200MHz Typical Experiments \[▶ 26\]](#)

300MHz Triple Resonance Solid Experiments

n°1	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346272	⁵⁹ Co	¹³ C W1346271	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346239
n°2	¹⁵ N	¹⁷ O	² H W1346273	²⁹ Si	⁵⁹ Co	¹³ C W1346271	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346239
n°3	¹⁵ N W1346278	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1346275	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346239
n°4	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346272	⁵⁹ Co	¹³ C	²⁷ Al W1346271	²³ Na	¹¹ B	³¹ P	¹ H W1346239
n°5	¹⁵ N	¹⁷ O W1346273	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346270	³¹ P	¹ H W1346239
n°6	¹⁵ N W1346278	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346270	³¹ P	¹ H W1346239
n°7	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na W1346271	¹¹ B	³¹ P W1346372	¹ H W1346239
n°8	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al W1346271	²³ Na	¹¹ B	³¹ P W1346372	¹ H W1346239
n°9	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1346271	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346372	¹ H W1346239
n°10	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co W1346271	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346372	¹ H W1346239
n°11	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346272	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346372	¹ H W1346239
n°12	¹⁵ N	¹⁷ O	² H W1346273	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346372	¹ H W1346239
n°13	¹⁵ N	¹⁷ O W1346273	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346372	¹ H W1346239

400MHz Triple Resonance Solid Experiments

n°1	^{15}N	^{17}O	^2H	^{29}Si W1346259	^{59}Co	^{13}C W1346488	^{27}Al	^{23}Na	^{11}B	^{31}P	^1H W1346240
n°2	^{15}N	^{17}O	^2H W1346260	^{29}Si	^{59}Co	^{13}C W1346440	^{27}Al	^{23}Na	^{11}B	^{31}P	^1H W1346240
n°3	^{15}N W1346628	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C W1346226	^{27}Al	^{23}Na	^{11}B	^{31}P	^1H W1346240
n°4	^{15}N	^{17}O	^2H	^{29}Si W1346259	^{59}Co	^{13}C	^{27}Al W1346488	^{23}Na	^{11}B	^{31}P	^1H W1346240
n°5	^{15}N	^{17}O W1346260	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B W1346257	^{31}P	^1H W1346240
n°6	^{15}N W1346261	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B W1346257	^{31}P	^1H W1346240
n°7	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na W1346488	^{11}B	^{31}P W1346349	^1H W1346240
n°8	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al W1346488	^{23}Na	^{11}B	^{31}P W1346349	^1H W1346240
n°9	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C W1346488	^{27}Al	^{23}Na	^{11}B	^{31}P W1346349	^1H W1346240
n°10	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co W1346258	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346349	^1H W1346240
n°11	^{15}N	^{17}O	^2H	^{29}Si W1346259	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346349	^1H W1346240
n°12	^{15}N	^{17}O	^2H W1346260	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346349	^1H W1346240
n°13	^{15}N	^{17}O W1346260	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346349	^1H W1346240

500MHz Triple Resonance Solid Experiments

n°1	^{15}N	^{17}O	^2H	^{29}Si W1346265	^{59}Co	^{13}C W1346489	^{27}Al	^{23}Na	^{11}B	^{31}P	^1H W1346223
n°2	^{15}N	^{17}O	^2H W1346266	^{29}Si	^{59}Co	^{13}C W1346460	^{27}Al	^{23}Na	^{11}B	^{31}P	^1H W1346223
n°3	^{15}N W1346629	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C W1346224	^{27}Al	^{23}Na	^{11}B	^{31}P	^1H W1346223
n°4	^{15}N	^{17}O	^2H	^{29}Si W1346265	^{59}Co	^{13}C	^{27}Al W1346489	^{23}Na	^{11}B	^{31}P	^1H W1346223
n°5	^{15}N	^{17}O W1346266	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B W1346263	^{31}P	^1H W1346223
n°6	^{15}N W1346267	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B W1346263	^{31}P	^1H W1346223
n°7	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na W1346489	^{11}B	^{31}P W1346351	^1H W1346223
n°8	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al W1346489	^{23}Na	^{11}B	^{31}P W1346351	^1H W1346223
n°9	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C W1346489	^{27}Al	^{23}Na	^{11}B	^{31}P W1346351	^1H W1346223
n°10	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co W1346264	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346351	^1H W1346223
n°11	^{15}N	^{17}O	^2H	^{29}Si W1346265	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346351	^1H W1346223
n°12	^{15}N	^{17}O	^2H W1346266	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346351	^1H W1346223
n°13	^{15}N	^{17}O W1346266	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346351	^1H W1346223

600MHz Triple Resonance Solid Experiments

n°1	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346281	⁵⁹ Co	¹³ C W1346490	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346241
n°2	¹⁵ N	¹⁷ O	² H W1346282	²⁹ Si	⁵⁹ Co	¹³ C W1346461	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346241
n°3	¹⁵ N W1346630	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1346230	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346241
n°4	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346281	⁵⁹ Co	¹³ C	²⁷ Al W1346490	²³ Na	¹¹ B	³¹ P	¹ H W1346241
n°5	¹⁵ N	¹⁷ O W1346282	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346280	³¹ P	¹ H W1346241
n°6	¹⁵ N W1346277	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346280	³¹ P	¹ H W1346241
n°7	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na W1346490	¹¹ B	³¹ P W1346367	¹ H W1346241
n°8	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al W1346490	²³ Na	¹¹ B	³¹ P W1346367	¹ H W1346241
n°9	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1346490	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346367	¹ H W1346241
n°10	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co W1346276	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346367	¹ H W1346241
n°11	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346281	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346367	¹ H W1346241
n°12	¹⁵ N	¹⁷ O	² H W1346282	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346367	¹ H W1346241
n°13	¹⁵ N	¹⁷ O W1346282	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346367	¹ H W1346241

700MHz Triple Resonance Solid Experiments

n°1	^{15}N	^{17}O	^2H	^{29}Si W1346320	^{59}Co	^{13}C W1346491	^{27}Al	^{23}Na	^{11}B	^{31}P	^1H W1346296
n°2	^{15}N	^{17}O	^2H W1346492	^{29}Si	^{59}Co	^{13}C W1346462	^{27}Al	^{23}Na	^{11}B	^{31}P	^1H W1346296
n°3	^{15}N W1346659	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C W1346316	^{27}Al	^{23}Na	^{11}B	^{31}P	^1H W1346296
n°4	^{15}N	^{17}O	^2H	^{29}Si W1346320	^{59}Co	^{13}C	^{27}Al W1346491	^{23}Na	^{11}B	^{31}P	^1H W1346296
n°5	^{15}N	^{17}O W1346492	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B W1346308	^{31}P	^1H W1346296
n°6	^{15}N W1346328	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B W1346308	^{31}P	^1H W1346296
n°7	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na W1346491	^{11}B	^{31}P W1346304	^1H W1346296
n°8	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al W1346491	^{23}Na	^{11}B	^{31}P W1346304	^1H W1346296
n°9	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C W1346491	^{27}Al	^{23}Na	^{11}B	^{31}P W1346304	^1H W1346296
n°10	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co W1346312	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346304	^1H W1346296
n°11	^{15}N	^{17}O	^2H	^{29}Si W1346320	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346304	^1H W1346296
n°12	^{15}N	^{17}O	^2H W1346492	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346304	^1H W1346296
n°13	^{15}N	^{17}O W1346492	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346304	^1H W1346296

750MHz Triple Resonance Solid Experiments

n°1	^{15}N	^{17}O	^2H	^{29}Si W1346321	^{59}Co	^{13}C W1346313	^{27}Al	^{23}Na	^{11}B	^{31}P	^1H W1346297
n°2	^{15}N	^{17}O	^2H W1346493	^{29}Si	^{59}Co	^{13}C W1346313	^{27}Al	^{23}Na	^{11}B	^{31}P	^1H W1346297
n°3	^{15}N W1346660	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C W1346317	^{27}Al	^{23}Na	^{11}B	^{31}P	^1H W1346297
n°4	^{15}N	^{17}O	^2H	^{29}Si W1346321	^{59}Co	^{13}C	^{27}Al W1346313	^{23}Na	^{11}B	^{31}P	^1H W1346297
n°5	^{15}N	^{17}O W1346493	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B W1346309	^{31}P	^1H W1346297
n°6	^{15}N W1346329	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B W1346309	^{31}P	^1H W1346297
n°7	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na W1346313	^{11}B	^{31}P W1346305	^1H W1346297
n°8	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al W1346313	^{23}Na	^{11}B	^{31}P W1346305	^1H W1346297
n°9	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C W1346313	^{27}Al	^{23}Na	^{11}B	^{31}P W1346305	^1H W1346297
n°10	^{15}N	^{17}O	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346305	^1H W1346297
n°11	^{15}N	^{17}O	^2H	^{29}Si W1346321	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346305	^1H W1346297
n°12	^{15}N	^{17}O	^2H W1346493	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346305	^1H W1346297
n°13	^{15}N	^{17}O W1346493	^2H	^{29}Si	^{59}Co	^{13}C	^{27}Al	^{23}Na	^{11}B	^{31}P W1346305	^1H W1346297

800Hz Triple Resonance Solid Experiments

n°1	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346322	⁵⁹ Co	¹³ C W1346314	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346298
n°2	¹⁵ N	¹⁷ O	² H W1346326	²⁹ Si	⁵⁹ Co	¹³ C W1346314	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346298
n°3	¹⁵ N W1346661	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1346318	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346298
n°4	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346322	⁵⁹ Co	¹³ C	²⁷ Al W1346314	²³ Na	¹¹ B	³¹ P	¹ H W1346298
n°5	¹⁵ N	¹⁷ O W1346326	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346310	³¹ P	¹ H W1346298
n°6	¹⁵ N W1346330	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346310	³¹ P	¹ H W1346298
n°7	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na W1346314	¹¹ B	³¹ P W1346306	¹ H W1346298
n°8	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al W1346314	²³ Na	¹¹ B	³¹ P W1346306	¹ H W1346298
n°9	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1346314	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346306	¹ H W1346298
n°10	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346306	¹ H W1346298
n°11	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346322	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346306	¹ H W1346298
n°12	¹⁵ N	¹⁷ O	² H W1346326	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346306	¹ H W1346298
n°13	¹⁵ N	¹⁷ O W1346326	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346306	¹ H W1346298

850MHz Triple Resonance Solid Experiments

n°1	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346718	⁵⁹ Co	¹³ C W1346710	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346655
n°2	¹⁵ N	¹⁷ O	² H W1346697	²⁹ Si	⁵⁹ Co	¹³ C W1346710	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346655
n°3	¹⁵ N W1346656	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1346657	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346655
n°4	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346718	⁵⁹ Co	¹³ C	²⁷ Al W1346710	²³ Na	¹¹ B	³¹ P	¹ H W1346655
n°5	¹⁵ N	¹⁷ O W1346697	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346714	³¹ P	¹ H W1346655
n°6	¹⁵ N W1346716	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346714	³¹ P	¹ H W1346655
n°7	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na W1346710	¹¹ B	³¹ P W1346715	¹ H W1346655
n°8	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al W1346710	²³ Na	¹¹ B	³¹ P W1346715	¹ H W1346655
n°9	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1346710	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346715	¹ H W1346655
n°10	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346715	¹ H W1346655
n°11	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346718	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346715	¹ H W1346655
n°12	¹⁵ N	¹⁷ O	² H W1346697	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346715	¹ H W1346655
n°13	¹⁵ N	¹⁷ O W1346697	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346715	¹ H W1346655

900MHz Triple Resonance Solid Experiments

n°1	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346323	⁵⁹ Co	¹³ C W1346315	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346299
n°2	¹⁵ N	¹⁷ O	² H W1346327	²⁹ Si	⁵⁹ Co	¹³ C W1346315	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346299
n°3	¹⁵ N W1346662	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1346319	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346299
n°4	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346323	⁵⁹ Co	¹³ C	²⁷ Al W1346315	²³ Na	¹¹ B	³¹ P	¹ H W1346299
n°5	¹⁵ N	¹⁷ O W1346327	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346311	³¹ P	¹ H W1346299
n°6	¹⁵ N W1346331	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346311	³¹ P	¹ H W1346299
n°7	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na W1346315	¹¹ B	³¹ P W1346307	¹ H W1346299
n°8	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al W1346315	²³ Na	¹¹ B	³¹ P W1346307	¹ H W1346299
n°9	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1346315	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346307	¹ H W1346299
n°10	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346307	¹ H W1346299
n°11	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346323	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346307	¹ H W1346299
n°12	¹⁵ N	¹⁷ O	² H W1346327	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346307	¹ H W1346299
n°13	¹⁵ N	¹⁷ O W1346327	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346307	¹ H W1346299

950MHz Triple Resonance Solid Experiments

n°1	¹⁵ N	¹⁷ O	² H	²⁹ Si W140203	⁵⁹ Co	¹³ C W140204	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346707
n°2	¹⁵ N	¹⁷ O	² H W139785	²⁹ Si	⁵⁹ Co	¹³ C W140204	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346707
n°3	¹⁵ N W143787	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na W143783	¹¹ B	³¹ P	¹ H W1346707
n°4	¹⁵ N	¹⁷ O	² H	²⁹ Si W140203	⁵⁹ Co	¹³ C	²⁷ Al W140204	²³ Na	¹¹ B	³¹ P	¹ H W1346707
n°5	¹⁵ N	¹⁷ O W139785	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W140205	³¹ P	¹ H W1346707
n°6	¹⁵ N W143779	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W140205	³¹ P	¹ H W1346707
n°7	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na W140204	¹¹ B	³¹ P W143776	¹ H W1346707
n°8	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al W140204	²³ Na	¹¹ B	³¹ P W143776	¹ H W1346707
n°9	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W140204	²⁷ Al	²³ Na	¹¹ B	³¹ P W143776	¹ H W1346707
n°10	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co W143783	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W143776	¹ H W1346707
n°11	¹⁵ N	¹⁷ O	² H	²⁹ Si W140203	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W143776	¹ H W1346707
n°12	¹⁵ N	¹⁷ O	² H W139785	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W143776	¹ H W1346707
n°13	¹⁵ N	¹⁷ O W139785	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W143776	¹ H W1346707

1000MHz Triple Resonance Solid Experiments

n°1	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346670	⁵⁹ Co	¹³ C W1346669	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346665
n°2	¹⁵ N	¹⁷ O	² H W1346671	²⁹ Si	⁵⁹ Co	¹³ C W1346669	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346665
n°3	¹⁵ N W1346672	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1346669	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1346665
n°4	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346670	⁵⁹ Co	¹³ C	²⁷ Al W1346669	²³ Na	¹¹ B	³¹ P	¹ H W1346665
n°5	¹⁵ N	¹⁷ O W1346671	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346668	³¹ P	¹ H W1346665
n°6	¹⁵ N W1346717	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346668	³¹ P	¹ H W1346665
n°7	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na W1346669	¹¹ B	³¹ P W1346667	¹ H W1346665
n°8	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al W1346669	²³ Na	¹¹ B	³¹ P W1346667	¹ H W1346665
n°9	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1346669	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346667	¹ H W1346665
n°10	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co W1346856	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346667	¹ H W1346665
n°11	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346670	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346667	¹ H W1346665
n°12	¹⁵ N	¹⁷ O	² H W1346671	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346667	¹ H W1346665
n°13	¹⁵ N	¹⁷ O W1346671	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346667	¹ H W1346665

1100MHz Triple Resonance Solid Experiments

n°1	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346863	⁵⁹ Co	¹³ C W1348165	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1342132
n°2	¹⁵ N	¹⁷ O	² H W1346865	²⁹ Si	⁵⁹ Co	¹³ C W1348165	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1342132
n°3	¹⁵ N W1342624	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1348165	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1342132
n°4	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346863	⁵⁹ Co	¹³ C	²⁷ Al W1348165	²³ Na	¹¹ B	³¹ P	¹ H W1342132
n°5	¹⁵ N	¹⁷ O W1346865	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346861	³¹ P	¹ H W1342132
n°6	¹⁵ N W1346867	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346861	³¹ P	¹ H W1342132
n°7	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na W1348165	¹¹ B	³¹ P W1346859	¹ H W1342132
n°8	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al W1348165	²³ Na	¹¹ B	³¹ P W1346859	¹ H W1342132
n°9	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W1348165	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346859	¹ H W1342132
n°10	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co W1346857	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346859	¹ H W1342132
n°11	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346863	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346859	¹ H W1342132
n°12	¹⁵ N	¹⁷ O	² H W1346865	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346859	¹ H W1342132
n°13	¹⁵ N	¹⁷ O W1346865	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346859	¹ H W1342132

1200MHz Triple Resonance Solid Experiments

n°1	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346864	⁵⁹ Co	¹³ C W135028	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1348174
n°2	¹⁵ N	¹⁷ O	² H W1346866	²⁹ Si	⁵⁹ Co	¹³ C W135028	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1348174
n°3	¹⁵ N W134926	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W135028	²⁷ Al	²³ Na	¹¹ B	³¹ P	¹ H W1348174
n°4	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346864	⁵⁹ Co	¹³ C	²⁷ Al W135028	²³ Na	¹¹ B	³¹ P	¹ H W1348174
n°5	¹⁵ N	¹⁷ O W1346866	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346862	³¹ P	¹ H W1348174
n°6	¹⁵ N W1346868	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B W1346862	³¹ P	¹ H W1348174
n°7	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na W135028	¹¹ B	³¹ P W1346860	¹ H W1348174
n°8	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al W135028	²³ Na	¹¹ B	³¹ P W1346860	¹ H W1348174
n°9	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co	¹³ C W135028	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346860	¹ H W1348174
n°10	¹⁵ N	¹⁷ O	² H	²⁹ Si	⁵⁹ Co W1346858	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346860	¹ H W1348174
n°11	¹⁵ N	¹⁷ O	² H	²⁹ Si W1346864	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346860	¹ H W1348174
n°12	¹⁵ N	¹⁷ O	² H W1346866	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346860	¹ H W1348174
n°13	¹⁵ N	¹⁷ O W1346866	² H	²⁹ Si	⁵⁹ Co	¹³ C	²⁷ Al	²³ Na	¹¹ B	³¹ P W1346860	¹ H W1348174

1.2.2 Triple Resonance $^1\text{H}/\text{X}/\text{Y}$



Figure 1.14: $^1\text{H}/\text{X}/\text{Y}$ Typical Configuration

This picture shows typical configuration of high power filters required for triple resonance solid-state NMR experiments ($^1\text{H}/\text{X}/\text{Y}$).

The two heavy band-pass filters, corresponding to X and Y nuclei studied, are not directly connected to the pre-amplifier (additional cables are provided). The ^1H pass filter is only connected to the pre-amplifier.

1.2.3 Triple resonance $^1\text{H}/^{19}\text{F}/\text{X}$




Figure 1.15: $^1\text{H}/^{19}\text{F}/\text{X}$ NMR Experiments Configuration

This picture shows typical configuration of high power filters required for $^1\text{H}/^{19}\text{F}$ NMR experiments. The X channel is not shown.

A schematic representation is attached here for the description of the different high power filters involved in the $^1\text{H}/^{19}\text{F}$ experiment.




Figure 1.16: Different High Power Filters for the $^1\text{H}/^{19}\text{F}$ Experiment

This diagram corresponds to ^{19}F observation with high power ^1H decoupling NMR experiment (^1H pre-amplifier is not required). The order of the different Pass/Stop high power filters is important.

In the case of ^1H observation with ^{19}F decoupling NMR experiment, the $^1\text{H}/^{19}\text{F}$ HP pre-amplifier is dedicated to ^1H channel (not required for ^{19}F decoupling).

1.3 AVANCE 300MHz

1.3.1 Lowpass Filter




Figure 1.17: Lowpass Filter HQ 300MHz Diagram

W1346642 - Filter HQ 300 0- ³¹ P LP (¹⁹ F- ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ³¹ P	< 121.5	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹⁹ F... ¹ H	282.4...300.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D050 [226] (mm)	144 x 63 x 34	
Weight (kg)	0.4	

Table 1.1: W1346642 - Filter HQ 300 0-³¹P LP (¹⁹F-¹H)

W135032 - Filter HQ 300 O - ¹⁵ N NR LP (² H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) :... ¹⁵ N	< 30.4	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)		
Frequency Stop (MHz) :... ² H	46.1	
Minimum Rejection (dB)		
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.2: W135032 - Filter HQ 300 O -¹⁵N NR LP (²H)

1.3.2 Bandstop Filter




Figure 1.18: Bandstop Filter HQ 300MHz Diagram

W1346253 - Filter HQ 300 ¹ H S		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ³¹ P, ¹⁹ F	< 121.5	282.4
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹ H	300.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.3: W1346253 - Filter HQ 300 ¹H S

W1346369 - Filter HQ 300 ¹⁹ F S		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ³¹ P, ¹ H	< 121.5	300.1
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹⁹ F	282.4	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.4: W1346369 - Filter HQ 300 ¹⁹F S

1.3.3 Bandpass Filter




Figure 1.19: Bandpass Filter HQ 300MHz Diagram

W1346239 - Filter HQ 300 ¹ H BP (⁻³¹ P, ¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H	300.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ³¹ P, ¹⁹ F	< 121.5	282.4
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R100 [▶ 217] (mm)	227 x 123 x 39	
Weight (kg)	1.2	

Table 1.5: W1346239 - Filter HQ 300 ¹H BP (⁻³¹P, ¹⁹F)

W1346235 - Filter HQ 300 ¹ H BP (¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H	300.1	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ³¹ P, ¹⁹ F	< 121.5	282.4
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 8R105P ▶ 219 (mm)	277 x 133 x 54	
Weight (kg)	2.6	

Table 1.6: W1346235 - Filter HQ 300 ¹H BP (¹⁹F)

W1346249 - Filter HQ 300 ¹⁹ F BP (⁻³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁹ F	282.4	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ³¹ P, ¹ H	< 121.5	300.1
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R100 ▶ 217 (mm)	227 x 123 x 39	
Weight (kg)	1.2	

Table 1.7: W1346249 - Filter HQ 300 ¹⁹F BP (⁻³¹P, ¹H)

W1346236 - Filter HQ 300 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁹ F	282.4	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ³¹ P, ¹ H	< 121.5	300.1
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 8R105P [▶ 219] (mm)	277 x 133 x 54	
Weight (kg)	2.6	

Table 1.8: W1346236 - Filter HQ 300 ¹⁹F BP (¹H)

W1346368 - Filter HQ 300 ¹⁹ F - ¹ H BP (- ³¹ P)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁹ F, ¹ H	282.4	300.1
Maximum Insertion Loss (dB)	0.8	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ³¹ P	< 121.5	
Minimum Rejection (dB)	110	
Mechanical Dimensions and Weight		
Case 5R100 [▶ 217] (mm)	227 x 123 x 39	
Weight (kg)	1.2	

Table 1.9: W1346368 - Filter HQ 300 ¹⁹F - ¹H BP (-³¹P)

W1346372 - Filter HQ 300 ⁷ Li - ³¹ P BP (- ¹³ C, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁷ Li... ³¹ P	116.6...121.5	
Maximum Insertion Loss (dB)	0.5	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ¹³ C, ¹⁹ F...	< 75.5	> 282.4
Minimum Rejection (dB)	90	80
Mechanical Dimensions and Weight		
Case 6H073 [▶ 222] (mm)	277 x 100 x 38	
Weight (kg)	1.6	

Table 1.10: W1346372 - Filter HQ 300 ⁷Li - ³¹P BP (-¹³C, ¹⁹F-)

W1346270 - Filter HQ 300 ⁷¹ Ga - ⁸⁷ Rb BP (- ¹³ C, ³¹ P-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁷¹ Ga... ⁸⁷ Rb	91.5...98.2	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ¹³ C, ³¹ P...	< 75.5	> 121.5
Minimum Rejection (dB)	90	85
Mechanical Dimensions and Weight		
Case 8H097 [▶ 223] (mm)	336 x 125 x 38	
Weight (kg)	2.5	

Table 1.11: W1346270 - Filter HQ 300 ⁷¹Ga - ⁸⁷Rb BP (-¹³C, ³¹P-)

W1346271 - Filter HQ 300 ⁵⁹ Co - ²³ Na BP (⁻²⁹ Si, ⁸⁷ Rb-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁵⁹ Co... ²³ Na	71.2...79.4	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ²⁹ Si, ⁸⁷ Rb..	< 59.6	> 80
Minimum Rejection (dB)	80	80
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.12: W1346271 - Filter HQ 300 ⁵⁹Co - ²³Na BP (⁻²⁹Si, ⁸⁷Rb-)

W1346275 - Filter HQ 300 ⁵⁹ Co - ²³ Na BP (⁻¹⁵ N)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁵⁹ Co... ²³ Na	71.2...79.4	
Maximum Insertion Loss (dB)	0.5	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹⁵ N	< 30.4	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 4H073 [▶ 221] (mm)	205 x 100 x 38	
Weight (kg)	1.1	

Table 1.13: W1346275 - Filter HQ 300 ⁵⁹Co - ²³Na BP (⁻¹⁵N)

W1346272 - Filter HQ 300 ²⁹ Si BP (- ² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ²⁹ Si	59.6	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ² H, ¹³ C...	< 46.0	> 75.5
Minimum Rejection (dB)	90	80
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.14: W1346272 - Filter HQ 300 ²⁹Si BP (-²H, ¹³C-)

W1346452 - Filter HQ 300 ¹⁷ O - ² H BP (¹¹ B-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁷ O... ² H	40.7...46.0	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹¹ B, ³¹ P...	> 96.3	
Minimum Rejection (dB)	75	80
Mechanical Dimensions and Weight		
Case 5H169 (mm)	241 x 197 x 38	
Weight (kg)	2.5	

Table 1.15: W1346452 - Filter HQ 300 ¹⁷O - ²H BP (¹¹B-)

W1346273 - Filter HQ 300 ¹³³ Cs - ² H BP (⁻¹⁵ N, ²⁹ Si-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : : ¹³³ Cs... ² H	39.4...46.0	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ¹⁵ N, ²⁹ Si...	< 30.4	> 59.6
Minimum Rejection (dB)	90	80
Mechanical Dimensions and Weight		
Case 8H169 [▶ 225] (mm)	336 x 197 x 38	
Weight (kg)	3.9	

Table 1.16: W1346273 - Filter HQ 300 ¹³³Cs - ²H BP (⁻¹⁵N, ²⁹Si-)

W1346274 - Filter HQ 300 ¹⁵ N BP (⁻¹⁴ N, ² H-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁵ N	30.4	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ¹⁴ N, ² H...	< 21.7	> 46.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H169 [▶ 225] (mm)	336 x 197 x 38	
Weight (kg)	3.9	

Table 1.17: W1346274 - Filter HQ 300 ¹⁵N BP (⁻¹⁴N, ²H-)

W1346278 - Filter HQ 300 ¹⁵ N BP (¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁵ N	30.4	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) ⁵⁹ Co, ¹³ C...:	71.2	> 75.5
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 5H169 (mm)	241 x 197 x 38	
Weight (kg)	2.5	

Table 1.18: W1346278 - Filter HQ 300 ¹⁵N BP (¹³C-)

1.3.4 Diplexer Filter




Figure 1.20: Diplexer Filter HQ 300MHz Diagram

W1346237 - Diplexer HQ 300 ¹ H / ¹⁹ F		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H, ¹⁹ F	300.1	282.4
Maximum Insertion Loss (dB)	0.6	0.6
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹⁹ F, ¹ H	282.4	300.1
Minimum Rejection (dB)	50	50
Mechanical Dimensions and Weight		
Case 4S063D [▶ 215] (mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.19: W1346237 - Diplexer HQ 300 ¹H / ¹⁹F

1.4 AVANCE 400MHz

1.4.1 Lowpass Filter




Figure 1.21: Lowpass Filter HQ 400MHz Diagram

W1346636 - Filter HQ 400 0- ³¹ P LP (¹⁹ F- ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ³¹ P	< 162.0	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹⁹ F... ¹ H	376.5...400.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D050 [▶ 226] (mm)	144 x 63 x 34	
Weight (kg)	0.4	

Table 1.20: W1346636 - Filter HQ 400 0-³¹P LP (¹⁹F-¹H)

W1346628 - Filter HQ 400 0- ¹⁵ N NR LP (¹³ C)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ¹⁵ N	< 40.6	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹³ C	100.6	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.21: W1346628 - Filter HQ 400 0-¹⁵N NR LP (¹³C)

1.4.2 Bandstop Filter




Figure 1.22: Bandstop Filter HQ 400MHz Diagram

W1346370 - Filter HQ 400 ¹⁹ F S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ³¹ P, ¹ H	< 162.0	400.1
Maximum Insertion Loss (dB)	0.4	0.4
Minimum Return Loss (dB)	15	20
Frequency Stop (MHz) : ¹⁹ F	376.5	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.22: W1346370 - Filter HQ 400 ¹⁹F S

W1346437 - Filter HQ 400 ² H P ⁷⁷ Se S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ² H	61.4	
Maximum Insertion Loss (dB)	0.4	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ⁷⁷ Se	76.3	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.23: W1346437 - Filter HQ 400 ²H P ⁷⁷Se S

1.4.3 Bandpass Filter




Figure 1.23: Bandpass Filter HQ 400MHz Diagram

W1346240 - Filter HQ 400 ¹ H BP (⁻³¹ P, ¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H	400.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹⁹ F	< 162.0	376.5
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R100 [▶ 217] (mm)	227 x 123 x 39	
Weight (kg)	1.2	

Table 1.24: W1346240 - Filter HQ 400 ¹H BP (⁻³¹P, ¹⁹F)

W1346204 - Filter HQ 400 ¹ H BP (¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H	400.1	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹⁹ F	< 162.0	376.5
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 8R105P [▶ 219] (mm)	277 x 133 x 54	
Weight (kg)	2.6	

Table 1.25: W1346204 - Filter HQ 400 ¹H BP (¹⁹F)

W1346250 - Filter HQ 400 ¹⁹ F BP (⁻³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁹ F	376.5	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹ H	< 162.0	400.1
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R100 [▶ 217] (mm)	227 x 123 x 39	
Weight (kg)	1.2	

Table 1.26: W1346250 - Filter HQ 400 ¹⁹F BP (⁻³¹P, ¹H)

W1346205 - Filter HQ 400 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁹ F	376.5	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹ H	< 162.0	400.1
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 8R105P ▶ 219 (mm)	277 x 133 x 54	
Weight (kg)	2.6	

Table 1.27: W1346205 - Filter HQ 400 ¹⁹F BP (¹H)

W1346373 - Filter HQ 400 ¹⁹ F - ¹ H BP (⁻³¹ P)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁹ F, ¹ H	376.5	400.1
Maximum Insertion Loss (dB)	0.8	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ... ³¹ P	< 162.0	
Minimum Rejection (dB)	110	
Mechanical Dimensions and Weight		
Case 5R100 ▶ 217 (mm)	227 x 123 x 39	
Weight (kg)	1.2	

Table 1.28: W1346373 - Filter HQ 400 ¹⁹F - ¹H BP (⁻³¹P)

W1346256 - Filter HQ 400 ⁷ Li - ³¹ P BP (⁻⁸⁷ Rb, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁷ Li... ³¹ P	155.5...162.0	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ⁸⁷ Rb, ¹⁹ F...	< 130.9	> 376.5
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.29: W1346256 - Filter HQ 400 ⁷Li - ³¹P BP (⁻⁸⁷Rb, ¹⁹F-)

W1346349 - Filter HQ 400 ⁷ Li - ³¹ P BP (⁻¹³ C, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁷ Li... ³¹ P	155.5...162.0	
Maximum Insertion Loss (dB)	0.5	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ¹³ C, ¹⁹ F...	< 100.6	> 376.5
Minimum Rejection (dB)	90	80
Mechanical Dimensions and Weight		
Case 6H073 [▶ 222] (mm)	277 x 100 x 38	
Weight (kg)	1.6	

Table 1.30: W1346349 - Filter HQ 400 ⁷Li - ³¹P BP (⁻¹³C, ¹⁹F-)

W1346257 - Filter HQ 400 ⁷¹ Ga - ⁸⁷ Rb BP (- ¹³ C, ³¹ P-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁷¹ Ga... ⁸⁷ Rb	122.0...130.9	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ¹³ C, ³¹ P...	< 100.6	> 162.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.31: W1346257 - Filter HQ 400 ⁷¹Ga - ⁸⁷Rb BP (-¹³C, ³¹P-)

W1346488 - Filter HQ 400 ¹³ C - ²³ Na BP (- ²⁹ Si, ⁷¹ Ga-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹³ C... ²³ Na	100.6...105.8	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ²⁹ Si, ⁷¹ Ga...	< 79.5	> 122.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.32: W1346488 - Filter HQ 400 ¹³C - ²³Na BP (-²⁹Si, ⁷¹Ga-)

W1346258 - Filter HQ 400 ⁵⁹ Co - ²³ Na BP (- ²⁹ Si, ⁸⁷ Rb-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁵⁹ Co... ²³ Na	94.9...105.8	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ²⁹ Si, ⁸⁷ Rb..	< 79.5	> 130.9
Minimum Rejection (dB)	80	80
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.33: W1346258 - Filter HQ 400 ⁵⁹Co - ²³Na BP (-²⁹Si, ⁸⁷Rb-)

W1346440 - Filter HQ 400 ⁵⁹ Co - ²³ Na BP (- ² H, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁵⁹ Co... ²³ Na	94.9...105.8	
Maximum Insertion Loss (dB)	0.6	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ² H, ¹⁹ F...	< 61.4	> 376.5
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6H073 [▶ 222] (mm)	277 x 100 x 38	
Weight (kg)	1.6	

Table 1.34: W1346440 - Filter HQ 400 ⁵⁹Co - ²³Na BP (-²H, ¹⁹F-)

W1346226 - Filter HQ 400 ⁵⁹ Co - ²³ Na BP (⁻¹⁵ N)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁵⁹ Co... ²³ Na	94.9...105.8	
Maximum Insertion Loss (dB)	0.5	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ¹⁵ N	< 40.6	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 4H073 [▶ 221] (mm)	205 x 100 x 38	
Weight (kg)	1.1	

Table 1.35: W1346226 - Filter HQ 400 ⁵⁹Co - ²³Na BP (⁻¹⁵N)

W1346259 - Filter HQ 400 ²⁹ Si BP (⁻² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ²⁹ Si	79.5	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ² H, ¹³ C...	< 61.4	> 100.6
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.36: W1346259 - Filter HQ 400 ²⁹Si BP (⁻²H, ¹³C-)

W1346442 - Filter HQ 400 ¹⁷ O - ² H BP (¹¹ B-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁷ O... ² H	54.2...61.4	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹¹ B, ³¹ P...	> 128.3	
Minimum Rejection (dB)	75	80
Mechanical Dimensions and Weight		
Case 5H169 (mm)	241 x 197 x 38	
Weight (kg)	2.5	

Table 1.37: W1346442 - Filter HQ 400 ¹⁷O - ²H BP (¹¹B-)

W1346260 - Filter HQ 400 ¹³³ Cs - ² H BP (¹⁵ N, ²⁹ Si-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹³³ Cs... ² H	52.5...61.4	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ¹⁵ N, ²⁹ Si...	< 40.6	> 79.5
Minimum Rejection (dB)	90	80
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.38: W1346260 - Filter HQ 400 ¹³³Cs - ²H BP (¹⁵N, ²⁹Si-)

W1346261 - Filter HQ 400 ¹⁵N BP (-¹⁴N, ²H-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁵ N	40.6	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ¹⁴ N, ² H...	< 28.9	> 61.4
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.39: W1346261 - Filter HQ 400 ¹⁵N BP (-¹⁴N, ²H-)

W1346439 - Filter HQ 400 ²³N ⁶⁵CU BP (-²⁹SI, ³¹P-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ²³ N... ⁶⁵ CU	105.8...113.7	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ²⁹ SI, ³¹ P...	< 79.5	> 162.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.40: W1346439 - Filter HQ 400 ²³N ⁶⁵CU BP (-²⁹SI, ³¹P-)

1.4.4 Diplexer Filter




Figure 1.24: Diplexer Filter HQ 400MHz Diagram

W1346345 - Diplexer HQ 400 ¹ H / ¹⁹ F		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H, ¹⁹ F	400.1	376.5
Maximum Insertion Loss (dB)	0.6	0.6
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹⁹ F, ¹ H	376.5	400.1
Minimum Rejection (dB)	70	70
Mechanical Dimensions and Weight		
Case 4S063D [▶ 215] (mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.41: W1346345 - Diplexer HQ 400 ¹H / ¹⁹F

1.5 AVANCE 500MHz

1.5.1 Lowpass Filter




Figure 1.25: Lowpass Filter HQ 500MHz Diagram

W1346637 - Filter HQ 500 0- ³¹ P LP (¹⁹ F- ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ³¹ P	< 202.5	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹⁹ F... ¹ H	470.6...500.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D050 [▶ 226] (mm)	144 x 63 x 34	
Weight (kg)	0.4	

Table 1.42: W1346637 - Filter HQ 500 0-³¹P LP (¹⁹F-¹H)

W1346629 - Filter HQ 500 0- ¹⁵ N NR LP (¹³ C)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ¹⁵ N	< 50.7	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹³ C	125.8	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.43: W1346629 - Filter HQ 500 0-¹⁵N NR LP (¹³C)

1.5.2 Bandstop Filter




Figure 1.26: Bandstop Filter HQ 500MHz Diagram

W1346232 - Filter HQ 500 ¹ H S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ³¹ P, ¹⁹ F	< 202.5	470.6
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹ H	500.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.44: W1346232 - Filter HQ 500 ¹H S

W1346358 - Filter HQ 500 ¹⁹ F S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ³¹ P, ¹ H	< 202.5	500.1
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹⁹ F	470.6	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.45: W1346358 - Filter HQ 500 ¹⁹F S

W1346430 - Filter HQ 500 ⁷ Li P ³¹ P S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁷ Li	194.4	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) ³¹ P	202.5	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.46: W1346430 - Filter HQ 500 ⁷Li P ³¹P S

W1346431 - Filter HQ 500 ³¹ P P ⁷ Li S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ³¹ P	202.5	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ⁷ Li	194.4	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.47: W1346431 - Filter HQ 500 ³¹P P ⁷Li S

W1346434 - Filter HQ 500 ⁶ Li P ² H S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁶ Li	73.6	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ² H	76.8	
Minimum Rejection (dB)	80	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.48: W1346434 - Filter HQ 500 ⁶Li P ²H S

W1346435 - Filter HQ 500 ² H P ⁶ Li S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ² H	76.8	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ⁶ Li	73.6	
Minimum Rejection (dB)	80	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.49: W1346435 - Filter HQ 500 ²H P ⁶Li S

W130896- Filter HQ 500 ¹³ C-P ²⁷ Al- S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹³ C	125.8	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : .. ²⁷ Al	130.3	
Minimum Rejection (dB)	80	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.50: W130896- Filter HQ 500 ¹³C-P ²⁷Al- S

W130897- Filter HQ 500 ²⁷ Al-P ¹³ C- S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ²⁷ Al	130.3	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : .. ¹³ C	125.8	
Minimum Rejection (dB)	80	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.51: W130897- Filter HQ 500 ²⁷Al-P ¹³C- S

1.5.3 Bandpass Filter




Figure 1.27: Band Pass Filter HQ 500MHz Diagram

W1346223 - Filter HQ 500 ¹ H BP (⁻³¹ P, ¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H	500.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹⁹ F	< 202.5	470.6
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R069 [▶ 217] (mm)	227 x 92 x 39	
Weight (kg)	0.9	

Table 1.52: W1346223 - Filter HQ 500 ¹H BP (⁻³¹P, ¹⁹F)

W1346247 - Filter HQ 500 ¹ H BP (¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H	500.1	
Maximum Insertion Loss (dB)	0.9	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹⁹ F	< 202.5	470.6
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R055P [▶ 219] (mm)	282 x 83 x 54	
Weight (kg)	1.6	

Table 1.53: W1346247 - Filter HQ 500 ¹H BP (¹⁹F)

W1346251 - Filter HQ 500 ¹⁹ F BP (⁻³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁹ F	470.6	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹ H	< 202.5	500.1
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R069 [▶ 217] (mm)	227 x 92 x 39	
Weight (kg)	0.9	

Table 1.54: W1346251 - Filter HQ 500 ¹⁹F BP (⁻³¹P, ¹H)

W1346248 - Filter HQ 500 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁹ F	470.6	
Maximum Insertion Loss (dB)	0.9	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹ H	< 202.5	500.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R055P [▶ 219] (mm)	282 x 83 x 54	
Weight (kg)	1.6	

Table 1.55: W1346248 - Filter HQ 500 ¹⁹F BP (¹H)

W1346374 - Filter HQ 500 ¹⁹ F - ¹ H BP (- ³¹ P)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁹ F, ¹ H	470.6	500.1
Maximum Insertion Loss (dB)	0.8	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ... ³¹ P	< 202.5	
Minimum Rejection (dB)	110	
Mechanical Dimensions and Weight		
Case 5R069 [▶ 217] (mm)	227 x 92 x 39	
Weight (kg)	0.9	

Table 1.56: W1346374 - Filter HQ 500 ¹⁹F - ¹H BP (-³¹P)

W1346262 - Filter HQ 500 ⁷ Li - ³¹ P BP (- ⁸⁷ Rb, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁷ Li... ³¹ P	194.4...202.5	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ⁸⁷ Rb, ¹⁹ F...	< 163.6	> 470.6
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.57: W1346262 - Filter HQ 500 ⁷Li - ³¹P BP (-⁸⁷Rb, ¹⁹F-)

W1346351 - Filter HQ 500 ⁷ Li - ³¹ P BP (- ¹³ C, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁷ Li... ³¹ P	194.4...202.5	
Maximum Insertion Loss (dB)	0.5	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ¹³ C, ¹⁹ F...	< 125.8	> 470.6
Minimum Rejection (dB)	80	70
Mechanical Dimensions and Weight		
Case 6H073 [▶ 222] (mm)	277 x 100 x 38	
Weight (kg)	1.6	

Table 1.58: W1346351 - Filter HQ 500 ⁷Li - ³¹P BP (-¹³C, ¹⁹F-)

W1346263 - Filter HQ 500 ⁷¹ Ga - ⁸⁷ Rb BP (- ¹³ C, ³¹ P-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁷¹ Ga... ⁸⁷ Rb	152.5...163.6	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ¹³ C, ³¹ P...	< 125.8	> 202.5
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H097 [▶ 223] (mm)	336 x 125 x 38	
Weight (kg)	2.5	

Table 1.59: W1346263 - Filter HQ 500 ⁷¹Ga - ⁸⁷Rb BP (-¹³C, ³¹P-)

W1346489 - Filter HQ 500 ¹³ C - ²³ Na BP (- ²⁹ Si, ⁷¹ Ga-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹³ C... ²³ Na	125.8...132.3	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ²⁹ Si, ⁷¹ Ga...	< 99.4	> 152.5
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.60: W1346489 - Filter HQ 500 ¹³C - ²³Na BP (-²⁹Si, ⁷¹Ga-)

W1346264 - Filter HQ 500 ⁵⁹ Co - ²³ Na BP (- ²⁹ Si, ⁸⁷ Rb-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁵⁹ Co... ²³ Na	118.7...132.3	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ²⁹ Si, ⁸⁷ Rb..	< 99.4	> 163.6
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H097 [▶ 223] (mm)	336 x 125 x 38	
Weight (kg)	2.5	

Table 1.61: W1346264 - Filter HQ 500 ⁵⁹Co - ²³Na BP (-²⁹Si, ⁸⁷Rb-)

W1346460 - Filter HQ 500 ⁵⁹ Co - ²³ Na BP (- ² H, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁵⁹ Co... ²³ Na	118.7...132.3	
Maximum Insertion Loss (dB)	0.6	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ² H, ¹⁹ F...	< 76.8	> 470.6
Minimum Rejection (dB)	90	80
Mechanical Dimensions and Weight		
Case 6H073 [▶ 222] (mm)	277 x 100 x 38	
Weight (kg)	1.6	

Table 1.62: W1346460 - Filter HQ 500 ⁵⁹Co - ²³Na BP (-²H, ¹⁹F-)

W1346224 - Filter HQ 500 ⁵⁹ Co - ²³ Na BP (- ¹⁵ N)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁵⁹ Co... ²³ Na	118.7...132.3	
Maximum Insertion Loss (dB)	0.5	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ¹⁵ N	< 50.7	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 4H073 [▶ 221] (mm)	205 x 100 x 38	
Weight (kg)	1.1	

Table 1.63: W1346224 - Filter HQ 500 ⁵⁹Co - ²³Na BP (-¹⁵N)

W1346265 - Filter HQ 500 ⁷⁷ Se - ²⁹ Si BP (- ² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁷⁷ Se..... ²⁹ Si	95.4.....99.4	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ² H, ¹³ C...	< 76.8	> 125.8
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.64: W1346265 - Filter HQ 500 ⁷⁷Se - ²⁹Si BP (-²H, ¹³C-)

W1346453 - Filter HQ 500 ¹⁷ O - ² H BP (¹¹ B-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁷ O... ² H	67.8...76.7	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹¹ B, ³¹ P...	> 160.4	
Minimum Rejection (dB)	75	85
Mechanical Dimensions and Weight		
Case 5H169 (mm)	241 x 197 x 38	
Weight (kg)	2.5	

Table 1.65: W1346453 - Filter HQ 500 ¹⁷O - ²H BP (¹¹B-)

W1346266 - Filter HQ 500 ¹³³ Cs - ² H BP (⁻¹⁵ N, ²⁹ Si-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹³³ Cs... ² H	65.6...76.8	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ¹⁵ N, ²⁹ Si...	< 50.7	> 99.4
Minimum Rejection (dB)	90	80
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.66: W1346266 - Filter HQ 500 ¹³³Cs - ²H BP (⁻¹⁵N, ²⁹Si-)

W1346267 - Filter HQ 500 ¹⁵ N BP (⁻¹⁴ N, ² H-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁵ N	50.7	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ¹⁴ N, ² H...	< 36.1	> 76.8
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.67: W1346267 - Filter HQ 500 ¹⁵N BP (⁻¹⁴N, ²H-)

W1346692 - Filter HQ 500 ¹ H BP (³ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H	500.1	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : .. ³ H	533.5	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 8R055P [▶ 219] (mm)	282 x 83 x 54	
Weight (kg)	1.6	

Table 1.68: W1346692 - Filter HQ 500 ¹H BP (³H)

W1346693 - Filter HQ 500 ¹ H BP (⁻³¹ P, ³ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H	500.1	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : .. ³¹ P, ³ H...	< 202.5	533.5
Minimum Rejection (dB)	110	50
Mechanical Dimensions and Weight		
Case 5R069 [▶ 217] (mm)	227 x 92 x 39	
Weight (kg)	0.9	

Table 1.69: W1346693 - Filter HQ 500 ¹H BP (⁻³¹P, ³H)

W1346694 - Filter HQ 500 ³ H BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ³ H	533.5	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹ H	500.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 8R055P [▶ 219] (mm)	282 x 83 x 54	
Weight (kg)	1.6	

Table 1.70: W1346694 - Filter HQ 500 ³H BP (¹H)

W1346695 - Filter HQ 500 ³ H BP (- ³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ³ H	533.5	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹ H...	< 202.5	500.1
Minimum Rejection (dB)	110	50
Mechanical Dimensions and Weight		
Case 5R069 [▶ 217] (mm)	227 x 92 x 39	
Weight (kg)	0.9	

Table 1.71: W1346695 - Filter HQ 500 ³H BP (-³¹P, ¹H)

W1346730 - Filter HQ 500 ²³ Na - ⁶⁵ Cu BP (- ²⁹ Si, ³¹ P-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ²³ Na, ⁶⁵ Cu...	132.3...142.1	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ²⁹ Si, ³¹ P...	< 199.4	> 202.5
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.72: W1346730 - Filter HQ 500 ²³Na -⁶⁵Cu BP (-²⁹Si, ³¹P-)

1.5.4 Diplexer Filter




Figure 1.28: Diplexer Filter HQ 500MHz Diagram

W1346346 - Diplexer HQ 500 ¹ H / ¹⁹ F		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H, ¹⁹ F	500.1	470.6
Maximum Insertion Loss (dB)	0.6	0.6
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹⁹ F, ¹ H	470.6	500.1
Minimum Rejection (dB)	70	70
Mechanical Dimensions and Weight		
Case 4S063D [▶ 215] (mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.73: W1346346 - Diplexer HQ 500 ¹H / ¹⁹F

W1346696 - Diplexer HQ 500 ¹ H / ³ H		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H, ³ H	500.1	533.5
Maximum Insertion Loss (dB)	0.6	0.6
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ³ H, ¹ H	533.5	500.1
Minimum Rejection (dB)	50	50
Mechanical Dimensions and Weight		
Case 4S063D [p 215] (mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.74: W1346696 - Diplexer HQ 500 ¹H / ³H

W136807- Diplexer HQ 500 ¹³ C / ²⁷ Al		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹³ C, ²⁷ Al	125.8	130.3
Maximum Insertion Loss (dB)	0.8	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ²⁷ Al, ¹³ C	130.3	125.8
Minimum Rejection (dB)	50	50
Mechanical Dimensions and Weight		
Case 4S063D [p 215] (mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.75: W136807- Diplexer HQ 500 ¹³C / ²⁷Al

1.6 AVANCE 600MHz

1.6.1 Lowpass Filter




Figure 1.29: Lowpass Filter HQ 600 MHz Diagram

W1346638 - Filter HQ 600 0- ³¹ P LP (¹⁹ F- ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ³¹ P	< 242.9	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹⁹ F... ¹ H	564.7...600.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D050 [▶ 226] (mm)	144 x 63 x 34	
Weight (kg)	0.4	

Table 1.76: W1346638 - Filter HQ 600 0-³¹P LP (¹⁹F-¹H)

W1346630 - Filter HQ 600 0- ¹⁵ N NR LP (¹³ C)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ¹⁵ N	< 60.8	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹³ C	150.9	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.77: W1346630 - Filter HQ 600 0-¹⁵N NR LP (¹³C)

W134543 - Filter HQ 600 0- ¹⁴ N NR (109Ag H ³) LP (¹³ C)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁴ N	< 43.3	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz): ¹³ C	150.9	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.78: W134543 - Filter HQ 600 0-¹⁴N NR (109Ag H³) LP (¹³C)

1.6.2 Bandstop Filter




Figure 1.30: Bandstop Filter HQ 600 MHz Diagram

W1346255 - Filter HQ 600 ¹ H S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ³¹ P, ¹⁹ F	< 242.9	564.7
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹ H	600.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.79: W1346255 - Filter HQ 600 ¹H S

W1346371 - Filter HQ 600 ¹⁹ F S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ... ³¹ P, ¹ H	< 242.9	600.1
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹⁹ F	564.7	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.80: W1346371 - Filter HQ 600 ¹⁹F S

W1346679 - Filter HQ 600 ⁶ Li- P ² H- S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁶ Li	88.3	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz): ² H	92.1	
Minimum Rejection (dB)	80	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.81: W1346679 - Filter HQ 600 ⁶Li- P ²H- S

W1346680 - Filter HQ 600 ² H- P ⁶ Li- S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ² H	92.1	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz): ⁶ Li	88.3	
Minimum Rejection (dB)	80	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.82: W1346680 - Filter HQ 600 ²H- P ⁶Li- S

1.6.3 Bandpass Filter




Figure 1.31: Bandpass Filter HQ 600 MHz Diagram

W1346241 - Filter HQ 600 ¹ H BP (⁻³¹ P, ¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H	600.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹⁹ F	< 242.9	564.7
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R069 [▶ 217] (mm)	227 x 92 x 39	
Weight (kg)	0.9	

Table 1.83: W1346241 - Filter HQ 600 ¹H BP (⁻³¹P, ¹⁹F)

W1346228 - Filter HQ 600 ¹ H BP (¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H	600.1	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹⁹ F	< 242.9	564.7
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R055P [▶ 219] (mm)	282 x 83 x 54	
Weight (kg)	1.6	

Table 1.84: W1346228 - Filter HQ 600 ¹H BP (¹⁹F)

W1346252 - Filter HQ 600 ¹⁹ F BP (⁻³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁹ F	564.7	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹ H	< 242.9	600.1
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R069 [▶ 217] (mm)	227 x 92 x 39	
Weight (kg)	0.9	

Table 1.85: W1346252 - Filter HQ 600 ¹⁹F BP (⁻³¹P, ¹H)

W1346229 - Filter HQ 600 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁹ F	564.7	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹ H	< 242.9	600.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R055P [▶ 219] (mm)	282 x 83 x 54	
Weight (kg)	1.6	

Table 1.86: W1346229 - Filter HQ 600 ¹⁹F BP (¹H)

W1346375 - Filter HQ 600 ¹⁹ F - ¹ H BP (³¹ P)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹⁹ F, ¹ H	564.7	600.1
Maximum Insertion Loss (dB)	0.8	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ... ³¹ P	< 242.9	
Minimum Rejection (dB)	110	
Mechanical Dimensions and Weight		
Case 5R069 [▶ 217] (mm)	227 x 92 x 39	
Weight (kg)	0.9	

Table 1.87: W1346375 - Filter HQ 600 ¹⁹F - ¹H BP (³¹P)

W1346279 - Filter HQ 600 ⁷ Li - ³¹ P BP (- ⁸⁷ Rb, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁷ Li... ³¹ P	233.2...242.9	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ⁸⁷ Rb, ¹⁹ F...	< 196.4	> 564.7
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.88: W1346279 - Filter HQ 600 ⁷Li - ³¹P BP (-⁸⁷Rb, ¹⁹F-)

W1346367 - Filter HQ 600 ⁷ Li - ³¹ P BP (- ¹³ C, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁷ Li... ³¹ P	233.2...242.9	
Maximum Insertion Loss (dB)	0.5	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ¹³ C, ¹⁹ F...	< 150.9	> 564.7
Minimum Rejection (dB)	90	80
Mechanical Dimensions and Weight		
Case 6H073 [▶ 222] (mm)	277 x 100 x 38	
Weight (kg)	1.6	

Table 1.89: W1346367 - Filter HQ 600 ⁷Li - ³¹P BP (-¹³C, ¹⁹F-)

W1346280 - Filter HQ 600 ⁷¹ Ga - ⁸⁷ Rb BP (- ¹³ C, ³¹ P-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ⁷¹ Ga... ⁸⁷ Rb	183.0...196.4	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ¹³ C, ³¹ P...	< 150.9	> 242.9
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.90: W1346280 - Filter HQ 600 ⁷¹Ga - ⁸⁷Rb BP (-¹³C, ³¹P-)

W1346490 - Filter HQ 600 ¹³ C - ²³ Na BP (- ²⁹ Si, ⁷¹ Ga)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) ¹³ C... ²³ Na	150.9...158.7	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ²⁹ Si, ⁷¹ Ga...	< 119.2	> 183.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.91: W1346490 - Filter HQ 600 ¹³C - ²³Na BP (-²⁹Si, ⁷¹Ga)

W1346276 - Filter HQ 600 ⁵⁹ Co - ²³ Na BP (- ²⁹ Si, ⁸⁷ Rb-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) ⁵⁹ Co... ²³ Na	142.4...158.7	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ²⁹ Si, ⁸⁷ Rb...	< 119.2	> 196.4
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H097 [▶ 223] (mm)	336 x 125 x 38	
Weight (kg)	2.5	

Table 1.92: W1346276 - Filter HQ 600 ⁵⁹Co - ²³Na BP (-²⁹Si, ⁸⁷Rb-)

W1346461 - Filter HQ 600 ⁵⁹ Co - ²³ Na BP (- ² H, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) ⁵⁹ Co... ²³ Na	142.4...158.7	
Maximum Insertion Loss (dB)	0.6	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ² H, ¹⁹ F...	< 92.1	> 564.7
Minimum Rejection (dB)	90	80
Mechanical Dimensions and Weight		
Case 6H073 [▶ 222] (mm)	277 x 100 x 38	
Weight (kg)	1.6	

Table 1.93: W1346461 - Filter HQ 600 ⁵⁹Co - ²³Na BP (-²H, ¹⁹F-)

W1346230 - Filter HQ 600 ⁵⁹ Co - ²³ Na BP (- ¹⁵ N)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) ⁵⁹ Co... ²³ Na	142.4...158.7	
Maximum Insertion Loss (dB)	0.5	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ¹⁵ N	< 60.8	
Minimum Rejection (dB)	80	
Mechanical Dimensions and Weight		
Case 4H073 [▶ 221] (mm)	205 x 100 x 38	
Weight (kg)	1.1	

Table 1.94: W1346230 - Filter HQ 600 ⁵⁹Co - ²³Na BP (-¹⁵N)

W1346281 - Filter HQ 600 ²⁹ Si BP (- ² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) ²⁹ Si	119.2	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ² H, ¹³ C...	< 92.1	> 150.9
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.95: W1346281 - Filter HQ 600 ²⁹Si BP (-²H, ¹³C-)

W1346454 - Filter HQ 600 ¹⁷O - ²H BP (¹¹B-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁷ O... ² H	81.3...92.1	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹¹ B...	>192.5	
Minimum Rejection (dB)	85	
Mechanical Dimensions and Weight		
Case 6H073 [▶ 222] (mm)	277 x 100 x 38	
Weight (kg)	1.6	

Table 1.96: W1346454 - Filter HQ 600 ¹⁷O - ²H BP (¹¹B-)

W1346282 - Filter HQ 600 ¹³³Cs - ²H BP (⁻¹⁵N, ²⁹Si-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹³³ Cs... ² H	78.7...92.1	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ¹⁵ N, ²⁹ Si.....	< 60.8	> 119.2
Minimum Rejection (dB)	90	75
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.97: W1346282 - Filter HQ 600 ¹³³Cs - ²H BP (⁻¹⁵N, ²⁹Si-)

W1346277 - Filter HQ 600 ¹⁵ N BP (- ¹⁴ N, ² H-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁵ N	60.8	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ... ¹⁴ N, ² H....	< 43.4	> 92.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.98: W1346277 - Filter HQ 600 ¹⁵N BP (-¹⁴N,²H-)

W1346283 - Filter HQ 600 ¹⁴ N BP (¹⁵ N-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁴ N	43.4	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁵ N....	> 60.8	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.99: W1346283 - Filter HQ 600 ¹⁴N BP (¹⁵N-)

W1346757 - Filter HQ 600 ²⁹ Si- ²⁰⁷ PB BP (- ² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ²⁹ Si.... ²⁰⁷ PB	119.2.....125.1	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz): ² H.... ¹³ C	< 92.1	> 150.9
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H073 [▶ 223] (mm)	336 x 197 x 38	
Weight (kg)	2.0	

Table 1.100: W1346757 - Filter HQ 600 ²⁹Si- ²⁰⁷PB BP (-²H, ¹³C-)

1.6.4 Diplexer Filter




Figure 1.32: Diplexer Filter HQ 600 MHz Diagram

W1346347 - Diplexer HQ 600 ¹ H / ¹⁹ F		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H, ¹⁹ F	600.1	564.7
Maximum Insertion Loss (dB)	0.6	0.6
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹⁹ F, ¹ H..	564.7	600.1
Minimum Rejection (dB)	70	70
Mechanical Dimensions and Weight		
Case 4S063D [▶ 215] (mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.101: W1346347 - Diplexer HQ 600 ¹H / ¹⁹F

1.7 AVANCE 700MHz

1.7.1 Lowpass Filter




Figure 1.33: Lowpass Filter HQ 700MHz Diagram

W1346639 - Filter HQ 700 0- ³¹ P LP (¹⁹ F- ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P	< 283.4	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹⁹ F... ¹ H..	658.8...700.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D050 [▶ 226] (mm)	144 x 63 x 34	
Weight (kg)	0.4	

Table 1.102: W1346639 - Filter HQ 700 0-³¹P LP (¹⁹F-¹H)

W1346659 - Filter HQ 700 0- ¹⁵ N NR LP (¹³ C)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁵ N	< 71.0	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ¹³ C	176.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.103: W1346659 - Filter HQ 700 0-¹⁵N NR LP (¹³C)

1.7.2 Bandstop Filter




Figure 1.34: Bandstop Filter HQ 700MHz Diagram

W1346284 - Filter HQ 700 ¹ H S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P, ¹⁹ F	< 283.4	658.8
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) :. ¹ H	700.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.104: W1346284 - Filter HQ 700 ¹H S

W1346387 - Filter HQ 700 ¹⁹ F S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P, ¹ H	< 283.4	700.1
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) :. ¹⁹ F	658.8	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.105: W1346387 - Filter HQ 700 ¹⁹F S

W1346394 - Filter HQ 700 ² H S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz):		
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)		
Frequency Stop (MHz) :. ² H	107.5	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.106: W1346394 - Filter HQ 700 ²H S

1.7.3 Bandpass Filter




Figure 1.35: Bandpass Filter HQ 700MHz Diagram

W1346296 - Filter HQ 700 ¹ H BP (- ³¹ P, ¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H	700.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :... ³¹ P, ¹⁹ F	< 283.4	658.8
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R045 [▶ 216] (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.107: W1346296 - Filter HQ 700 ¹H BP (-³¹P, ¹⁹F)

W1346288 - Filter HQ 700 ¹ H BP (¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H	700.1	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹⁹ F..	< 283.4	658.8
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R055P [▶ 219] (mm)	282 x 83 x 54	
Weight (kg)	1.6	

Table 1.108: W1346288 - Filter HQ 700 ¹H BP (¹⁹F)

W1346300 - Filter HQ 700 ¹⁹ F BP (- ³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 2% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F	658.8	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :.. ³¹ P, ¹ H	< 283.4	700.1
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R045 [▶ 216] (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.109: W1346300 - Filter HQ 700 ¹⁹F BP (-³¹P, ¹H)

W1346292 - Filter HQ 700 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F	658.8	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :... ³¹ P, ¹ H..	< 283.4	700.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R055P ▶ 219 (mm)	282 x 83 x 54	
Weight (kg)	1.6	

Table 1.110: W1346292 - Filter HQ 700 ¹⁹F BP (¹H)

W1346304 - Filter HQ 700 ⁷ Li - ³¹ P BP (⁻⁸⁷ Rb, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷ Li... ³¹ P	272.1...283.4	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. . ⁸⁷ Rb, ¹⁹ F....	< 229.1	> 658.8
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 6R125 ▶ 218 (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.111: W1346304 - Filter HQ 700 ⁷Li - ³¹P BP (⁻⁸⁷Rb, ¹⁹F-)

W1346308 - Filter HQ 700 ⁷¹ Ga - ⁸⁷ Rb BP (- ¹³ C, ⁷ Li-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷¹ Ga... ⁸⁷ Rb	213.5...229.1	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ¹³ C, ⁷ Li.....	< 176.0	> 272.1
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.112: W1346308 - Filter HQ 700 ⁷¹Ga - ⁸⁷Rb BP (-¹³C, ⁷Li-)

W1346491 - Filter HQ 700 ¹³ C - ²³ Na BP (- ²⁹ Si, ⁷¹ Ga-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹³ C... ²³ Na	176.0...185.2	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ²⁹ Si, ⁷¹ Ga...	< 139.1	> 213.5
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.113: W1346491 - Filter HQ 700 ¹³C - ²³Na BP (-²⁹Si, ⁷¹Ga-)

W1346312 - Filter HQ 700 ⁵⁹ Co - ²³ Na BP (- ²⁹ Si, ⁸⁷ Rb-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁵⁹ Co... ²³ Na	166.1...185.2	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ²⁹ Si, ⁸⁷ Rb...	< 139.1	> 229.1
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 8H073 [▶ 223] (mm)	336 x 100 x 38	
Weight (kg)	2.0	

Table 1.114: W1346312 - Filter HQ 700 ⁵⁹Co - ²³Na BP (-²⁹Si, ⁸⁷Rb-)

W1346462 - Filter HQ 700 ⁵⁹ Co - ²³ Na BP (- ² H, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁵⁹ Co... ²³ Na	166.1...185.2	
Maximum Insertion Loss (dB)	0.6	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : .. ² H, ¹⁹ F...	< 107.5	> 658.8
Minimum Rejection (dB)	90	80
Mechanical Dimensions and Weight		
Case 6H073 [▶ 222] (mm)	277 x 100 x 38	
Weight (kg)	1.6	

Table 1.115: W1346462 - Filter HQ 700 ⁵⁹Co - ²³Na BP (-²H, ¹⁹F-)

W1346316 - Filter HQ 700 ⁵⁹ Co - ²³ Na BP (- ¹⁵ N)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁵⁹ Co... ²³ Na	166.1...185.2	
Maximum Insertion Loss (dB)	0.5	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ... ¹⁵ N	< 71.0	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 4H073 [▶ 221] (mm)	205 x 100 x 38	
Weight (kg)	1.1	

Table 1.116: W1346316 - Filter HQ 700 ⁵⁹Co - ²³Na BP (-¹⁵N)

W1346320 - Filter HQ 700 ⁷⁷ Se - ²⁹ Si BP (- ² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷⁷ Se ... ²⁹ Si	133.5....139.1	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :... .. ² H, ¹³ C.....	< 107.5	> 176.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.117: W1346320 - Filter HQ 700 ⁷⁷Se -²⁹Si BP (-²H, ¹³C-)

W1346455 - Filter HQ 700 ¹⁷ O - ² H BP (¹¹ B-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁷ O... ² H	94.9...107.5	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹¹ B, ³¹ P.....	> 224.5	
Minimum Rejection (dB)	75	85
Mechanical Dimensions and Weight		
Case 6H073 [▶ 222] (mm)	277 x 100 x 38	
Weight (kg)	1.6	

Table 1.118: W1346455 - Filter HQ 700 ¹⁷O - ²H BP (¹¹B-)

W1346492 - Filter HQ 700 ¹⁷ O - ² H BP (- ¹⁵ N, ²⁹ Si-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁷ O... ² H	94.9...107.5	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. .. ¹⁵ N, ²⁹ Si....	< 71.0	> 139.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.119: W1346492 - Filter HQ 700 ¹⁷O - ²H BP (-¹⁵N, ²⁹Si-)

W1346324 - Filter HQ 700 ¹³³ Cs - ² H BP (- ¹⁵ N, ²⁹ Si-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹³³ Cs... ² H	91.8...107.5	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ¹⁵ N, ²⁹ Si....	< 71.0	> 139.1
Minimum Rejection (dB)	90	80
Mechanical Dimensions and Weight		
Case 8H097 [▶ 223] (mm)	336 x 125 x 38	
Weight (kg)	2.5	

Table 1.120: W1346324 - Filter HQ 700 ¹³³Cs - ²H BP (-¹⁵N, ²⁹Si-)

W1346328 - Filter HQ 700 ¹⁵ N BP (- ¹⁴ N, ² H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁵ N	71.0	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	12	
Frequency Stop (MHz) :. ... ¹⁴ N, ² H.....	< 50.6	> 107.5
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.121: W1346328 - Filter HQ 700 ¹⁵N BP (-¹⁴N, ²H)

W1346336 - Filter HQ 700 ¹⁴ N BP (¹⁵ N-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁴ N	50.6	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁵ N	> 71.0	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.122: W1346336 - Filter HQ 700 ¹⁴N BP (¹⁵N-)

1.7.4 Diplexer Filter




Figure 1.36: Diplexer Filter HQ 700MHz Diagram

W1346421 - Diplexer HQ 700 ¹ H / ¹⁹ F		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H, ¹⁹ F	700.1	658.80
Maximum Insertion Loss (dB)	0.6	0.6
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹⁹ F, ¹ H..	658.8	700.1
Minimum Rejection (dB)	60	60
Mechanical Dimensions and Weight		
Case 4S063D [▶ 215] mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.123: W1346421 - Diplexer HQ 700 ¹H / ¹⁹F

1.8 AVANCE 750MHz

1.8.1 Lowpass Filter




Figure 1.37: Lowpass Filter HQ 750MHz Diagram

W1346644 - Filter HQ 750 0- ³¹ P LP (¹⁹ F- ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P	< 303.7	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁹ F... ¹ H..	705.8...750.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D050 [▶ 226] (mm)	144 x 63 x 34	
Weight (kg)	0.4	

Table 1.124: W1346644 - Filter HQ 750 0-³¹P LP (¹⁹F-¹H)

W1346660 - Filter HQ 750 0- ¹⁵ N NR LP (¹³ C)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁵ N	< 76.0	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹³ C..	188.6	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.125: W1346660 - Filter HQ 750 0-¹⁵N NR LP (¹³C)

1.8.2 Bandstop Filter




Figure 1.38: Bandstop Filter HQ 750MHz Diagram

W1346285 Filter HQ 750 ¹ H S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P, ¹⁹ F	< 303.7	705.8
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) :. ¹ H	750.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.126: W1346285 Filter HQ 750 ¹H S

W1346385 - Filter HQ 750 ¹⁹ F S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P, ¹ H	< 303.7	750.1
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) :. ¹⁹ F...	705.8	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.127: W1346385 - Filter HQ 750 ¹⁹F S

1.8.3 Bandpass Filter




Figure 1.39: Bandpass Filter HQ 750MHz Diagram

W1346297 - Filter HQ 750 ¹ H BP (- ³¹ P, ¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 2% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H	750.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹⁹ F....	< 303.7	705.8
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R045 [▶ 216] (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.128: W1346297 - Filter HQ 750 ¹H BP (-³¹P, ¹⁹F)

W1346289 - Filter HQ 750 ¹ H BP (¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H	750.1	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : .. ³¹ P, ¹⁹ F...	< 303.7	705.8
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R055P [▶ 219] (mm)	282 x 83 x 54	
Weight (kg)	1.6	

Table 1.129: W1346289 - Filter HQ 750 ¹H BP (¹⁹F)

W1346301 - Filter HQ 750 ¹⁹ F BP (⁻³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F	705.8	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹ H....	< 303.7	750.1
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R045 [▶ 216] (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.130: W1346301 - Filter HQ 750 ¹⁹F BP (⁻³¹P, ¹H)

W1346293 - Filter HQ 750 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F	705.8	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ... ³¹ P, ¹ H...	< 303.7	750.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R055P ▶ 219 (mm)	282 x 83 x 54	
Weight (kg)	1.6	

Table 1.131: W1346293 - Filter HQ 750 ¹⁹F BP (¹H)

W1346305 - Filter HQ 750 ⁷ Li - ³¹ P BP (- ⁸⁷ Rb, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷ Li... ³¹ P	291.5...303.7	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. . ⁸⁷ Rb, ¹⁹ F....	< 245.4	> 705.8
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 6R125 ▶ 218 (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.132: W1346305 - Filter HQ 750 ⁷Li - ³¹P BP (-⁸⁷Rb, ¹⁹F-)

W1346309 - Filter HQ 750 ⁷¹ Ga - ⁸⁷ Rb BP (- ¹³ C, ⁷ Li-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷¹ Ga... ⁸⁷ Rb	228.8...245.4	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ¹³ C, ⁷ Li...	< 188.6	> 291.5
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.133: W1346309 - Filter HQ 750 ⁷¹Ga - ⁸⁷Rb BP (-¹³C, ⁷Li-)

W1346313 - Filter HQ 750 ¹³ C - ²³ Na BP (- ²⁹ Si, ⁷¹ Ga-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹³ C... ²³ Na	188.6...198.4	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ²⁹ Si, ⁷¹ Ga..	< 149.0	> 228.8
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.134: W1346313 - Filter HQ 750 ¹³C - ²³Na BP (-²⁹Si, ⁷¹Ga-)

W1346317 - Filter HQ 750 ⁵⁹ Co - ²³ Na BP (- ¹⁵ N)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁵⁹ Co... ²³ Na	178.0...198.4	
Maximum Insertion Loss (dB)	0.5	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ¹⁵ N	< 76.0	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 4H073 [▶ 221] (mm)	205 x 100 x 38	
Weight (kg)	1.1	

Table 1.135: W1346317 - Filter HQ 750 ⁵⁹Co - ²³Na BP (-¹⁵N)

W1346321 - Filter HQ 750 ⁷⁷ Se - ²⁹ Si BP (- ² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ⁷⁷ Se, ²⁹ Si..	143.1...149.0	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ² H, ¹³ C.....	< 115.2	> 188.6
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.136: W1346321 - Filter HQ 750 ⁷⁷Se -²⁹Si BP (-²H, ¹³C-)

W1346493 - Filter HQ 750 ¹⁷ O - ² H BP (- ¹⁵ N, ²⁹ Si-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁷ O... ² H	101.7...115.2	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : . ¹⁵ N, ²⁹ Si..	< 76.0	> 149.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.137: W1346493 - Filter HQ 750 ¹⁷O - ²H BP (-¹⁵N, ²⁹Si-)

W1346325 - Filter HQ 750 ¹³³ Cs - ² H BP (- ¹⁵ N, ²⁹ Si-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹³³ Cs... ² H	98.4...115.2	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : . ¹⁵ N, ²⁹ Si...	< 76.0	> 149.0
Minimum Rejection (dB)	80	80
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.138: W1346325 - Filter HQ 750 ¹³³Cs - ²H BP (-¹⁵N, ²⁹Si-)

W1346329 - Filter HQ 750 ¹⁵N BP (-¹⁴N, ²H-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁵ N	76.0	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ¹⁴ N, ² H...	< 54.2	> 115.2
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H169 [▶ 225] (mm)	336 x 197 x 38	
Weight (kg)	3.9	

Table 1.139: W1346329 - Filter HQ 750 ¹⁵N BP (-¹⁴N, ²H-)

1.8.4 Diplexer Filter




Figure 1.40: Diplexer Filter HQ 750MHz Diagram

W1346378 - Diplexer HQ 750 ¹ H / ¹⁹ F		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H, ¹⁹ F	750.1	705.8
Maximum Insertion Loss (dB)	0.6	0.6
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) ∴ ¹⁹ F, ¹ H...	705.8	750.1
Minimum Rejection (dB)	60	60
Mechanical Dimensions and Weight		
Case 4S063D [▶ 215] (mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.140: W1346378 - Diplexer HQ 750 ¹H / ¹⁹F

1.9 AVANCE 800MHz

1.9.1 Lowpass Filter




Figure 1.41: Lowpass Filter HQ 800MHz Diagram

W1346640 - Filter HQ 800 0- ³¹ P LP (¹⁹ F- ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P	< 323.9	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁹ F... ¹ H...	752.9...800.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D050 [▶ 226] (mm)	144 x 63 x 34	
Weight (kg)	0.4	

Table 1.141: W1346640 - Filter HQ 800 0-³¹P LP (¹⁹F-¹H)

W1346661 - Filter HQ 800 0- ¹⁵ N NR LP (¹³ C)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁵ N	< 81.1	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹³ C	201.2	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.142: W1346661 - Filter HQ 800 0-¹⁵N NR LP (¹³C)

1.9.2 Bandstop Filter




Figure 1.42: Bandstop Filter HQ 800MHz Diagram

W1346286 - Filter HQ 800 ¹ H S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P, ¹⁹ F	< 323.9	752.9
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹ H	800.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.143: W1346286 - Filter HQ 800 ¹H S

W1346386 - Filter HQ 800 ¹⁹ F S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P, ¹ H	< 323.9	800.1
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) :. ¹⁹ F	752.9	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.144: W1346386 - Filter HQ 800 ¹⁹F S

1.9.3 Bandpass Filter




Figure 1.43: Bandpass Filter HQ 800MHz Diagram

W1346298 - Filter HQ 800 ¹ H BP (⁻³¹ P, ¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 2% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H	800.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :... ³¹ P, ¹⁹ F..	< 323.9	752.9
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R045 [▶ 216] (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.145: W1346298 - Filter HQ 800 ¹H BP (⁻³¹P, ¹⁹F)

W1346290 - Filter HQ 800 ¹ H BP (¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H	800.1	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹⁹ F....	< 323.9	752.9
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R041P [▶ 220] (mm)	282 x 69 x 54	
Weight (kg)	0.0	

Table 1.146: W1346290 - Filter HQ 800 ¹H BP (¹⁹F)

W1346302 - Filter HQ 800 ¹⁹ F BP (⁻³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 2% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F	752.9	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P, ¹ H..	< 323.9	800.1
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R045 [▶ 216] (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.147: W1346302 - Filter HQ 800 ¹⁹F BP (⁻³¹P, ¹H)

W1346294 - Filter HQ 800 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F	752.9	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : .. ³¹ P, ¹ H...	< 323.9	800.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R041P ▶ 220 (mm)	282 x 69 x 54	
Weight (kg)	0.0	

Table 1.148: W1346294 - Filter HQ 800 ¹⁹F BP (¹H)

W1346306 - Filter HQ 800 ⁷ Li - ³¹ P BP (⁻⁸⁷ Rb, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷ Li... ³¹ P	311.0...323.9	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ⁸⁷ Rb, ¹⁹ F.....	< 261.8	> 752.9
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 ▶ 218 (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.149: W1346306 - Filter HQ 800 ⁷Li - ³¹P BP (⁻⁸⁷Rb, ¹⁹F-)

W1346310 - Filter HQ 800 ⁷¹ Ga - ⁸⁷ Rb BP (- ¹³ C, ⁷ Li-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷¹ Ga... ⁸⁷ Rb	244.0...261.8	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ¹³ C, ⁷ Li...	< 201.2	> 311.0
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.150: W1346310 - Filter HQ 800 ⁷¹Ga - ⁸⁷Rb BP (-¹³C, ⁷Li-)

W1346314 - Filter HQ 800 ¹³ C - ²³ Na BP (- ²⁹ Si, ⁷¹ Ga-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹³ C... ²³ Na	201.2...211.7	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : . ²⁹ Si, ⁷¹ Ga....	< 159.0	> 244.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.151: W1346314 - Filter HQ 800 ¹³C - ²³Na BP (-²⁹Si, ⁷¹Ga-)

W1346318 - Filter HQ 800 ⁵⁹ Co - ²³ Na BP (- ¹⁵ N)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁵⁹ Co... ²³ Na	189.8...211.7	
Maximum Insertion Loss (dB)	0.5	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ¹⁵ N	< 81.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 4H073 [▶ 221] (mm)	205 x 100 x 38	
Weight (kg)	1.1	

Table 1.152: W1346318 - Filter HQ 800 ⁵⁹Co - ²³Na BP (-¹⁵N)

W1346322 - Filter HQ 800 ⁷⁷ Se - ²⁹ Si BP (- ² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz):... ⁷⁷ Se, ²⁹ Si...	152.6...159.0	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ... ² H, ¹³ C..	< 122.8	> 201.2
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.153: W1346322 - Filter HQ 800 ⁷⁷Se - ²⁹Si BP (-²H, ¹³C-)

W1346326 - Filter HQ 800 ¹⁷ O - ² H BP (- ¹⁵ N, ²⁹ Si-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁷ O... ² H	108.5...122.8	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. . ¹⁵ N, ²⁹ Si...	< 81.1	> 159.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.154: W1346326 - Filter HQ 800 ¹⁷O - ²H BP (-¹⁵N, ²⁹Si-)

W1346330 - Filter HQ 800 ¹⁵ N BP (- ¹⁴ N, ² H-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁵ N	81.1	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. . ¹⁴ N, ² H....	< 57.8	> 122.8
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H169 [▶ 225] (mm)	336 x 197 x 38	
Weight (kg)	3.9	

Table 1.155: W1346330 - Filter HQ 800 ¹⁵N BP (-¹⁴N, ²H-)

1.10 AVANCE 850MHz

1.10.1 Lowpass Filter




Figure 1.44: Lowpass Filter HQ 850MHz Diagram

W1346677 - Filter HQ 850 0- ³¹ P LP (¹⁹ F- ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P	< 344.1	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁹ F... ¹ H...	799.9...850.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D050 [▶ 226] (mm)	144 x 63 x 34	
Weight (kg)	0.4	

Table 1.156: W1346677 - Filter HQ 850 0-³¹P LP (¹⁹F-¹H)

W1346656 - Filter HQ 850 0- ¹⁵ N NR LP (¹³ C)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁵ N	< 86.2	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹³ C	213.8	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.157: W1346656 - Filter HQ 850 0-¹⁵N NR LP (¹³C)

W128330 - Filter HQ 850 0- ³³ S NR LP (²⁵ Mg, H3)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³³ S	< 65.3	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ¹³ C	213.8	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.158: W128330 - Filter HQ 850 0-³³S NR LP (²⁵Mg, H3)

W128329 - Filter HQ 850 0- ⁸⁹ Y NR LP (¹⁸³ W, H3)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ⁸⁹ Y	< 41.7	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹³ C	213.8	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.159: W128329 - Filter HQ 850 0-⁸⁹Y NR LP (¹⁸³W, H3)

W128331 - Filter HQ 850 0- ¹⁵ N NR LP (⁹¹ Zr, H3)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁵ N	< 86.2	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹³ C	213.8	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.160: W128331 - Filter HQ 850 0-¹⁵N NR LP (⁹¹Zr, H3)

W131538 - Filter HQ 850 0- ³⁷ Cl NR LP (³³ S, H3)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³⁷ Cl	69.3	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹³ C	213.8	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.161: W131538 - Filter HQ 850 0-³⁷Cl NR LP (³³S, H3)

1.10.2 Bandpass Filter




Figure 1.45: Bandpass Filter HQ 850MHz Diagram

W119950 Filter HQ 850 ¹⁴ N BP		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁴ N	61.4	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁴ N	> 86.2	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 8H139 [▶ 224] (mm)	336 x 167 x 38	
Weight (kg)	3.4	

Table 1.162: W119950 Filter HQ 850 ¹⁴N BP

W1346655 - Filter HQ 850 ¹ H BP (⁻³¹ P, ¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 2% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H	850.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ³¹ P, ¹⁹ F....	< 344.1	799.9
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R045 [▶ 216] (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.163: W1346655 - Filter HQ 850 ¹H BP (⁻³¹P, ¹⁹F)

W1346674- Filter HQ 850 ¹ H BP (¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹ H	850.1	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ³¹ P, ¹⁹ F...	< 344.1	799.9
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R041P [▶ 220] (mm)	282 x 69 x 54	
Weight (kg)	0.0	

Table 1.164: W1346674- Filter HQ 850 ¹H BP (¹⁹F)

W1346675 - Filter HQ 850 ¹⁹ F BP (⁻³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 2% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁹ F	799.9	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ³¹ P, ¹ H...	< 344.1	850.1
Minimum Rejection (dB)	120	50
Mechanical Dimensions and Weight		
Case 5R045 [▶ 216] (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.165: W1346675 - Filter HQ 850 ¹⁹F BP (⁻³¹P, ¹H)

W1346676 - Filter HQ 850 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): .. ¹⁹ F	799.9	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ³¹ P, ¹ H	< 344.1	850.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R041P [▶ 220] (mm)	282 x 69 x 54	
Weight (kg)	0.0	

Table 1.166: W1346676 - Filter HQ 850 ¹⁹F BP (¹H)

W1346715 - Filter HQ 850 ⁷ Li - ³¹ P BP (- ⁸⁷ Rb, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ⁷ Li, ³¹ P	330.4...344.1	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ⁸⁷ R, ¹⁹ F...	< 278.2	> 799.9
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.167: W1346715 - Filter HQ 850 ⁷Li -³¹P BP (-⁸⁷Rb, ¹⁹F-)

W1346714 - Filter HQ 850 ⁷¹ Ga- ⁸⁷ Rb BP (- ¹³ C, ⁷ Li-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ⁷¹ Ga, ⁸⁷ Rb	259.3...278.2	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ¹³ C, ⁷ Li...	< 213.8	> 330.4
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.168: W1346714 - Filter HQ 850 ⁷¹Ga-⁸⁷Rb BP (-¹³C, ⁷Li-)

W1346710 - Filter HQ 850 ¹³C-²³Na BP (-²⁹Si, ⁷¹Ga-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹³ C, ²³ Na	213.8...224.8	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ²⁹ Si, ⁷¹ Ga...	< 168.9	> 259.3
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.169: W1346710 - Filter HQ 850 ¹³C-²³Na BP (-²⁹Si, ⁷¹Ga-)

W1346657 - Filter HQ 850 ⁵⁹Co-²³Na BP (-¹⁵N)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ⁵⁹ Co, ²³ Na	201.7...224.9	
Maximum Insertion Loss (dB)	0.5	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁵ N...	< 86.2	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 4H073 [▶ 221] (mm)	205 x 100 x 38	
Weight (kg)	1.1	

Table 1.170: W1346657 - Filter HQ 850 ⁵⁹Co-²³Na BP (-¹⁵N)

W1346718 - Filter HQ 850 ⁷⁷ Se - ²⁹ Si BP (- ² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ⁷⁷ Se, ²⁹ Si...	162.1...168.9	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ² H, ¹³ C...	< 130.5	> 213.8
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.171: W1346718 - Filter HQ 850 ⁷⁷Se - ²⁹Si BP (-²H, ¹³C-)

W1346697 - Filter HQ 850 ¹⁷ O- ² H BP (- ¹⁵ N, ²⁹ Si-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁷ O, ² H	115.2...130.5	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ¹⁵ N, ²⁹ Si...	< 86.2	> 168.9
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.172: W1346697 - Filter HQ 850 ¹⁷O-²H BP (-¹⁵N, ²⁹Si-)

W1346716- Filter HQ 850 ¹⁵ N BP (- ¹⁴ N, ² H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁵ N	82.6	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁴ N, ² H...	< 61.4	> 130.5
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H169 [▶ 225] (mm)	336 x 197 x 38	
Weight (kg)	3.9	

Table 1.173: W1346716- Filter HQ 850 ¹⁵N BP (-¹⁴N, ²H)

1.10.3 Diplexer Filter




Figure 1.46: Diplexer Filter HQ 850MHz Diagram

W116514 - Diplexer HQ 850 ¹ H / ¹⁹ F		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): .. ¹ H, ¹⁹ F...	850.1	799.9
Maximum Insertion Loss (dB)	0.6	0.6
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) :. ¹⁹ F, ¹ H...	799.9	850.1
Minimum Rejection (dB)	50	50
Mechanical Dimensions and Weight		
Case 4S063D [▶ 215] (mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.174: W116514 - Diplexer HQ 850 ¹H / ¹⁹F

1.11 AVANCE 900MHz

1.11.1 Lowpass Filter




Figure 1.47: Lowpass Filter HQ 900MHz Diagram

W1346641 - Filter HQ 900 0- ³¹ P LP (¹⁹ F- ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P	< 364.4	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁹ F... ¹ H...	847.0...900.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D050 [▶ 226] (mm)	144 x 63 x 34	
Weight (kg)	0.4	

Table 1.175: W1346641 - Filter HQ 900 0-³¹P LP (¹⁹F-¹H)

W1346662 - Filter HQ 900 0- ¹⁵ N NR LP (¹³ C)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁵ N	< 91.2	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹³ C	226.3	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.176: W1346662 - Filter HQ 900 0-¹⁵N NR LP (¹³C)

1.11.2 Bandstop Filter




Figure 1.48: Bandstop Filter HQ 900MHz Diagram

W1346287 - Filter HQ 900 ¹ H S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P, ¹⁹ F	< 364.4	847.0
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) :. ¹ H	900.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.177: W1346287 - Filter HQ 900 ¹H S

1.11.3 Bandpass Filter




Figure 1.49: Bandpass Filter HQ 900MHz Diagram

W1346299 - Filter HQ 900 ¹ H BP (⁻³¹ P, ¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 2% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H	900.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ... ³¹ P, ¹⁹ F...	< 364.4	847.0
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R045 [216] (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.178: W1346299 - Filter HQ 900 ¹H BP (⁻³¹P, ¹⁹F)

W1346291 - Filter HQ 900 ¹ H BP (¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H	900.1	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :: ... ³¹ P, ¹⁹ F...	< 364.4	847.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R041P ▶ 220 (mm)	282 x 69 x 54	
Weight (kg)	0.0	

Table 1.179: W1346291 - Filter HQ 900 ¹H BP (¹⁹F)

W1346303 - Filter HQ 900 ¹⁹ F BP (⁻³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 2% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F	847.0	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :: ... ³¹ P, ¹ H...	< 364.4	900.1
Minimum Rejection (dB)	120	50
Mechanical Dimensions and Weight		
Case 5R045 ▶ 216 (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.180: W1346303 - Filter HQ 900 ¹⁹F BP (⁻³¹P, ¹H)

W1346295 - Filter HQ 900 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F	847.0	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ... ³¹ P, ¹ H...	< 364.4	900.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8R041P [▶ 220] (mm)	282 x 69 x 54	
Weight (kg)	0.0	

Table 1.181: W1346295 - Filter HQ 900 ¹⁹F BP (¹H)

W1346307 - Filter HQ 900 ⁷ Li - ³¹ P BP (⁻⁸⁷ Rb, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷ Li... ³¹ P	349.8...364.4	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :... ⁸⁷ Rb, ¹⁹ F..	< 294.5	> 847.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.182: W1346307 - Filter HQ 900 ⁷Li - ³¹P BP (⁻⁸⁷Rb, ¹⁹F-)

W1346311 - Filter HQ 900 ⁷¹ Ga - ⁸⁷ Rb BP (⁻¹³ C, ⁷ Li-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷¹ Ga... ⁸⁷ Rb	274.5...294.5	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ... ¹³ C, ⁷ Li.....	< 226.3	> 349.8
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.183: W1346311 - Filter HQ 900 ⁷¹Ga - ⁸⁷Rb BP (⁻¹³C, ⁷Li-)

W1346315 - Filter HQ 900 ¹³ C - ²³ Na BP (⁻²⁹ Si, ⁷¹ Ga-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹³ C... ²³ Na	226.3...238.1	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. . ²⁹ Si, ⁷¹ Ga..	< 178.8	> 274.5
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.184: W1346315 - Filter HQ 900 ¹³C - ²³Na BP (⁻²⁹Si, ⁷¹Ga-)

W1346319 - Filter HQ 900 ⁵⁹ Co - ²³ Na BP (- ¹⁵ N)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁵⁹ Co... ²³ Na	213.6...238.1	
Maximum Insertion Loss (dB)	0.5	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁵ N	< 91.2	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 4H073 [▶ 221] (mm)	205 x 100 x 38	
Weight (kg)	1.1	

Table 1.185: W1346319 - Filter HQ 900 ⁵⁹Co - ²³Na BP (-¹⁵N)

W1346323 - Filter HQ 900 ⁷⁷ Se - ²⁹ Si BP (- ² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷⁷ Se,... ²⁹ Si	171.1...178.8	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ² H, ¹³ C...	< 138.2	> 226.3
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.186: W1346323 - Filter HQ 900 ⁷⁷Se - ²⁹Si BP (-²H, ¹³C-)

W1346458 - Filter HQ 900 ¹⁷ O - ² H BP (¹¹ B-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁷ O... ² H	122.0...138.2	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ¹¹ B, ³¹ P...	>288.8	
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.187: W1346458 - Filter HQ 900 ¹⁷O - ²H BP (¹¹B-)

W1346327 - Filter HQ 900 ¹⁷ O - ² H BP (¹⁵ N, ²⁹ Si-)		
RF Power Specification		
Maximum Power Operating	1kW @100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁷ O... ² H	122.0...138.2	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ¹⁵ N, ²⁹ Si...	< 91.2	> 178.8
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.188: W1346327 - Filter HQ 900 ¹⁷O - ²H BP (¹⁵N, ²⁹Si-)

W1346331 - Filter HQ 900 ¹⁵ N BP (⁻¹⁴ N, ² H-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁵ N	91.2	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : .. ¹⁴ N, ² H..	< 65.0	> 118.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 8H169 [▶ 225] (mm)	336 x 197 x 38	
Weight (kg)	3.9	

Table 1.189: W1346331 - Filter HQ 900 ¹⁵N BP (⁻¹⁴N, ²H-)

W1346335 - Filter HQ 900 ¹⁵ N BP (¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁵ N	91.2	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : .. ¹³ C	213.6	> 226.3
Minimum Rejection (dB)	85	90
Mechanical Dimensions and Weight		
Case 5H169 (mm)	241 x 197 x 38	
Weight (kg)	2.5	

Table 1.190: W1346335 - Filter HQ 900 ¹⁵N BP (¹³C-)

1.11.4 Diplexer Filter




Figure 1.50: Diplexer Filter HQ 900MHz Diagram

W1346423 - Diplexer HQ 900 ¹ H / ¹⁹ F		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H,... ¹⁹ F	900.1	847.0
Maximum Insertion Loss (dB)	0.6	0.6
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹⁹ F,... ¹ H	847.0	900.1
Minimum Rejection (dB)	50	50
Mechanical Dimensions and Weight		
Case 4S063D [▶ 215] (mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.191: W1346423 - Diplexer HQ 900 ¹H / ¹⁹F

1.12 AVANCE 950MHz

1.12.1 Lowpass Filter




Figure 1.51: Lowpass Filter HQ 950MHz Diagram

W143786 - Filter HQ 950 0- ³¹ P LP (¹⁹ F- ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P	< 384.6	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁹ F... ¹ H...	894.0...950.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D050 [▶ 226] (mm)	144 x 63 x 34	
Weight (kg)	0.4	

Table 1.192: W143786 - Filter HQ 950 0-³¹P LP (¹⁹F-¹H)

W143787 - Filter HQ 950 0- ¹⁵ N NR LP (¹³ C)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁵ N	< 96.3	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹³ C	238.9	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.193: W143787 - Filter HQ 950 0-¹⁵N NR LP (¹³C)

1.12.2 Bandstop Filter




Figure 1.52: Bandstop Filter HQ 950MHz Diagram

W143773 - Filter HQ 950 ¹ H S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P, ¹⁹ F	< 384.6	894.0
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) :. ¹ H	950.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.194: W143773 - Filter HQ 950 ¹H S

1.12.3 Bandpass Filter




Figure 1.53: Bandpass Filter HQ 950MHz Diagram

W1346707 - Filter HQ 950 ¹ H BP (- ³¹ P, ¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 2% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H	950.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ... ³¹ P, ¹⁹ F...	< 384.6	894.0
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R045 [▶ 216] (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.195: W1346707 - Filter HQ 950 ¹H BP (-³¹P, ¹⁹F)

W143775 - Filter HQ 950 ¹⁹ F BP (- ³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 2% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F	894.0	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ... ³¹ P, ¹ H...	< 384.6	950.1
Minimum Rejection (dB)	110	45
Mechanical Dimensions and Weight		
Case 5R045 [▶ 216] (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.196: W143775 - Filter HQ 950 ¹⁹F BP (-³¹P, ¹H)

W143774 - Filter HQ 950 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F	894.0	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ... ³¹ P, ¹ H...	< 384.6	950.1
Minimum Rejection (dB)	80	90
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.197: W143774 - Filter HQ 950 ¹⁹F BP (¹H)

W143776 - Filter HQ 950 ⁷ Li - ³¹ P BP (- ⁸⁷ Rb, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷ Li... ³¹ P	369.2...384.6	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :... ⁸⁷ Rb, ¹⁹ F..	< 310.9	> 894.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.198: W143776 - Filter HQ 950 ⁷Li - ³¹P BP (-⁸⁷Rb, ¹⁹F-)

W140205 - Filter HQ 950 ⁷¹ Ga - ⁸⁷ Rb BP (- ¹³ C, ⁷ Li-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷¹ Ga... ⁸⁷ Rb	289.8...310.9	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ... ¹³ C, ⁷ Li...	< 238.9	> 369.3
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.199: W140205 - Filter HQ 950 ⁷¹Ga - ⁸⁷Rb BP (-¹³C, ⁷Li-)

W140204 - Filter HQ 950 ¹³ C - ²³ Na BP (⁻²⁹ Si, ⁷¹ Ga-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹³ C... ²³ Na	238.9...251.3	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ⁻²⁹ Si, ⁷¹ Ga..	< 178.8	> 289.8
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.200: W140204 - Filter HQ 950 ¹³C - ²³Na BP (⁻²⁹Si, ⁷¹Ga-)

W143783 - Filter HQ 950 ⁵⁹ Co - ²³ Na BP (⁻² H, ³¹ P-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁵⁹ Co... ²³ Na	225.4...251.3	
Maximum Insertion Loss (dB)	0.6	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ⁻² H, ³¹ P...	< 145.9	> 384.6
Minimum Rejection (dB)	90	80
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.201: W143783 - Filter HQ 950 ⁵⁹Co - ²³Na BP (⁻²H, ³¹P-)

W139785 - Filter HQ 950 ¹⁷ O - ² H BP (- ¹⁵ N, ²⁹ Si-)		
RF Power Specification		
Maximum Power Operating	1kW @100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁷ O... ² H	128.8...145.9	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ¹⁵ N, ²⁹ Si..	< 96.3	> 188.8
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.202: W139785- Filter HQ 950 ¹⁷O - ²H BP (-¹⁵N, ²⁹Si-)

W143778 - Filter HQ 950 ¹⁵ N BP (- ¹⁴ N, ² H-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁵ N	96.3	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : .. ¹⁴ N, ² H..	< 68.7	> 145.9
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.203: W143778 - Filter HQ 950 ¹⁵N BP (-¹⁴N, ²H-)

W143779 - Filter HQ 950 ¹⁵ N BP (¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁵ N	96.3	
Maximum Insertion Loss (dB)	0.8	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ⁵⁹ Co, ¹³ C..	225.4	> 238.9
Minimum Rejection (dB)	85	90
Mechanical Dimensions and Weight		
Case 5H169 (mm)	241 x 197 x 38	
Weight (kg)	2.5	

Table 1.204: W143779 - Filter HQ 950 ¹⁵N BP (¹³C-)

W143780 - Filter HQ 950 ¹⁴ N BP (¹⁵ N-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁴ N	68.7	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : ¹⁵ N	> 96.3	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case *(mm)	L x W x H	
Weight (kg)	0.0	

Table 1.205: W143780 - Filter HQ 950 ¹⁴N BP (¹⁵N-)

W1346869 - Filter HQ 950 ¹ H BP (¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H	950.1	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ³¹ P, ¹⁹ F...	<384.6	894.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case *(mm)	L x W x H	
Weight (kg)	0.0	

Table 1.206: W1346869 - Filter HQ 950 ¹H BP (¹⁹F)

W140203- Filter HQ 950 ⁷⁷ Se - ²⁹ Si BP (² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) ⁷⁷ Se... ²⁹ Si	181.2...188.8	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)		
Frequency Stop (MHz) ² H... ¹³ C	< 145.9	> 238.9
Minimum Rejection (dB)		
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.207: W140203- Filter HQ 950 ⁷⁷Se -²⁹Si BP (²H, ¹³C-)

1.12.4 Diplexer Filter




Figure 1.54: Diplexer Filter HQ 950MHz Diagram

W134781 - Diplexer HQ 950 ¹ H / ¹⁹ F		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H, ... ¹⁹ F	950.1	894.0
Maximum Insertion Loss (dB)	0.6	0.6
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹⁹ F, ¹ H	894.0	950.1
Minimum Rejection (dB)	50	50
Mechanical Dimensions and Weight		
Case 4S063D [▶ 215] (mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.208: W134781 - Diplexer HQ 950 ¹H / ¹⁹F

1.13 AVANCE 1000MHz

1.13.1 Lowpass Filter




Figure 1.55: Lowpass Filter HQ 1000MHz Diagram

W1346678- Filter HQ 1000 0- ³¹ P LP (¹⁹ F- ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P	< 404.9	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁹ F... ¹ H...	941.1...1000.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D050 [▶ 226] (mm)	144 x 63 x 34	
Weight (kg)	0.4	

Table 1.209: W1346678- Filter HQ 1000 0-³¹P LP (¹⁹F-¹H)

W1346672- Filter HQ 1000 0- ¹⁵ N NR LP (¹³ C)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁵ N	< 101.4	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹³ C.	251.5	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.210: W1346672- Filter HQ 1000 0-¹⁵N NR LP (¹³C)

1.13.2 Bandstop Filter




Figure 1.56: Bandstop Filter HQ 1000MHz Diagram

W1346663 - Filter HQ 1000 ¹ H S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P, ¹⁹ F	< 404.9	941.1
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) :. ¹ H...	1000.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.211: W1346663 - Filter HQ 1000 ¹H S

W1346664 - Filter HQ 1000 ¹⁹ F S		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P, ¹ H	< 404.9	1000.1
Maximum Insertion Loss (dB)	0.4	0.8
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) :. ¹⁹ F...	941.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 3S063 [▶ 215] (mm)	161 x 106 x 50	
Weight (kg)	1.1	

Table 1.212: W1346664 - Filter HQ 1000 ¹⁹F S

1.13.3 Bandpass Filter




Figure 1.57: Bandpass Filter HQ 1000MHz Diagram

W1346665 - Filter HQ 1000 ¹ H BP (⁻³¹ P, ¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 2% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H	1000.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ³¹ P, ¹⁹ F....	< 404.9	941.1
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R045 [▶ 216] (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.213: W1346665 - Filter HQ 1000 ¹H BP (⁻³¹P, ¹⁹F)

W1346666 - Filter HQ 1000 ¹⁹ F BP (- ³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F	941.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁹ F	< 404.9	1000.1
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R045 [▶ 216] (mm)	227 x 68 x 39	
Weight (kg)	0.7	

Table 1.214: W1346666 - Filter HQ 1000 ¹⁹F BP (-³¹P, ¹H)

W1346667 - Filter HQ 1000 ⁷ Li - ³¹ P BP (- ⁸⁷ Rb, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷ Li... ³¹ P	388.7...404.9	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ⁸⁷ Rb, ¹⁹ F	< 327.2	> 941.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.215: W1346667 - Filter HQ 1000 ⁷Li -³¹P BP (-⁸⁷Rb, ¹⁹F-)

W1346668- Filter HQ 1000 ⁷¹Ga- ⁸⁷Rb BP (-¹³C, ⁷Li-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁷¹ Ga... ⁸⁷ Rb	305.0...327.2	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ¹³ C, ⁷ Li...	< 251.5	> 388.7
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.216: W1346668- Filter HQ 1000 ⁷¹Ga- ⁸⁷Rb BP (-¹³C, ⁷Li-)

W1346669 - Filter HQ 1000 ¹³C-²³Na BP (-²⁹Si, ⁷¹Ga-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹³ C... ²³ Na	251.5...264.5	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ²⁹ Si, ⁷¹ Ga...	< 198.7	> 305.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.217: W1346669 - Filter HQ 1000 ¹³C-²³Na BP (-²⁹Si, ⁷¹Ga-)

W1346670 - Filter HQ 1000 ²⁹ Si BP (- ² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): .. ²⁹ Si	198.7	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :.. ² H, ¹³ C	< 153.5	> 251.5
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	0.0	

Table 1.218: W1346670 - Filter HQ 1000 ²⁹Si BP (-²H, ¹³C-)

W1346671 - Filter HQ 1000 ¹⁷ O - ² H BP (- ¹⁵ N, ²⁹ Si-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁷ O, ² H	135.6...153.5	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) :. ¹⁵ N, ²⁹ Si...	< 101.4	> 198.7
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.219: W1346671 - Filter HQ 1000 ¹⁷O -²H BP (-¹⁵N, ²⁹Si-)

W1346717 - Filter HQ 1000 ¹⁵ N BP (- ¹⁴ N, ² H-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁵ N	101.4	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁴ N, ² H...	< 72.3	> 153.5
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.220: W1346717 - Filter HQ 1000 ¹⁵N BP (-¹⁴N, ²H-)

W1346856 - Filter HQ 1000 ⁵⁹ Co - ²³ Na BP (- ² H, ³¹ P-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ⁵⁹ Co... ²³ Na	273.3...264.5	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)		
Frequency Stop (MHz) :. ² H, ³¹ P...	< 153.5	> 404.9
Minimum Rejection (dB)		
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.221: W1346856 - Filter HQ 1000 ⁵⁹Co -²³Na BP (-²H, ³¹P-)

W1346870 - Filter HQ 1000 ¹ H BP (¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H...	1000.1	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ³¹ P, ¹⁹ F...	< 404.9	941.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.222: W1346870 - Filter HQ 1000 ¹H BP (¹⁹F)

W1346871 - Filter HQ 1000 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F...	941.0	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ³¹ P, ¹ H...	< 404.9	1000.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.223: W1346871 - Filter HQ 1000 ¹⁹F BP (¹H)

1.13.4 Diplexer Filter




Figure 1.58: Diplexer Filter HQ 1000MHz Diagram

W1346876 - Diplexer HQ 1000 ¹ H / ¹⁹ F		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H,... ¹⁹ F	1000.1	941.1
Maximum Insertion Loss (dB)	0.6	0.6
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹⁹ F, ¹ H...	941.1	1000.1
Minimum Rejection (dB)	50	50
Mechanical Dimensions and Weight		
Case 4S063D [▶ 215] (mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.224: W1346876 - Diplexer HQ 1000 ¹H / ¹⁹F

1.14 AVANCE 1100MHz

1.14.1 Lowpass Filter




Figure 1.59: Lowpass Filter HQ 1100MHz Diagram

W1342624 - Filter HQ 1100 0- ¹⁵ N NR LP (¹³ C)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁵ N	<111.5	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹³ C	276.6	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.225: W1342624 - Filter HQ 1100 0-¹⁵N NR LP (¹³C)

W1342373 - Filter HQ 1100 0- ³¹ P LP (¹⁹ F- ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P	<445.3	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁹ F,... ¹ H	1035.1...1100.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D050 [▶ 226] (mm)	144 x 63 x 34	
Weight (kg)	0.4	

Table 1.226: W1342373 - Filter HQ 1100 0-³¹P LP (¹⁹F-¹H)

1.14.2 Bandpass Filter




Figure 1.60: Bandpass Filter HQ 1100MHz Diagram

W1342132- Filter HQ 1100 ¹ H BP (- ³¹ P, ¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ¹ H	1100.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : : . ³¹ P, ¹⁹ F	< 445.3	1035.1
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R032 [▶ 216] (mm)	227 x 55 x 39	
Weight (kg)	0.7	

Table 1.227: W1342132- Filter HQ 1100 ¹H BP (-³¹P, ¹⁹F)

W1346872- Filter HQ 1100 ¹ H BP (¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ¹ H	1100.1	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : . : ³¹ P, ¹⁹ F...	< 445.3	1035.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.228: W1346872- Filter HQ 1100 ¹H BP (¹⁹F)

W1346859 - Filter HQ 1100 ⁷ Li - ³¹ P BP (- ⁸⁷ Rb, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ⁷ Li... ³¹ P	427.6...445.3	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : . : ⁸⁷ Rb, ¹⁹ F...	< 360.0	> 1035.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.229: W1346859 - Filter HQ 1100 ⁷Li-³¹P BP (-⁸⁷Rb, ¹⁹F-)

W1346861- Filter HQ 1100 ⁷¹ Ga – ⁸⁷ Rb BP (⁻¹³ C, ⁷ Li-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ⁷¹ Ga... ⁸⁷ Rb	335.5...360.0	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : : . ¹³ C, ⁷ Li...	< 276.6	> 427.5
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.230: W1346861- Filter HQ 1100 ⁷¹Ga –⁸⁷Rb BP (⁻¹³C, ⁷Li-)

W1346863 - Filter HQ 1100 ⁷⁷ Se – ²⁹ Si BP (⁻² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ⁷⁷ Se... ²⁹ Si	209.8...218.6	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)		
Frequency Stop (MHz) : : . ² H, ¹³ C...	< 168.9	> 276.6
Minimum Rejection (dB)		
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.231: W1346863 - Filter HQ 1100 ⁷⁷Se –²⁹Si BP (⁻²H, ¹³C-)

W1346865 - Filter HQ 1100 ¹⁷ O – ² H BP (- ¹⁵ N, ²⁹ Si)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ¹⁷ O... ² H	149.1...168.9	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : . : ¹⁵ N, ²⁹ Si...	< 111.5	> 218.6
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.232: W1346865 - Filter HQ 1100 ¹⁷O-²H BP (-¹⁵N, ²⁹Si)

W1346867- Filter HQ 1100 ¹⁵ N BP (- ¹⁴ N, ² H-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ¹⁵ N	111.5	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : . : ¹⁴ N, ² H...	< 79.5	> 168.9
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.233: W1346867- Filter HQ 1100 ¹⁵N BP (-¹⁴N, ²H-)

W1346873- Filter HQ 1100 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ¹⁹ F	1035.1	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : : . ³¹ P, ¹ H...	< 445.3	1100.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.234: W1346873- Filter HQ 1100 ¹⁹F BP (¹H)

W1348165- Filter HQ 1100 ¹³ C- ²³ Na BP (- ²⁹ Si, ⁷¹ Ga-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ¹³ C... ²³ Na	276.6...291.0	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : : . ²⁹ Si,... ⁷¹ Ga	< 218.6	> 335.5
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.235: W1348165- Filter HQ 1100 ¹³C-²³Na BP (-²⁹Si, ⁷¹Ga-)

W1346473- Filter HQ 1100 ¹⁹ F BP (- ³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 2% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ¹⁹ F	1035.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : . ³¹ P, . ¹ H	< 445.3	1100.1
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R032 [▶ 216] (mm)	227 x 55 x 39	
Weight (kg)	0.0	

Table 1.236: W1346473- Filter HQ 1100 ¹⁹F BP (-³¹P, ¹H)

W1346857 - Filter HQ 1100 ⁵⁹ Co – ²³ Na BP (- ² H, ³¹ P-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ⁵⁹ Co... ²³ Na	261.0...291.0	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)		
Frequency Stop (MHz) : . . ² H, ³¹ P...	< 168.9	> 445.3
Minimum Rejection (dB)		
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.237: W1346857 - Filter HQ 1100 ⁵⁹Co-²³Na BP (-²H, ³¹P-)

1.14.3 Diplexer Filter




Figure 1.61: Diplexer Filter HQ 1100MHz Diagram

W1346877 - Diplexer HQ 1100 ¹ H / ¹⁹ F		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹ H, ... ¹⁹ F	1100.1	1035.1
Maximum Insertion Loss (dB)	0.6	0.6
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹⁹ F, ¹ H	1035.1	1100.1
Minimum Rejection (dB)	50	50
Mechanical Dimensions and Weight		
Case 4S063D [▶ 215] (mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.238: W1346877- Diplexer HQ 1100 ¹H / ¹⁹F

1.15 AVANCE 1200MHz

1.15.1 Lowpass Filter




Figure 1.62: Lowpass Filter HQ 1200MHz Diagram

W134926 - Filter HQ 1200 0- ¹⁵ N NR LP (¹³ C)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ¹⁵ N	<121.7	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹³ C	301.8	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D095 [▶ 226] (mm)	240 x 55 x 37	
Weight (kg)	0.7	

Table 1.239: W134926 - Filter HQ 1200 0-¹⁵N NR LP (¹³C)

W134925 - Filter HQ 1200 0- ³¹ P LP (¹⁹ F- ¹ H)		
RF Power Specification		
Maximum Power Operating	1.5kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ... ³¹ P	<485.9	
Maximum Insertion Loss (dB)	0.2	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) :. ¹⁹ F, ¹ H..	1129.2...1200.1	
Minimum Rejection (dB)	90	
Mechanical Dimensions and Weight		
Case 2D050 [▶ 226] (mm)	144 x 63 x 34	
Weight (kg)	0.4	

Table 1.240: W134925 - Filter HQ 1200 0-³¹P LP (¹⁹F-¹H)

1.15.2 Bandpass Filter




Figure 1.63: Bandpass Filter HQ 1200MHz Diagram

W1348174 Filter HQ 1200 ¹ H BP (- ³¹ P, ¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H	1200.1	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ... ³¹ P,.. ¹⁹ F	< 485.8	1129.2
Minimum Rejection (dB)	120	45
Mechanical Dimensions and Weight		
Case 5R032 [▶ 216] (mm)	227 x 55 x 39	
Weight (kg)	0.7	

Table 1.241: W1348174 Filter HQ 1200 ¹H BP (-³¹P, ¹⁹F)

W1346874- Filter HQ 1200 ¹ H BP (¹⁹ F)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ¹ H	1200.1	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ³¹ P, ¹⁹ F...	< 485.8	1129.2
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.242: W1346874- Filter HQ 1200 ¹H BP (¹⁹F)

W1346860 - Filter HQ 1200 ⁷ Li - ³¹ P BP (- ⁸⁷ Rb, ¹⁹ F-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ⁷ Li... ³¹ P	466.4...485.8	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : : ⁸⁷ Rb, ¹⁹ F...	< 392.7	> 1129.2
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.243: W1346860 - Filter HQ 1200 ⁷Li-³¹P BP (-⁸⁷Rb, ¹⁹F-)

W1346862- Filter HQ 1200 ⁷¹ Ga – ⁸⁷ Rb BP (- ¹³ C, ⁷ Li-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ⁷¹ Ga... ⁸⁷ Rb	366.0...392.7	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : . : ¹³ C, ⁷ Li...	< 301.8	> 466.4
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.244: W1346862- Filter HQ 1200 ⁷¹Ga-⁸⁷Rb BP (-¹³C, ⁷Li-)

W1346864 - Filter HQ 1200 ⁷⁷ Se – ²⁹ Si BP (- ² H, ¹³ C-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ⁷⁷ Se... ²⁹ Si	228.9...283.4	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)		
Frequency Stop (MHz) : . : ² H, ¹³ C...	< 184.2	> 301.8
Minimum Rejection (dB)		
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.245: W1346864 - Filter HQ 1200 ⁷⁷Se-²⁹Si BP (-²H, ¹³C-)

W1346866 - Filter HQ 1200 ¹⁷ O – ² H BP (⁻¹⁵ N, ²⁹ Si)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ¹⁷ O... ² H	162.7...184.2	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : : . ¹⁵ N, ²⁹ Si...	< 121.7	> 238.4
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.246: W1346866 - Filter HQ 1200 ¹⁷O-²H BP (⁻¹⁵N, ²⁹Si)

W1346868- Filter HQ 1200 ¹⁵ N BP (⁻¹⁴ N, ² H-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ¹⁵ N	121.7	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : : . ¹⁴ N, ² H...	< 86.7	> 184.2
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 7R155 [▶ 218] (mm)	336 x 194 x 52	
Weight (kg)	3.4	

Table 1.247: W1346868- Filter HQ 1200 ¹⁵N BP (⁻¹⁴N, ²H-)

W1346879- Filter HQ 1200 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ¹⁹ F	1129.2	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)	20	
Frequency Stop (MHz) : ³¹ P, ¹ H	< 485.8	1200.1
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.248: W1346879- Filter HQ 1200 ¹⁹F BP (¹H)

W135028- Filter HQ 1200 ¹³ C- ²³ Na BP (⁻²⁹ Si, ⁷¹ Ga-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ¹³ C... ²³ Na	301.8...317.5	
Maximum Insertion Loss (dB)	1	
Minimum Return Loss (dB)	15	
Frequency Stop (MHz) : . ²⁹ Si, ⁷¹ Ga	< 238.4	> 366.0
Minimum Rejection (dB)	90	90
Mechanical Dimensions and Weight		
Case 6R125 [▶ 218] (mm)	250 x 153 x 39	
Weight (kg)	1.5	

Table 1.249: W135028- Filter HQ 1200 ¹³C-²³Na BP (⁻²⁹Si, ⁷¹Ga-)

W135027 Filter HQ 1200 ¹⁹ F BP (- ³¹ P, ¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F	1129.2	
Maximum Insertion Loss (dB)	0.7	
Minimum Return Loss (dB)	20	
Frequency Stop (MHz): ... ³¹ P, ¹ H	< 485.8	1200.1
Minimum Rejection (dB)		
Mechanical Dimensions and Weight		
Case 5R032 [▶ 216] (mm)	227 x 55 x 39	
Weight (kg)	0.7	

Table 1.250: W135027 Filter HQ 1200 ¹⁹F BP (-³¹P, ¹H)

W1346879 Filter HQ 1200 ¹⁹ F BP (¹ H)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): ¹⁹ F	1129.2	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)		
Frequency Stop (MHz)	< 485.8	1200.1
Minimum Rejection (dB): ... ³¹ P, .. ¹ H		
Mechanical Dimensions and Weight		
Case 8R041P [▶ 220] (mm)	282 x 69 x 54	
Weight (kg)	0.0	

Table 1.251: W1346879 Filter HQ 1200 ¹⁹F BP (¹H)

W1346858 - Filter HQ 1200 ⁵⁹ Co – ²³ Na BP (² H, ³¹ P-)		
RF Power Specification		
Maximum Power Operating	1kW @ 100 ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz): . ⁵⁹ Co... ²³ Na	284.8...317.5	
Maximum Insertion Loss (dB)		
Minimum Return Loss (dB)		
Frequency Stop (MHz) : . ² H, ³¹ P...	< 184.2	> 485.8
Minimum Rejection (dB)		
Mechanical Dimensions and Weight		
Case * (mm)	L x W x H	
Weight (kg)	0.0	

Table 1.252: W1346858 - Filter HQ 1200 ⁵⁹Co-²³Na BP (²H, ³¹P-)

1.15.3 Diplexer Filter




Figure 1.64: Diplexer Filter HQ 1200MHz Diagram

W1346878 - Diplexer HQ 1200 ¹ H / ¹⁹ F		
RF Power Specification		
Maximum Power Operating	1kW @ 100ms 5% duty cycle	
Operating Temperature Range (°C)	15...50	
RF Low Level Specification		
Frequency Pass (MHz) : ¹ H,... ¹⁹ F	1200.1	1129.2
Maximum Insertion Loss (dB)	0.6	0.6
Minimum Return Loss (dB)	20	20
Frequency Stop (MHz) : ¹⁹ F, ¹ H	1129.2	1200.1
Minimum Rejection (dB)	50	50
Mechanical Dimensions and Weight		
Case 4S063D [▶ 215] (mm)	188 x 106 x 69	
Weight (kg)	1.5	

Table 1.253: W1346878- Diplexer HQ 1200 ¹H / ¹⁹F

1.16 Filters Category

1.16.1 Lowpass Filter

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346642	2	300	<121.5	282.4...300.1	0.2	20	90	2D050
W1346636	2	400	<162.0	376.5...400.1	0.2	20	90	2D050
W1346637	2	500	<202.5	470.6...500.1	0.2	20	90	2D050
W1346638	2	600	<242.9	564.7...600.1	0.2	20	90	2D050
W1346639	2	700	<283.4	658.8...700.1	0.2	20	90	2D050
W1346644	2	750	<303.7	705.8...750.1	0.2	20	90	2D050
W1346640	2	800	<323.9	752.9...800.1	0.2	20	90	2D050
W1346677	2	850	<344.1	799.9...850.1	0.2	20	90	2D050
W1346641	2	900	<364.4	847.0...900.1	0.2	20	90	2D050
W143786	2	950	<384.6	894.0...950.1	0.2	20	90	2D050
W1346678	2	1000	<404.9	941.1...1000.1	0.2	20	90	2D050
W1342373	2	1100	<445.3	1035.1...1100.1	0.2	20	90	2D050
W134925	2	1200	<485.9	1129.2...1200.1	0.2	20	90	2D050

Table 1.254: Filters HQ 0-31P LP (¹⁹F-¹H)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346628	2	400	<40.6	100.6	0.2	20	90	2D095
W1346629	2	500	<50.7	125.8	0.2	20	90	2D095
W1346630	2	600	<60.8	150.9	0.2	20	90	2D095
W1346659	2	700	<71.0	176.1	0.2	20	90	2D095
W1346660	2	750	<76.0	188.6	0.2	20	90	2D095
W1346661	2	800	<81.1	201.2	0.2	20	90	2D095
W1346656	2	850	<86.2	213.8	0.2	20	90	2D095
W1346662	2	900	<91.2	226.3	0.2	20	90	2D095
W143787	2	950	<96.3	238.9	0.2	20	90	2D095
W1346672	2	1000	<101.4	251.5	0.2	20	90	2D095
W1342624	2	1100	<111.5	276.6	0.2	20	90	2D095
W134926	2	1200	<121.7	301.8	0.2	20	90	2D095

Table 1.255: Filters HQ 0-¹⁵N NR LP (¹³C)

1.16.2 Bandstop Filter

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346253	3	300	<121.5 / 282.4	300.1	0.4 / 0.8	20 / 20	90	3S063
W1346232	3	500	<202.5 / 470.6	500.1	0.4 / 0.8	20 / 20	90	3S063
W1346255	3	600	<242.9 / 564.7	600.1	0.4 / 0.8	20 / 20	90	3S063
W1346284	3	700	<283.4 / 658.8	700.1	0.4 / 0.8	20 / 20	90	3S063
W1346285	3	750	<303.7 / 705.8	750.1	0.4 / 0.8	20 / 20	90	3S063
W1346286	3	800	<323.9 / 752.9	800.1	0.4 / 0.8	20 / 20	90	3S063
W1346287	3	900	<364.4 / 847.0	900.1	0.4 / 0.8	20 / 20	90	3S063
W143773	3	950	<384.6 / 894.0	950.1	0.4 / 0.8	20 / 20	90	3S063
W1346663	3	1000	<404.9 / 941.1	1000.1	0.4 / 0.8	20 / 20	90	3S063

Table 1.256: Filters HQ ¹H S

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346369	3	300	<121.5 / 300.1	282.4	0.4 / 0.8	20 / 20	90	3S063
W1346370	3	400	<162.0 / 400.1	376.5	0.4 / 0.4	15 / 20	90	3S063
W1346358	3	500	<202.5 / 500.1	470.6	0.4 / 0.8	20 / 20	90	3S063
W1346371	3	600	<242.9 / 600.1	564.7	0.4 / 0.8	20 / 20	90	3S063
W1346387	3	700	<283.4 / 700.1	658.8	0.4 / 0.8	20 / 20	90	3S063
W1346385	3	750	<303.7 / 750.1	705.8	0.4 / 0.8	20 / 20	90	3S063
W1346386	3	800	<323.9 / 800.1	752.9	0.4 / 0.8	20 / 20	90	3S063
W1346664	3	1000	<404.9 / 1000.1	941.1	0.4 / 0.8	20 / 20	90	3S063

Table 1.257: Filters HQ ¹⁹F S

NMR Filters

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346430	3	500	194.4	202.5	0.8	20	90	3S063

Table 1.258: Filters HQ $^7\text{Li P} - ^{31}\text{P S}$

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346431	3	500	202.5	194.4	0.8	20	90	3S063

Table 1.259: Filters HQ $^{31}\text{P P} - ^7\text{Li S}$

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346437	3	400	61.4	76.3	0.4	15	90	3S063

Table 1.260: Filters HQ $^2\text{H P} - ^{77}\text{Se S}$

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346394	3	700		107.5			90	3S063

Table 1.261: Filters HQ $^2\text{H S}$

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346434	3	500	73.6	76.8	0.8	20	80	3S063

Table 1.262: Filters HQ $^6\text{Li P} - ^2\text{H S}$

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346435	3	500	76.8	73.6	0.8	20	80	3S063

Table 1.263: Filters HQ ^2H P - ^6Li S

1.16.3 Bandpass Filters

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max.Insertion (dB)	Min.Return Loss (dB)	Min. Rejection (dB)	Case
W1346239	5	300	300.1	<121.5 / 282.4	0.7	20	120 / 45	5R100
W1346240	5	400	400.1	<162.0 / 376.5	0.7	20	120 / 45	5R100
W1346223	5	500	500.1	<202.5 / 470.6	0.7	20	120 / 45	5R069
W1346241	5	600	600.1	<242.9 / 564.7	0.7	20	120 / 45	5R069
W1346296	5	700	700.1	<283.4 / 658.8	0.7	20	120 / 45	5R045
W1346297	5	750	750.1	<303.7 / 705.8	0.7	20	120 / 45	5R045
W1346298	5	800	800.1	<323.9 / 752.9	0.7	20	120 / 45	5R045
W1346655	5	850	850.1	<344.1 / 799.9	0.7	20	120 / 45	5R045
W1346299	5	900	900.1	<364.4 / 847.0	0.7	20	120 / 45	5R045
W1346707	5	950	950.1	<384.6 / 894.0	0.7	20	120 / 45	5R045
W1346665	5	1000	1000.0	<404.9 / 941.1	0.7	20	120 / 45	5R045
W1342132	5	1100	1100.1	<445.3 / 1035.1	0.7	20	120 / 45	5R032
W1348174	5	1200	1200.1	<485.8 / 1129.2	0.7	20	120 / 45	5R032

Table 1.264: Filters HQ ¹H BP (-³¹P, ¹⁹F)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346235	8	300	300.1	<121.5 / 282.4	0.8	20	80 / 90	8R105P
W1346204	8	400	400.1	<162.0 / 376.5	0.8	20	80 / 90	8R105P
W1346247	8	500	500.1	<202.5 / 470.6	0.9	20	90 / 90	8R055P
W1346228	8	600	600.1	<242.9 / 564.7	0.8	20	90 / 90	8R055P
W1346288	8	700	700.1	<283.4 / 658.8	0.8	20	90 / 90	8R055P
W1346289	8	750	750.1	<303.7 / 705.8	0.8	20	90 / 90	8R055P
W1346290	8	800	800.1	<323.9 / 752.9	0.8	20	90 / 90	8R041P
W1346674	8	850	850.1	<344.1 / 799.9	0.8	20	90 / 90	8R041P
W1346291	8	900	900.1	<364.4 / 847.0	1	20	90 / 90	8R041P
W1346869	8	950	950.1	<384.6 / 894.0		20	90 / 90	
W1346870	8	1000	1000.1	<404.9 / 941.1		20	90 / 90	
W1346872	8	1100	1100.1	<445.3 / 1035.1		20	90 / 90	
W1346874	8	1200	1200.1	<485.8 / 1129.2		20	90 / 90	

Table 1.265: Filters HQ ¹H BP (¹⁹F)

NMR Filters

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346249	5	300	282.4	<121.5 / 300.1	0.7	20	120 / 45	5R100
W1346250	5	400	376.5	<162.0 / 400.1	0.7	20	110 / 45	5R100
W1346251	5	500	470.6	<202.5 / 500.1	0.7	20	120 / 45	5R069
W1346252	5	600	564.7	<242.9 / 600.1	0.7	20	120 / 45	5R069
W1346300	5	700	658.8	<283.4 / 700.1	0.7	20	120 / 45	5R045
W1346301	5	750	705.8	<303.7 / 750.1	0.7	20	120 / 45	5R045
W1346302	5	800	752.9	<323.9 / 800.1	0.7	20	120 / 45	5R045
W1346675	5	850	799.9	<344.1 / 850.1	0.7	20	120 / 50	5R045
W1346303	5	900	847.0	<364.4 / 900.1	0.7	20	120 / 50	5R045
W143775	5	950	894.0	<384.6 / 950.1	0.7	20	110 / 45	5R045
W1346666	5	1000	941.1	<404.9 / 1000.1	0.7	20	120 / 45	5R045
W1346473	5	1100	1035.1	<445.3 / 1100.1	0.7	20	120 / 45	5R032
W135027	5	1200	1129.2	<485.8 / 1200.1	0.7	20	120 / 45	5R032

Table 1.266: Filters HQ ¹⁹F BP (-³¹P, ¹H)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346236	8	300	282.4	<121.5 / 300.1	0.8	20	80 / 90	8R105P
W1346205	8	400	376.5	<162.0 / 400.1	0.8	20	80 / 90	8R105P
W1346248	8	500	470.6	<202.5 / 500.1	0.9	20	90 / 90	8R055P
W1346229	8	600	564.7	<242.9 / 600.1	0.8	20	90 / 90	8R055P
W1346292	8	700	658.8	<283.4 / 700.1	0.8	20	90 / 90	8R055P
W1346293	8	750	705.8	<303.7 / 750.1	0.8	20	90 / 90	8R055P
W1346294	8	800	752.9	<323.9 / 800.1	0.8	20	90 / 90	8R041P
W1346676	8	850	799.9	<344.1 / 850.1	0.8	20	90 / 90	8R041P
W1346295	8	900	847.0	<364.4 / 900.1	1	20	90 / 90	8R041P
W143774	8	950	894.0	<384.6 / 950.1		20	80 / 90	
W1346871	8	1000	941.1	<404.9 / 1000.1		20	90 / 90	
W1346873	8	1100	1035.1	<445.3 / 1100.1		20	90 / 90	
W1346879	8	1200	1129.2	<485.8 / 1200.1		20	90 / 90	

Table 1.267: Filters HQ ¹⁹F BP (1H)

NMR Filters

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346368	5	300	282.4 / 300.1	<121.5	0.8 / 0.8	20 / 20	110	5R100
W1346373	5	400	376.5 / 400.1	<162.0	0.8 / 0.8	20 / 20	110	5R100
W1346374	5	500	470.6 / 500.1	<202.5	0.8 / 0.8	20 / 20	110	5R069
W1346375	5	600	564.7 / 600.1	<242.9	0.8 / 0.8	20 / 20	110	5R069

Table 1.268: Filters HQ ¹⁹F - ¹H BP (-³¹P)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346256	6	400	155.5...162.0	<130.9 / >376.5	1	15	90 / 90	6R125
W1346262	6	500	194.4...202.5	<163.6 / >470.6	1	15	80 / 90	6R125
W1346279	6	600	233.2...242.9	<196.4 / >564.7	1	15	80 / 90	6R125
W1346304	6	700	272.1...283.4	<229.1 / >658.8	1	15	80 / 90	6R125
W1346305	6	750	291.5...303.7	<245.4 / >705.8	1	15	80 / 90	6R125
W1346306	6	800	311.0...323.9	<261.8 / >752.9	1	15	90 / 90	6R125
W1346715	6	850	330.4...344.1	<278.2 / >799.9	1	15	90 / 90	6R125
W1346307	6	900	349.8...364.4	<294.5 / >847.0	1	15	90 / 90	6R125
W143776	6	950	369.2...384.6	<310.9 / >894.0	1	15	90 / 90	6R125
W1346667	6	1000	388.7...404.9	<327.2 / >941.1	1	15	90 / 90	6R125
W1346859	6	1100	427.6...445.3	<360.0 / >1035.1	1	15	90 / 90	6R125
W1346860	6	1200	466.4...485.8	<392.7 / >1129.2	1	15	90 / 90	6R125

Table 1.269: Filters HQ ⁷Li - ³¹P BP (-⁸⁷Rb, ¹⁹F-)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346372	6	300	116.6...121.5	<75.5 / >282.4	0.5	20	90 / 80	6H073
W1346349	6	400	155.5...162.0	<100.6 / >376.5	0.5	20	90 / 80	6H073
W1346351	6	500	194.4...202.5	<125.8 / >470.6	0.5	20	80 / 70	6H073
W1346367	6	600	233.2...242.9	<150.9 / >564.7	0.5	20	90 / 80	6H073

Table 1.270: Filters HQ ⁷Li - ³¹P BP (-¹³C, ¹⁹F-)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346270	8	300	91.5...98.2	<75.5 / >121.5	1	15	90 / 85	8H097
W1346257	7	400	122.0...130.9	<100.6 / >162.0	1	15	90 / 90	7R155
W1346263	8	500	152.5...163.6	<125.8 / >202.5	1	15	90 / 90	8H097
W1346280	6	600	183.0...196.4	<150.9 / >242.9	1	15	80 / 90	6R125

Table 1.271: Filters HQ ⁷¹Ga - ⁸⁷Rb BP (-¹³C, ³¹P-)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346308	6	700	213.5...229.1	<176.0 / >272.1	1	15	80 / 90	6R125
W1346309	6	750	228.8...245.4	<188.6 / >291.5	1	15	80 / 90	6R125
W1346310	6	800	244.0...261.8	<201.2 / >311.0	1	15	80 / 90	6R125
W1346714	6	850	259.3...278.2	<213.8 / >330.4	1	15	90 / 90	6R125
W1346311	6	900	274.5...294.5	<226.3 / >349.8	1	15	90 / 90	6R125
W140205	6	950	289.8...310.9	<238.9 / >369.3	1	15	90 / 90	6R125
W1346668	6	1000	305.0...327.2	<251.5 / >388.7	1	15	90 / 90	6R125
W1346861	6	1100	355.5...360.0	<276.6 / >427.5	1	15	90 / 90	6R125
W1346862	6	1200	366.0...392.7	<301.8 / >466.4	1	15	90 / 90	6R125

Table 1.272: Filters HQ ^{71}Ga - ^{87}Rb BP (^{-13}C , ^7Li -)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346488	7	400	100.6...105.8	<79.5 / >122.0	1	15	90 / 90	7R155
W1346489	7	500	125.8...132.3	<99.4 / >152.5	1	15	90 / 90	7R155
W1346490	6	600	150.9...158.7	<119.2 / >183.0	1	15	90 / 90	6R125
W1346491	6	700	176.0...185.2	<139.1 / >213.5	1	15	90 / 90	6R125
W1346313	6	750	188.6...198.4	<149.0 / >228.8	1	15	90 / 90	6R125
W1346314	6	800	201.2...211.7	<159.0 / >244.0	1	15	90 / 90	6R125
W1346710	6	850	213.8...224.8	<168.9 / >259.3	1	15	90 / 90	6R125
W1346315	6	900	226.3...238.1	<178.8 / >274.5	1	15	90 / 90	6R125
W140204	6	950	238.9...251.3	<178.8 / >289.8	1	15	90 / 90	6R125
W1346669	6	1000	251.5...264.5	<198.7 / >305.0	1	15	90 / 90	6R125
W1348165	6	1100	276.6...291.0	<218.6 / >335.5	1	15	90 / 90	6R125
W135028	6	1200	301.8...317.5	<238.4 / >366.0	1	15	90 / 90	6R125

Table 1.273: Filters HQ ^{13}C - ^{23}Na BP (^{-29}Si , ^{71}Ga -)

NMR Filters

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346271	8	300	71.2...79.4	<59.6 / >80.0	1	15	80 / 80	8H139
W1346258	8	400	94.9...105.8	<79.5 / >130.9	1	15	80 / 80	8H139
W1346264	8	500	118.7...132.3	<99.4 / >163.6	1	15	90 / 90	8H097
W1346276	8	600	142.4...158.7	<119.2 / >196.4	1	15	90 / 90	8H097
W1346312	8	700	166.1...185.2	<139.1 / >229.1	1	15	80 / 90	8H073

Table 1.274: Filters HQ ⁵⁹Co - ²³Na BP (-²⁹Si, ⁸⁷Rb-)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346440	6	400	94.9...105.8	<61.4 / >376.5	0.6	20	90 / 90	6H073
W1346460	6	500	118.7...132.3	<76.8 / >470.6	0.6	20	90 / 80	6H073
W1346461	6	600	142.4...158.7	<92.1 / >564.7	0.6	20	90 / 80	6H073
W1346462	6	700	166.1...185.2	<107.5 / >658.8	0.6	20	90 / 80	6H073

Table 1.275: Filters HQ ⁵⁹Co - ²³Na BP (-²H, ¹⁹F-)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346275	4	300	71.2...79.4	<30.4	0.5	20	90	4H073
W1346226	4	400	94.9...105.8	<40.6	0.5	20	90	4H073
W1346224	4	500	118.7...132.3	<50.7	0.5	20	90	4H073
W1346230	4	600	142.4...158.7	<60.8	0.5	20	80	4H073
W1346316	4	700	166.1...185.2	<71.0	0.5	20	90	4H073
W1346317	4	750	178.0...198.4	<76.0	0.5	20	90	4H073
W1346318	4	800	189.8...211.7	<81.1	0.5	20	90	4H073
W1346657	4	850	201.7...224.9	<86.2	0.5	20	90	4H073
W1346319	4	900	213.6...238.1	<91.2	0.5	20	90	4H073

Table 1.276: Filters HQ ^{59}Co - ^{23}Na BP (^{-15}N)

NMR Filters

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346265	7	500	95.4...99.4	<76.8 / >125.8	1	15	90 / 90	7R155
W1346281	7	600	119.2	<92.1 / >150.9	1	15	90 / 90	7R155
W1346320	7	700	133.5...139.1	<107.5 / >176.0	1	15	90 / 90	7R155
W1346321	7	750	143.1...149.0	<115.2 / >188.6	1	15	90 / 90	7R155
W1346322	6	800	152.6...159.0	<122.8 / >201.2	1	15	90 / 90	6R125
W1346718	6	850	162.1...168.9	<130.5 / >213.8	1	15	90 / 90	6R125
W1346323	6	900	171.1...178.8	<138.2 / >226.3	1	15	90 / 90	6R125
W140203	6	950						
W1346670	6	1000	198.7	<153.5 / >251.5	1	15	90 / 90	6R125
W1346863	6	1100	209.8...218.6	<168.9 / >276.6				
W1346864	6	1200	228.9...283.4	<184.2 / >301.8				

Table 1.277: Filters HQ ⁷⁷Se-²⁹Si BP (-²H, ¹³C-)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346452	5	300	40.7...46.0	>96.3	0.8	20	75 / 80	5H169
W1346442	5	400	54.2...61.4	>128.3	0.8	20	75 / 80	5H169
W1346453	6	500	67.8...76.7	>160.4	0.8	20	75 / 85	5H169
W1346454	6	600	81.3...92.1	>192.5	0.8	20	85	6H073
W1346455	6	700	94.9...107.5	>224.5	0.8	20	75 / 85	6H073
W1346458	6	900	122.0...138.2	>288.8	0.8	15	90 / 90	7R155

Table 1.278: Filters HQ ¹⁷O - ²H BP (¹B-)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346492	7	700	94.9...107.5	<71.0 / >139.1	1	15	90 / 90	7R155
W1346493	7	750	101.7...115.2	<76.0 / >149.0	1	15	90 / 90	7R155
W1346326	7	800	108.5...122.8	<81.1 / >159.0	1	15	90 / 90	7R155
W1346697	7	850	115.2...130.5	<86.2 / >168.9	1	15	90 / 90	7R155
W1346327	7	900	122.0...138.2	<91.2 / >178.8	1	15	90 / 90	7R155
W139785	7	950	128.8...145.9	<96.3 / >188.8	1	15	90 / 90	7R155
W1346671	7	1000	135.6...153.5	<101.4 / >198.7	1	15	90 / 90	7R155
W1346865		1100	149.1...168.9	<111.5 / >218.6	1	15	90 / 90	6R125
W1346866		1200	162.7...184.2	<121.7 / >238.4	1	15	90 / 90	

Table 1.279: Filters HQ ¹⁷O - ²H BP (⁻¹⁵N, ²⁹Si-)

NMR Filters

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346273	8	300	39.4...46.0	<30.4 / >59.6	1	15	90 / 80	8H169
W1346260	8	400	52.5...61.4	<40.6 / >79.5	1	15	90 / 80	8H139
W1346266	8	500	65.6...76.8	<50.7 / >99.4	1	15	90 / 80	8H139
W1346282	8	600	78.7...92.1	<60.8 / >119.2	1	15	90 / 75	8H139
W1346324	8	700	91.8...107.5	<71.0 / >139.1	1	15	90 / 80	8H097
W1346325	7	750	98.4...115.2	<76.0 / >149.0	1	15	80 / 80	7R155

Table 1.280: Filters HQ ¹³³Cs - ²H BP (⁻¹⁵N, ²⁹Si-)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346274	8	300	30.4	<21.7 / >46.0	1	15	90 / 90	8H169
W1346261	8	400	40.6	<28.9 / >61.4	1	20	90 / 90	8H139
W1346267	8	500	50.7	<36.1 / >76.8	1	20	90 / 90	8H139
W1346277	8	600	60.8	<43.4 / >92.1	1	20	90 / 90	8H139
W1346328	8	700	71.0	<50.6 / >107.5	1	12	90 / 90	8H139
W1346329	8	750	76.0	<54.2 / >115.2	1	20	90 / 90	8H169
W1346330	8	800	81.1	<57.8 / >122.8	1	20	90 / 90	8H169
W1346716	8	850	86.2	<61.4 / >130.5	1	20	90 / 90	8H169
W1346331	8	900	91.2	<65.0 / >118.0	1	20	90 / 90	8H169
W143778	8	950	96.3	<68.7 / >145.9	1	20	90 / 90	*
W1346717	7	1000	101.4	<72.3 / >153.5	1	20	90 / 90	7R155
W1346867	7	1100	111.5	<79.5 / >168.9	1	20	90 / 90	7R155
W1346868	8	1200	121.7	<86.7 / >184.2	1	20	90 / 90	7R155

Table 1.281: Filters HQ ¹⁵N BP (⁻¹⁴N, ²H-)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346278	5	300	30.4	71.2 / >75.5	0.8	20	90 / 90	5H169
W1346335	5	900	91.2	213.6 / >226.3	0.8	20	85 / 90	5H169
W143779	5	950	96.3	225.4 / >238.9	0.8	20	85 / 90	5H169

Table 1.282: Filters HQ ¹⁵N BP (¹³C)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346283	8	600	43.4	>60.8	1	20	90	8H139
W1346336	8	700	50.6	>71.0	1	20	90	8H139
W143780	8	950	68.7	>96.3	1	15	90	8H139

Table 1.283: Filters HQ ¹⁴N BP (¹⁵N-)

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W143783		950	225.4...251.3	<145.9 / >384.6	0.6	20	90 / 80	
W1346856		1000	273.3...264.5	<153.5 / >404.9	0.6	20		
W1346857		1100	261.0...291.0	<168.9 / >445.3	0.6	20		
W1346858		1200	284.8...317.5	<184.2 / >485.8	0.6	20		

Table 1.284: Filters HQ ⁵⁹Co -²³Na BP (²H, ³¹P)

1.16.4 Diplexer Filter

Part Number	Number of Cells	Filter (MHz)	Frequency Pass (MHz)	Frequency Stop (MHz)	Max. Insertion (dB)	Min. Return Loss (dB)	Min. Rejection (dB)	Case
W1346237	4	300	300.1 / 282.4	282.4 / 300.1	0.6 / 0.6	20 / 20	50 / 50	4S063D
W1346345	4	400	400.1 / 376.5	376.5 / 400.1	0.6 / 0.6	20 / 20	70 / 70	4S063D
W1346346	4	500	500.1 / 470.6	470.6 / 500.1	0.6 / 0.6	20 / 20	70 / 70	4S063D
W1346347	4	600	600.1 / 564.7	564.7 / 600.1	0.6 / 0.6	20 / 20	70 / 70	4S063D
W1346378	4	750	750.1 / 705.8	705.8 / 750.1	0.6 / 0.6	20 / 20	60 / 60	4S063D
W116514	4	850	850.1 / 799.9	799.9 / 850.1	0.6 / 0.6	20 / 20	50 / 50	4S063D
W1346423	4	900	900.1 / 847.0	847.0 / 900.1	0.6 / 0.6	20 / 20	50 / 50	4S063D
W134781	4	950	950.1 / 894.0	894.0 / 950.1	0.6 / 0.6	20 / 20	50 / 50	4S063D
W1346876	4	1000	1000.1 / 941.1	941.1 / 1000.1	0.6 / 0.6	20 / 20	50 / 50	4S063D
W1346877	4	1100	1100.0 / 1035.1	1035.1 / 1100.1	0.6 / 0.6	20 / 20	50 / 50	4S063D
W1346878	4	1200	1200.1 / 1129.2	1129.2 / 1200.1	0.6 / 0.6	20 / 20	50 / 50	4S063D

Table 1.285: Diplexer HQ ¹H / ¹⁹F

2 Cases and Accessories

2.1 Notch Filters

2.1.1 Case 3S063

Mechanical dimensions (mm): 161 x 106 x 50

Weight (kg): 1,1



Figure 2.1: 3S063 View

2.1.2 Case 4S063D

Mechanical dimensions (mm): 242 x 107 x 68

Weight (kg): 1,5



Figure 2.2: 4S063D View

2.2 Bandpass Filters

2.2.1 Case 5R032

Mechanical dimensions (mm): 227 x 55 x 39

Weight (kg): 0,7




Figure 2.3: 5R032 View

2.2.2 Case 5R045

Mechanical dimensions (mm): 227 x 68 x 39

Weight (kg): 0,7



Figure 2.4: 5R045 View

2.2.3 Case 5R069

Mechanical dimensions (mm) : 227 x 92 x 39

Weight (kg): 0,9



Figure 2.5: 5R069 View

2.2.4 Case 5R100

Mechanical dimensions (mm) : 227 x 123 x 39

Weight (kg): 1,2



Figure 2.6: 5R100 View

2.2.5 Case 6R125

Mechanical dimensions (mm) : 250 x 153 x 39

Weight (kg): 1,5



Figure 2.7: 6R125 View

2.2.6 Case 7R155

Mechanical dimensions (mm) : 336 x 194 x 52

Weight (kg): 3,4



Figure 2.8: 7R155 View

2.2.7 Case 8R055P

Mechanical dimensions (mm): 282 x 83 x 54

Weight (kg): 1,6



Figure 2.9: 8R055P View

2.2.8 Case 8R105P

Mechanical dimensions (mm) : 277 x 133 x 54

Weight (kg) : 2,6



Figure 2.10: 8R105P View

2.2.9 Case 8R041P

Mechanical dimensions (mm) : 282 x 69 x 54

Weight (kg) : 0,0

image

8R041P View

2.3 Helical Filters

2.3.1 Case 4H073

Mechanical dimensions (mm): 205 x 100 x 38

Weight (kg): 1,1



Figure 2.11: 4H073 View

2.3.2 Case 5H133

Mechanical dimensions (mm): 241 x 157 x 38

Weight (kg): 1,9



Figure 2.12: 5H133 View

2.3.3 Case 5H069

Mechanical dimensions (mm): 241 x 197 x 38

Weight (kg): 2,5



Figure 2.13: 5H069 View

2.3.4 Case 6H073

Mechanical dimensions (mm) : 277 x 100 x 38

Weight (kg): 1,6



Figure 2.14: 6H073 View

2.3.5 Case 8H073

Mechanical dimensions (mm): 336 x 100 x 38

Weight (kg): 2.0



Figure 2.15: 8H073 View

2.3.6 Case 8H097

Mechanical dimensions (mm): 336 x 125 x 38

Weight (kg): 2.5



Figure 2.16: 8H097 View

2.3.7 Case 8H139

Mechanical dimensions (mm): 336 x 167 x 38

Weight (kg): 3,4



Figure 2.17: 8H139 View

2.3.8 Case 8H169

Mechanical dimensions (mm): 336 x 197 x 38

Weight (kg): 3,9



8H169 View

2.4 Low Pass Filters

2.4.1 Case 2D050

Mechanical dimensions (mm) : 144 x 58 x 34

Weight (kg) : 0,4



Figure 2.18: 2D050 View

2.4.2 Case 2D095

Mechanical dimensions (mm): 240 x 55 x 37

Weight (kg): 0,7



Figure 2.19: 2D095 View

2.5 Accessories



Figure 2.20: W1303715 - N / Female-Female 0,5m Coaxial Cable



Figure 2.21: 30649 - N / Female-Female Adapter



Figure 2.22: 4241 - N / Male-Male Adapter



Figure 2.23: 30209 - N / Male-Female Right Angle Adapter

3 Typical Curves

3.1 Bandpass Filters Curves




Figure 3.1: W1346239 - Filter HQ 300 ¹H BP (-³¹P, ¹⁹F)

Typical Curves




Figure 3.2: W1346249 - Filter HQ 300 ¹⁹F BP (-³¹P, ¹H)




Figure 3.3: W1346240 - Filter HQ 400 ¹H BP (-³¹P, ¹⁹F)




Figure 3.4: W1346250 - Filter HQ 400 ¹⁹F BP (-³¹P, ¹H)




Figure 3.5: W1346223 - Filter HQ 500 ¹H BP(-³¹P, ¹⁹F)

Typical Curves




Figure 3.6: W1346251 - Filter HQ 500 ¹⁹F BP (-³¹P, ¹H)




Figure 3.7: W1346241 - Filter HQ 600 ¹H BP (-³¹P, ¹⁹F)



Figure 3.8: W1346296 - Filter HQ 700 ¹H BP (⁻³¹P, ¹⁹F)




Figure 3.9: W1346300 - Filter HQ 700 ¹⁹F BP (⁻³¹P, ¹H)

Typical Curves



Figure 3.10: W1346297 - Filter HQ 750 ¹H BP (-³¹P, ¹⁹F)




Figure 3.11: W1346298 - Filter HQ 800 ¹H BP (-³¹P, ¹⁹F)




Figure 3.12: W1346299- Filter HQ 900 ¹H BP (-³¹P, ¹⁹F)




Figure 3.13: Filter HQ 1000 ¹H BP (-³¹P, ¹⁹F)

Typical Curves




Figure 3.14: Filter HQ 1100 ¹H BP (-³¹P, ¹⁹F)




Figure 3.15: Filter HQ 1200 ¹H BP (-³¹P, ¹⁹F)

4 Chemicool

4.1 Periodic Table

	1																18	
	1a	2											13	14	15	16	17	VIII
1	H	11a											11a	11a	11a	11a	11a	He
2	Li	Be	3	4	5	6	7	8	9	10	11	12	B	C	N	O	F	Ne
3	Na	Mg	IIIb	IVb	Vb	VIb	VIIb	--- VIIIb ---		IIb			Al	Si	P	S	Cl	Ar
4	K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br	Kr
5	Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I	Xe
6	Cs	Ba	*	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po	At	Rn
7	Fr	Ra	**	Unq	Unp	Unh	Ns	Hs	Mt	Uun	Uuu							
			*	La	Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu
			**	Ac	Th	Pa	U	Np	Pu	Am	Cm	Bk	Cf	Es	Fm	Md	No	Lr



Alkali metals	Transition metals	Halogens	Lanthanoids
Alkali earth metals	Other metals	Noble gases	Actinoids

Figure 4.1: Periodic Table of Elements

4.2 NMR Applications Frequencies

4.2.1 NMR Frequency Table

See also

-  [NMR Frequency Table \(1\) \[▶ 239\]](#)
-  [NMR Frequency Table \(2\) \[▶ 240\]](#)

Isotope	NMR Field (Tesla)												
	7,049	9,398	11,747	14,095	16,444	17,619	18,793	19,967	21,142	22,316	23,490	25,840	28,189
¹ H	300,130	400,130	500,130	600,130	700,130	750,130	800,130	850,130	900,130	950,130	1000,130	1100,130	1200,130
² H	46,072	61,423	76,773	92,124	107,475	115,150	122,825	130,501	138,176	145,851	153,527	168,877	184,228
³ H	320,131	426,795	533,459	640,123	746,787	800,120	853,452	906,784	960,116	1013,448	1066,780	1173,444	1280,108
³ He	228,636	304,815	380,994	457,173	533,352	571,441	609,531	647,620	685,710	723,799	761,889	838,068	914,247
⁶ Li	44,167	58,883	73,599	88,315	103,031	110,389	117,747	125,105	132,463	139,821	147,179	161,895	176,611
⁷ Li	116,642	155,506	194,370	233,233	272,097	291,529	310,961	330,393	349,825	369,257	388,689	427,553	466,416
⁹ Be	42,174	56,226	70,278	84,330	98,382	105,408	112,434	119,460	126,485	133,511	140,537	154,589	168,641
¹⁰ B	32,245	42,989	53,732	64,476	75,220	80,592	85,963	91,335	96,707	102,079	107,451	118,194	128,938
¹¹ B	96,294	128,378	160,462	192,546	224,630	240,672	256,714	272,757	288,799	304,841	320,883	352,967	385,051
¹³ C	75,468	100,613	125,758	150,903	176,048	188,621	201,194	213,766	226,339	238,911	251,484	276,629	301,774
¹⁴ N	21,688	28,914	36,140	43,367	50,593	54,206	57,819	61,432	65,045	68,658	72,271	79,498	86,724
¹⁵ N	30,423	40,560	50,696	60,833	70,969	76,038	81,106	86,174	91,243	96,311	101,379	111,516	121,652
¹⁷ O	40,687	54,243	67,800	81,356	94,913	101,691	108,469	115,248	122,026	128,804	135,582	149,139	162,695
¹⁹ F	282,404	376,498	470,592	564,686	658,780	705,827	752,873	799,920	846,967	894,014	941,061	1035,155	1129,249
²¹ Ne	23,693	31,587	39,481	47,376	55,270	59,217	63,164	67,111	71,058	75,006	78,953	86,847	94,741
²³ Na	79,390	105,842	132,294	158,746	185,197	198,423	211,649	224,875	238,101	251,327	264,553	291,005	317,457
²⁵ Mg	18,373	24,495	30,616	36,738	42,860	45,921	48,981	52,042	55,103	58,164	61,225	67,346	73,468
²⁷ Al	78,204	104,261	130,317	156,374	182,431	195,459	208,488	221,516	234,544	247,573	260,601	286,658	312,714
²⁹ Si	59,627	79,494	99,361	119,228	139,095	149,029	158,962	168,896	178,829	188,763	198,696	218,563	238,431
³¹ P	121,495	161,976	202,457	242,937	283,418	303,659	323,899	344,139	364,380	384,620	404,861	445,341	485,822
³³ S	23,038	30,714	38,390	46,066	53,742	57,580	61,418	65,256	69,094	72,932	76,770	84,446	92,122
³⁵ Cl	29,406	39,204	49,002	58,799	68,597	73,496	78,395	83,294	88,193	93,091	97,990	107,788	117,586
³⁷ Cl	24,478	32,634	40,790	48,945	57,101	61,179	65,257	69,335	73,413	77,491	81,569	89,724	97,880
³⁹ K	14,005	18,671	23,338	28,004	32,670	35,003	37,337	39,670	42,003	44,336	46,669	51,335	56,002
⁴¹ K	7,687	10,248	12,809	15,371	17,932	19,213	20,493	21,774	23,054	24,335	25,616	28,177	30,738
⁴³ Ca	20,199	26,929	33,659	40,389	47,119	50,484	53,849	57,214	60,580	63,945	67,310	74,040	80,770
⁴⁵ Sc	72,907	97,199	121,491	145,782	170,074	182,220	194,366	206,512	218,658	230,804	242,950	267,241	291,533
⁴⁷ Ti	16,920	22,558	28,195	33,833	39,470	42,289	45,108	47,927	50,745	53,564	56,383	62,020	67,658
⁴⁹ Ti	16,924	22,563	28,202	33,841	39,480	42,299	45,118	47,938	50,757	53,577	56,396	62,035	67,674
⁵⁰ V	29,924	39,894	49,865	59,835	69,805	74,791	79,776	84,761	89,746	94,731	99,716	109,687	119,657
⁵¹ V	78,943	105,246	131,549	157,852	184,155	197,306	210,458	223,609	236,761	249,912	263,064	289,366	315,669
⁵³ Cr	16,965	22,618	28,270	33,923	39,575	42,401	45,228	48,054	50,880	53,707	56,533	62,185	67,838
⁵⁵ Mn	74,400	99,189	123,979	148,768	173,557	185,952	198,346	210,741	223,136	235,530	247,925	272,714	297,503
⁵⁷ Fe	9,718	12,956	16,194	19,432	22,670	24,289	25,908	27,527	29,146	30,765	32,384	35,621	38,859
⁵⁹ Co	71,212	94,939	118,666	142,393	166,120	177,984	189,847	201,711	213,574	225,438	237,301	261,028	284,755
⁶¹ Ni	26,820	35,756	44,692	53,628	62,565	67,033	71,501	75,969	80,437	84,905	89,373	98,309	107,245
⁶³ Cu	79,581	106,097	132,612	159,128	185,643	198,901	212,159	225,416	238,674	251,932	265,190	291,705	318,221
⁶⁵ Cu	85,248	113,652	142,055	170,459	198,863	213,065	227,266	241,468	255,670	269,872	284,074	312,478	340,881
⁶⁷ Zn	18,779	25,036	31,293	37,550	43,807	46,935	50,064	53,192	56,321	59,449	62,578	68,835	75,092
⁶⁹ Ga	72,035	96,036	120,038	144,039	168,040	180,041	192,041	204,042	216,043	228,043	240,044	264,045	288,046
⁷¹ Ga	91,530	122,027	152,524	183,020	213,517	228,766	244,014	259,262	274,511	289,759	305,007	335,504	366,001
⁷³ Ge	10,469	13,957	17,445	20,933	24,422	26,166	27,910	29,654	31,398	33,142	34,886	38,374	41,862
⁷⁵ As	51,390	68,513	85,635	102,758	119,880	128,442	137,003	145,564	154,125	162,687	171,248	188,371	205,493
⁷⁷ Se	57,239	76,310	95,382	114,453	133,525	143,060	152,596	162,132	171,667	181,203	190,739	209,810	228,882
⁷⁹ Br	75,195	100,249	125,303	150,357	175,412	187,939	200,466	212,993	225,520	238,047	250,574	275,628	300,682
⁸¹ Br	81,055	108,062	135,068	162,075	189,082	202,585	216,088	229,591	243,095	256,598	270,101	297,108	324,115
⁸⁵ Rb	28,977	38,632	48,287	57,941	67,596	72,424	77,251	82,078	86,906	91,733	96,561	106,216	115,870
⁸⁷ Rb	98,204	130,924	163,645	196,365	229,086	245,446	261,806	278,167	294,527	310,887	327,247	359,968	392,688
⁸⁷ Sr	13,007	17,341	21,675	26,008	30,342	32,509	34,676	36,843	39,010	41,177	43,344	47,677	52,011

Isotope	NMR Field (Tesla)												
	7,049	9,398	11,747	14,095	16,444	17,619	18,793	19,967	21,142	22,316	23,490	25,840	28,189
⁸⁹ Y	14,707	19,607	24,507	29,408	34,308	36,758	39,208	41,658	44,108	46,558	49,008	53,909	58,809
⁹¹ Zr	27,901	37,197	46,494	55,790	65,086	69,734	74,383	79,031	83,679	88,327	92,975	102,271	111,568
⁹³ Nb	73,460	97,936	122,412	146,888	171,364	183,602	195,840	208,078	220,316	232,554	244,792	269,268	293,745
⁹⁵ Mo	19,559	26,076	32,593	39,110	45,626	48,885	52,143	55,402	58,660	61,918	65,177	71,694	78,211
⁹⁷ Mo	19,970	26,624	33,278	39,931	46,585	49,912	53,239	56,566	59,893	63,220	66,546	73,200	79,854
⁹⁹ Ru	13,821	18,426	23,031	27,636	32,241	34,544	36,846	39,149	41,451	43,754	46,056	50,661	55,266
¹⁰¹ Ru	15,491	20,652	25,814	30,975	36,137	38,717	41,298	43,879	46,460	49,040	51,621	56,782	61,944
¹⁰³ Rh	9,563	12,749	15,936	19,122	22,308	23,901	25,494	27,088	28,681	30,274	31,867	35,053	38,240
¹⁰⁵ Pd	13,734	18,310	22,886	27,462	32,038	34,326	36,614	38,902	41,190	43,478	45,766	50,342	54,918
¹⁰⁷ Ag	12,149	16,197	20,245	24,293	28,341	30,365	32,389	34,413	36,436	38,460	40,484	44,532	48,580
¹⁰⁹ Ag	13,967	18,621	23,274	27,928	32,582	34,908	37,235	39,562	41,889	44,216	46,543	51,196	55,850
¹¹¹ Cd	63,674	84,889	106,105	127,320	148,536	159,144	169,751	180,359	190,967	201,575	212,182	233,398	254,613
¹¹³ Cd	66,608	88,801	110,994	133,187	155,380	166,477	177,573	188,670	199,766	210,863	221,959	244,152	266,345
¹¹³ In	65,626	87,492	109,358	131,224	153,089	164,022	174,955	185,888	196,821	207,754	218,687	240,553	262,419
¹¹⁵ In	65,766	87,679	109,591	131,504	153,416	164,372	175,329	186,285	197,241	208,197	219,154	241,066	262,979
¹¹⁷ Sn	106,943	142,575	178,207	213,840	249,472	267,288	285,104	302,920	320,736	338,552	356,369	392,001	427,633
¹¹⁹ Sn	111,920	149,211	186,501	223,792	261,082	279,727	298,373	317,018	335,663	354,308	372,954	410,244	447,535
¹²¹ Sb	71,823	95,754	119,684	143,615	167,546	179,511	191,476	203,441	215,407	227,372	239,337	263,268	287,199
¹²³ Sb	38,894	51,853	64,812	77,771	90,730	97,210	103,689	110,169	116,648	123,128	129,607	142,566	155,525
¹²³ Te	78,543	104,713	130,882	157,052	183,222	196,306	209,391	222,476	235,561	248,646	261,731	287,900	314,070
¹²⁵ Te	94,690	126,240	157,789	189,339	220,889	236,663	252,438	268,213	283,988	299,763	315,538	347,087	378,637
¹²⁷ I	60,048	80,055	100,063	120,070	140,077	150,081	160,085	170,088	180,092	190,096	200,099	220,107	240,114
¹²⁹ Xe	83,467	111,277	139,088	166,898	194,708	208,613	222,518	236,424	250,329	264,234	278,139	305,949	333,760
¹³¹ Xe	24,742	32,986	41,230	49,473	57,717	61,839	65,961	70,083	74,205	78,326	82,448	90,692	98,936
¹³³ Cs	39,365	52,481	65,597	78,713	91,829	98,387	104,945	111,503	118,061	124,619	131,177	144,293	157,409
¹³⁵ Ba	29,816	39,750	49,685	59,619	69,553	74,521	79,488	84,455	89,422	94,389	99,357	109,291	119,225
¹³⁷ Ba	33,353	44,466	55,579	66,692	77,804	83,361	88,917	94,474	100,030	105,587	111,143	122,256	133,369
¹³⁸ La	39,600	52,794	65,989	79,183	92,377	98,974	105,571	112,169	118,766	125,363	131,960	145,154	158,349
¹³⁹ La	42,395	56,521	70,646	84,772	98,897	105,960	113,023	120,086	127,148	134,211	141,274	155,399	169,525
¹⁴¹ Pr	91,890	122,507	153,123	183,740	214,357	229,665	244,974	260,282	275,590	290,899	306,207	336,824	367,441
¹⁴³ Nd	16,350	21,798	27,245	32,693	38,141	40,864	43,588	46,312	49,036	51,760	54,483	59,931	65,379
¹⁷⁷ Hf	12,180	16,238	20,296	24,355	28,413	30,442	32,471	34,500	36,529	38,559	40,588	44,646	48,704
¹⁷⁹ Hf	7,650	10,199	12,748	15,297	17,846	19,120	20,394	21,669	22,943	24,218	25,492	28,041	30,590
¹⁸¹ Ta	35,984	47,973	59,963	71,952	83,942	89,937	95,931	101,926	107,921	113,916	119,910	131,900	143,889
¹⁸³ W	12,505	16,672	20,838	25,005	29,171	31,254	33,338	35,421	37,504	39,587	41,671	45,837	50,004
¹⁸⁵ Re	67,603	90,128	112,652	135,177	157,701	168,964	180,226	191,488	202,750	214,013	225,275	247,800	270,324
¹⁸⁷ Re	68,284	91,035	113,787	136,538	159,290	170,666	182,041	193,417	204,793	216,169	227,544	250,296	273,047
¹⁸⁹ Os	23,306	31,071	38,837	46,602	54,367	58,250	62,133	66,015	69,898	73,780	77,663	85,428	93,194
¹⁹¹ Ir	5,400	7,199	8,998	10,798	12,597	13,496	14,396	15,296	16,195	17,095	17,995	19,794	21,593
¹⁹³ Ir	5,860	7,812	9,765	11,717	13,670	14,646	15,622	16,599	17,575	18,551	19,527	21,480	23,432
¹⁹⁵ Pt	64,518	86,015	107,511	129,008	150,505	161,253	172,001	182,750	193,498	204,246	214,995	236,491	257,988
¹⁹⁷ Au	5,310	7,079	8,848	10,618	12,387	13,272	14,156	15,041	15,925	16,810	17,695	19,464	21,233
¹⁹⁹ Hg	53,756	71,667	89,578	107,489	125,400	134,355	143,311	152,266	161,221	170,177	179,132	197,043	214,954
²⁰¹ Hg	19,843	26,454	33,066	39,677	46,289	49,595	52,900	56,206	59,512	62,818	66,123	72,735	79,346
²⁰³ Ti	171,444	228,567	285,690	342,814	399,937	428,499	457,060	485,622	514,183	542,745	571,307	628,430	685,553
²⁰⁵ Ti	173,127	230,811	288,495	346,179	403,863	432,705	461,547	490,389	519,231	548,073	576,915	634,599	692,283
²⁰⁷ Pb	62,789	83,710	104,630	125,551	146,471	156,932	167,392	177,852	188,313	198,773	209,233	230,154	251,074
²⁰⁹ Bi	48,229	64,298	80,368	96,437	112,506	120,541	128,576	136,611	144,645	152,680	160,715	176,784	192,853
²³⁵ U	5,527	7,369	9,210	11,052	12,893	13,814	14,735	15,655	16,576	17,497	18,418	20,259	22,101

4.2.2 With increasing NMR Frequency

See also

- 📄 [With increasing NMR Frequency \(1\) \[▶ 242\]](#)
- 📄 [With increasing NMR Frequency \(2\) \[▶ 243\]](#)

Isotope	NMR Field (Tesla)												
	7,049	9,398	11,747	14,095	16,444	17,619	18,793	19,967	21,142	22,316	23,490	25,840	28,189
¹⁹⁷ Au	5,310	7,079	8,848	10,618	12,387	13,272	14,156	15,041	15,925	16,810	17,695	19,464	21,233
¹⁹¹ Ir	5,400	7,199	8,998	10,798	12,597	13,496	14,396	15,296	16,195	17,095	17,995	19,794	21,593
²³⁵ U	5,527	7,368	9,210	11,052	12,893	13,814	14,735	15,655	16,576	17,497	18,418	20,259	22,101
¹⁹³ Ir	5,860	7,812	9,765	11,717	13,670	14,646	15,622	16,599	17,575	18,551	19,527	21,480	23,432
¹⁷⁹ Hf	7,650	10,199	12,748	15,297	17,846	19,120	20,394	21,669	22,943	24,218	25,492	28,041	30,590
⁴¹ K	7,687	10,248	12,809	15,371	17,932	19,213	20,493	21,774	23,054	24,335	25,616	28,177	30,738
¹⁰³ Rh	9,563	12,749	15,936	19,122	22,308	23,901	25,494	27,088	28,681	30,274	31,867	35,053	38,240
⁵⁷ Fe	9,718	12,956	16,194	19,432	22,670	24,289	25,908	27,527	29,146	30,765	32,384	35,621	38,859
⁷³ Ge	10,469	13,957	17,445	20,933	24,422	26,166	27,910	29,654	31,398	33,142	34,886	38,374	41,862
¹⁰⁷ Ag	12,149	16,196	20,245	24,293	28,341	30,365	32,389	34,413	36,436	38,460	40,484	44,532	48,580
¹⁷⁷ Hf	12,180	16,238	20,296	24,355	28,413	30,442	32,471	34,500	36,529	38,559	40,588	44,646	48,704
¹⁸³ W	12,505	16,671	20,838	25,005	29,171	31,254	33,338	35,421	37,504	39,587	41,671	45,837	50,004
⁸⁷ Sr	13,007	17,340	21,675	26,008	30,342	32,509	34,676	36,843	39,010	41,177	43,344	47,677	52,011
¹⁰⁵ Pd	13,734	18,309	22,886	27,462	32,038	34,326	36,614	38,902	41,190	43,478	45,766	50,342	54,918
⁹⁹ Ru	13,821	18,425	23,031	27,636	32,241	34,544	36,846	39,149	41,451	43,754	46,056	50,661	55,266
¹⁰⁹ Ag	13,967	18,620	23,274	27,928	32,582	34,908	37,235	39,562	41,889	44,216	46,543	51,196	55,850
³⁹ K	14,005	18,671	23,338	28,004	32,670	35,003	37,337	39,670	42,003	44,336	46,669	51,335	56,002
⁸⁹ Y	14,707	19,607	24,507	29,408	34,308	36,758	39,208	41,658	44,108	46,558	49,008	53,909	58,809
¹⁰¹ Ru	15,491	20,652	25,814	30,975	36,137	38,717	41,298	43,879	46,460	49,040	51,621	56,782	61,944
¹⁴³ Nd	16,350	21,797	27,245	32,693	38,141	40,864	43,588	46,312	49,036	51,760	54,483	59,931	65,379
⁴⁷ Ti	16,920	22,557	28,195	33,833	39,470	42,289	45,108	47,927	50,745	53,564	56,383	62,020	67,658
⁴⁹ Ti	16,924	22,562	28,202	33,841	39,480	42,299	45,118	47,938	50,757	53,577	56,396	62,035	67,674
⁵³ Cr	16,965	22,617	28,270	33,923	39,575	42,401	45,228	48,054	50,880	53,707	56,533	62,185	67,838
²⁵ Mg	18,373	24,494	30,616	36,738	42,860	45,921	48,981	52,042	55,103	58,164	61,225	67,346	73,468
⁶⁷ Zn	18,779	25,035	31,293	37,550	43,807	46,935	50,064	53,192	56,321	59,449	62,578	68,835	75,092
⁹⁵ Mo	19,559	26,075	32,593	39,110	45,626	48,885	52,143	55,402	58,660	61,918	65,177	71,694	78,211
²⁰¹ Hg	19,843	26,454	33,066	39,677	46,289	49,595	52,900	56,206	59,512	62,818	66,123	72,735	79,346
⁹⁷ Mo	19,970	26,623	33,278	39,931	46,585	49,912	53,239	56,566	59,893	63,220	66,546	73,200	79,854
⁴³ Ca	20,199	26,928	33,659	40,389	47,119	50,484	53,849	57,214	60,580	63,945	67,310	74,040	80,770
¹⁴ N	21,688	28,913	36,140	43,367	50,593	54,206	57,819	61,432	65,045	68,658	72,271	79,498	86,724
³³ S	23,038	30,713	38,390	46,066	53,742	57,580	61,418	65,256	69,094	72,932	76,770	84,446	92,122
¹⁸⁹ Os	23,306	31,070	38,837	46,602	54,367	58,250	62,133	66,015	69,898	73,780	77,663	85,428	93,194
²¹ Ne	23,693	31,586	39,481	47,376	55,270	59,217	63,164	67,111	71,058	75,006	78,953	86,847	94,741
³⁷ Cl	24,478	32,633	40,790	48,945	57,101	61,179	65,257	69,335	73,413	77,491	81,569	89,724	97,880
¹³¹ Xe	24,742	32,985	41,230	49,473	57,717	61,839	65,961	70,083	74,205	78,326	82,448	90,692	98,936
⁶¹ Ni	26,820	35,755	44,692	53,628	62,565	67,033	71,501	75,969	80,437	84,905	89,373	98,309	107,245
⁹¹ Zr	27,901	37,196	46,494	55,790	65,086	69,734	74,383	79,031	83,679	88,327	92,975	102,271	111,568
⁸⁵ Rb	28,977	38,631	48,287	57,941	67,596	72,424	77,251	82,078	86,906	91,733	96,561	106,216	115,870
³⁵ Cl	29,406	39,202	49,002	58,799	68,597	73,496	78,395	83,294	88,193	93,091	97,990	107,788	117,586
¹³⁵ Ba	29,816	39,749	49,685	59,619	69,553	74,521	79,488	84,455	89,422	94,389	99,357	109,291	119,225
⁵⁰ V	29,924	39,893	49,865	59,835	69,805	74,791	79,776	84,761	89,746	94,731	99,716	109,687	119,657
¹⁵ N	30,423	40,558	50,696	60,833	70,969	76,038	81,106	86,174	91,243	96,311	101,379	111,516	121,652
¹⁰ B	32,245	42,987	53,732	64,476	75,220	80,592	85,963	91,335	96,707	102,079	107,451	118,194	128,938
¹³⁷ Ba	33,353	44,464	55,579	66,692	77,804	83,361	88,917	94,474	100,030	105,587	111,143	122,256	133,369
¹⁸¹ Ta	35,984	47,972	59,963	71,952	83,942	89,937	95,931	101,926	107,921	113,916	119,910	131,900	143,889
¹²³ Sb	38,894	51,851	64,812	77,771	90,730	97,210	103,689	110,169	116,648	123,128	129,607	142,566	155,525
¹³³ Cs	39,365	52,479	65,597	78,713	91,829	98,387	104,945	111,503	118,061	124,619	131,177	144,293	157,409
¹³⁸ La	39,600	52,793	65,989	79,183	92,377	98,974	105,571	112,169	118,766	125,363	131,960	145,154	158,349
¹⁷ O	40,687	54,242	67,800	81,356	94,913	101,691	108,469	115,248	122,026	128,804	135,582	149,139	162,695

Isotope	NMR Field (Tesla)												
	7,049	9,398	11,747	14,095	16,444	17,619	18,793	19,967	21,142	22,316	23,490	25,840	28,189
⁹ Be	42,174	56,226	70,278	84,330	98,382	105,408	112,434	119,460	126,485	133,511	140,537	154,589	168,641
¹³⁹ La	42,395	56,521	70,646	84,772	98,897	105,960	113,023	120,086	127,148	134,211	141,274	155,399	169,525
⁶ Li	44,167	58,883	73,599	88,315	103,031	110,389	117,747	125,105	132,463	139,821	147,179	161,895	176,611
² H	46,072	61,423	76,773	92,124	107,475	115,150	122,825	130,501	138,176	145,851	153,527	168,877	184,228
²⁰⁹ Bi	48,229	64,298	80,368	96,437	112,506	120,541	128,576	136,611	144,645	152,680	160,715	176,784	192,853
⁷⁵ As	51,390	68,513	85,635	102,758	119,880	128,442	137,003	145,564	154,125	162,687	171,248	188,371	205,493
¹⁹⁹ Hg	53,756	71,667	89,578	107,489	125,400	134,355	143,311	152,266	161,221	170,177	179,132	197,043	214,954
⁷⁷ Se	57,239	76,310	95,382	114,453	133,525	143,060	152,596	162,132	171,667	181,203	190,739	209,810	228,882
²⁹ Si	59,627	79,494	99,361	119,228	139,095	149,029	158,962	168,896	178,829	188,763	198,696	218,563	238,431
¹²⁷ I	60,048	80,055	100,063	120,070	140,077	150,081	160,085	170,088	180,092	190,096	200,099	220,107	240,114
²⁰⁷ Pb	62,789	83,710	104,630	125,551	146,471	156,932	167,392	177,852	188,313	198,773	209,233	230,154	251,074
¹¹¹ Cd	63,674	84,889	106,105	127,320	148,536	159,144	169,751	180,359	190,967	201,575	212,182	233,398	254,613
¹⁹⁵ Pt	64,518	86,015	107,511	129,008	150,505	161,253	172,001	182,750	193,498	204,246	214,995	236,491	257,988
¹¹³ In	65,626	87,492	109,358	131,224	153,089	164,022	174,955	185,888	196,821	207,754	218,687	240,553	262,419
¹¹⁵ In	65,766	87,679	109,591	131,504	153,416	164,372	175,329	186,285	197,241	208,197	219,154	241,066	262,979
¹¹³ Cd	66,608	88,801	110,994	133,187	155,380	166,477	177,573	188,670	199,766	210,863	221,959	244,152	266,345
¹⁸⁵ Re	67,603	90,128	112,652	135,177	157,701	168,964	180,226	191,488	202,750	214,013	225,275	247,800	270,324
¹⁸⁷ Re	68,284	91,035	113,787	136,538	159,290	170,666	182,041	193,417	204,793	216,169	227,544	250,296	273,047
⁵⁹ Co	71,212	94,939	118,666	142,393	166,120	177,984	189,847	201,711	213,574	225,438	237,301	261,028	284,755
¹²¹ Sb	71,823	95,754	119,684	143,615	167,546	179,511	191,476	203,441	215,407	227,372	239,337	263,268	287,199
⁶⁹ Ga	72,035	96,036	120,038	144,039	168,040	180,041	192,041	204,042	216,043	228,043	240,044	264,045	288,046
⁴⁵ Sc	72,907	97,199	121,491	145,782	170,074	182,220	194,366	206,512	218,658	230,804	242,950	267,241	291,533
⁹³ Nb	73,460	97,936	122,412	146,888	171,364	183,602	195,840	208,078	220,316	232,554	244,792	269,268	293,745
⁵⁵ Mn	74,400	99,189	123,979	148,768	173,557	185,952	198,346	210,741	223,136	235,530	247,925	272,714	297,503
⁷⁹ Br	75,195	100,249	125,303	150,357	175,412	187,939	200,466	212,993	225,520	238,047	250,574	275,628	300,682
¹³ C	75,468	100,613	125,758	150,903	176,048	188,621	201,194	213,766	226,339	238,911	251,484	276,629	301,774
²⁷ Al	78,204	104,261	130,317	156,374	182,431	195,459	208,488	221,516	234,544	247,573	260,601	286,658	312,714
¹²² Te	78,543	104,713	130,882	157,052	183,222	196,306	209,391	222,476	235,561	248,646	261,731	287,900	314,070
⁵¹ V	78,943	105,246	131,549	157,852	184,155	197,306	210,458	223,609	236,761	249,912	263,064	289,366	315,669
²³ Na	79,390	105,842	132,294	158,746	185,197	198,423	211,649	224,875	238,101	251,327	264,553	291,005	317,457
⁶³ Cu	79,581	106,097	132,612	159,128	185,643	198,901	212,159	225,416	238,674	251,932	265,190	291,705	318,221
⁸¹ Br	81,055	108,062	135,068	162,075	189,082	202,585	216,088	229,591	243,095	256,598	270,101	297,108	324,115
¹²⁹ Xe	83,467	111,277	139,088	166,898	194,708	208,613	222,518	236,424	250,329	264,234	278,139	305,949	333,760
⁶⁵ Cu	85,248	113,652	142,055	170,459	198,863	213,065	227,266	241,468	255,670	269,872	284,074	312,478	340,881
⁷¹ Ga	91,530	122,027	152,524	183,020	213,517	228,766	244,014	259,262	274,511	289,759	305,007	335,504	366,001
¹⁴¹ Pr	91,890	122,507	153,123	183,740	214,357	229,665	244,974	260,282	275,590	290,899	306,207	336,824	367,441
¹²⁵ Te	94,690	126,240	157,789	189,339	220,889	236,663	252,438	268,213	283,988	299,763	315,538	347,087	378,637
¹¹ B	96,294	128,378	160,462	192,546	224,630	240,672	256,714	272,757	288,799	304,841	320,883	352,967	385,051
⁸⁷ Rb	98,204	130,924	163,645	196,365	229,086	245,446	261,806	278,167	294,527	310,887	327,247	359,968	392,688
¹¹⁷ Sn	106,943	142,575	178,207	213,840	249,472	267,288	285,104	302,920	320,736	338,552	356,369	392,001	427,633
¹¹⁹ Sn	111,920	149,211	186,501	223,792	261,082	279,727	298,373	317,018	335,663	354,308	372,954	410,244	447,535
⁷ Li	116,642	155,506	194,370	233,233	272,097	291,529	310,961	330,393	349,825	369,257	388,689	427,553	466,416
³¹ P	121,495	161,976	202,457	242,937	283,418	303,659	323,899	344,139	364,380	384,620	404,861	445,341	485,822
²⁰³ Tl	171,444	228,567	285,690	342,814	399,937	428,499	457,060	485,622	514,183	542,745	571,307	628,430	685,553
²⁰⁵ Tl	173,127	230,811	288,495	346,179	403,863	432,705	461,547	490,389	519,231	548,073	576,915	634,599	692,283
³ He	228,636	304,815	380,994	457,173	533,352	571,441	609,531	647,620	685,710	723,799	761,889	838,068	914,247
¹⁹ F	282,404	376,498	470,592	564,686	658,780	705,827	752,873	799,920	846,967	894,014	941,061	1035,155	1129,249
¹ H	300,130	400,130	500,130	600,130	700,130	750,130	800,130	850,130	900,130	950,130	1000,130	1100,130	1200,130
³ H	320,131	426,795	533,459	640,123	746,787	800,120	853,452	906,784	960,116	1013,448	1066,780	1173,444	1280,108

5 Contact

Manufacturer:

Bruker BioSpin SAS
34, rue de l'Industrie
67166 WISSEMBOURG Cedex
France
Phone: + 33 3 88 06 60 00
Fax: + 33 3 88 06 68 20
<http://www.bruker.com>
WEEE FR311020911

NMR Hotlines:

Contact our NMR service centers.

Bruker BioSpin NMR provides dedicated hotlines and service centers, so that our specialists can respond as quickly as possible to all your service requests, applications questions, software or technical needs.

Please select the NMR service center or hotline you wish to contact from our list available at:

<http://www.bruker.com/service/information-communication/helpdesk.html>

List of Figures

Figure 1.1:	AVANCE 300MHz [▶ 30]	8
Figure 1.2:	AVANCE 400MHz [▶ 43]	8
Figure 1.3:	AVANCE 500MHz [▶ 57]	8
Figure 1.4:	AVANCE 600MHz [▶ 77]	8
Figure 1.5:	AVANCE 700MHz [▶ 93]	9
Figure 1.6:	AVANCE 750MHz [▶ 107]	9
Figure 1.7:	AVANCE 800MHz [▶ 119]	9
Figure 1.8:	AVANCE 850MHz [▶ 129]	9
Figure 1.9:	AVANCE 900MHz [▶ 141]	10
Figure 1.10:	AVANCE 950MHz [▶ 152]	10
Figure 1.11:	AVANCE 1000MHz [▶ 163]	10
Figure 1.12:	AVANCE 1100MHz [▶ 174]	10
Figure 1.13:	AVANCE 1200MHz [▶ 183]	11
Figure 1.14:	¹ H/ ¹⁹ F Typical Configuration	27
Figure 1.15:	¹ H/ ¹⁹ F/X NMR Experiments Configuration.....	28
Figure 1.16:	Different High Power Filters for the ¹ H/ ¹⁹ F Experiment.....	28
Figure 1.17:	Lowpass Filter HQ 300MHz Diagram	30
Figure 1.18:	Bandstop Filter HQ 300MHz Diagram	32
Figure 1.19:	Bandpass Filter HQ 300MHz Diagram	34
Figure 1.20:	Diplexer Filter HQ 300MHz Diagram	42
Figure 1.21:	Lowpass Filter HQ 400MHz Diagram	43
Figure 1.22:	Bandstop Filter HQ 400MHz Diagram	45
Figure 1.23:	Bandpass Filter HQ 400MHz Diagram	47
Figure 1.24:	Diplexer Filter HQ 400MHz Diagram	56
Figure 1.25:	Lowpass Filter HQ 500MHz Diagram	57
Figure 1.26:	Bandstop Filter HQ 500MHz Diagram	59
Figure 1.27:	Band Pass Filter HQ 500MHz Diagram	64
Figure 1.28:	Diplexer Filter HQ 500MHz Diagram	75
Figure 1.29:	Lowpass Filter HQ 600 MHz Diagram	77
Figure 1.30:	Bandstop Filter HQ 600 MHz Diagram	79
Figure 1.31:	Bandpass Filter HQ 600 MHz Diagram	82
Figure 1.32:	Diplexer Filter HQ 600 MHz Diagram	92
Figure 1.33:	Lowpass Filter HQ 700MHz Diagram	93
Figure 1.34:	Bandstop Filter HQ 700MHz Diagram	95
Figure 1.35:	Bandpass Filter HQ 700MHz Diagram	97
Figure 1.36:	Diplexer Filter HQ 700MHz Diagram	106
Figure 1.37:	Lowpass Filter HQ 750MHz Diagram	107
Figure 1.38:	Bandstop Filter HQ 750MHz Diagram	109
Figure 1.39:	Bandpass Filter HQ 750MHz Diagram	111

List of Figures

Figure 1.40:	Diplexer Filter HQ 750MHz Diagram	118
Figure 1.41:	Lowpass Filter HQ 800MHz Diagram	119
Figure 1.42:	Bandstop Filter HQ 800MHz Diagram	121
Figure 1.43:	Bandpass Filter HQ 800MHz Diagram	123
Figure 1.44:	Lowpass Filter HQ 850MHz Diagram	129
Figure 1.45:	Bandpass Filter HQ 850MHz Diagram	133
Figure 1.46:	Diplexer Filter HQ 850MHz Diagram	140
Figure 1.47:	Lowpass Filter HQ 900MHz Diagram	141
Figure 1.48:	Bandstop Filter HQ 900MHz Diagram	143
Figure 1.49:	Bandpass Filter HQ 900MHz Diagram	144
Figure 1.50:	Diplexer Filter HQ 900MHz Diagram	151
Figure 1.51:	Lowpass Filter HQ 950MHz Diagram	152
Figure 1.52:	Bandstop Filter HQ 950MHz Diagram	154
Figure 1.53:	Bandpass Filter HQ 950MHz Diagram	155
Figure 1.54:	Diplexer Filter HQ 950MHz Diagram	162
Figure 1.55:	Lowpass Filter HQ 1000MHz Diagram	163
Figure 1.56:	Bandstop Filter HQ 1000MHz Diagram	165
Figure 1.57:	Bandpass Filter HQ 1000MHz Diagram	167
Figure 1.58:	Diplexer Filter HQ 1000MHz Diagram	173
Figure 1.59:	Lowpass Filter HQ 1100MHz Diagram	174
Figure 1.60:	Bandpass Filter HQ 1100MHz Diagram	176
Figure 1.61:	Diplexer Filter HQ 1100MHz Diagram	182
Figure 1.62:	Lowpass Filter HQ 1200MHz Diagram	183
Figure 1.63:	Bandpass Filter HQ 1200MHz Diagram	185
Figure 1.64:	Diplexer Filter HQ 1200MHz Diagram	192
Figure 2.1:	3S063 View	215
Figure 2.2:	4S063D View	215
Figure 2.3:	5R032 View	216
Figure 2.4:	5R045 View	216
Figure 2.5:	5R069 View	217
Figure 2.6:	5R100 View	217
Figure 2.7:	6R125 View	218
Figure 2.8:	7R155 View	218
Figure 2.9:	8R055P View	219
Figure 2.10:	8R105P View	219
Figure 2.11:	4H073 View	221
Figure 2.12:	5H133 View	221
Figure 2.13:	5H069 View	222
Figure 2.14:	6H073 View	222
Figure 2.15:	8H073 View	223
Figure 2.16:	8H097 View	223
Figure 2.17:	8H139 View	224
Figure 2.18:	2D050 View	226

Figure 2.19:	2D095 View	226
Figure 2.20:	W1303715 - N / Female-Female 0,5m Coaxial Cable	227
Figure 2.21:	30649 - N / Female-Female Adapter	227
Figure 2.22:	4241 - N / Male-Male Adapter.....	227
Figure 2.23:	30209 - N / Male-Female Right Angle Adapter.....	227
Figure 3.1:	W1346239 - Filter HQ 300 ¹ H BP (- ³¹ P, ¹⁹ F)	229
Figure 3.2:	W1346249 - Filter HQ 300 ¹⁹ F BP (- ³¹ P, ¹ H)	230
Figure 3.3:	W1346240 - Filter HQ 400 ¹ H BP (- ³¹ P, ¹⁹ F)	230
Figure 3.4:	W1346250 - Filter HQ 400 ¹⁹ F BP (- ³¹ P, ¹ H)	231
Figure 3.5:	W1346223 - Filter HQ 500 ¹ H BP(- ³¹ P, ¹⁹ F)	231
Figure 3.6:	W1346251 - Filter HQ 500 ¹⁹ F BP (- ³¹ P, ¹ H)	232
Figure 3.7:	W1346241 - Filter HQ 600 ¹ H BP (- ³¹ P, ¹⁹ F)	232
Figure 3.8:	W1346296 - Filter HQ 700 ¹ H BP (- ³¹ P, ¹⁹ F)	233
Figure 3.9:	W1346300 - Filter HQ 700 ¹⁹ F BP (- ³¹ P, ¹ H)	233
Figure 3.10:	W1346297 - Filter HQ 750 ¹ H BP (- ³¹ P, ¹⁹ F)	234
Figure 3.11:	W1346298 - Filter HQ 800 ¹ H BP (- ³¹ P, ¹⁹ F)	234
Figure 3.12:	W1346299- Filter HQ 900 ¹ H BP (- ³¹ P, ¹⁹ F)	235
Figure 3.13:	Filter HQ 1000 ¹ H BP (- ³¹ P, ¹⁹ F).....	235
Figure 3.14:	Filter HQ 1100 ¹ H BP (- ³¹ P, ¹⁹ F).....	236
Figure 3.15:	Filter HQ 1200 ¹ H BP (- ³¹ P, ¹⁹ F).....	236
Figure 4.1:	Periodic Table of Elements.....	237

List of Tables

Table 1.1:	W1346642 - Filter HQ 300 0- ³¹ P LP (¹⁹ F- ¹ H).....	30
Table 1.2:	W135032 - Filter HQ 300 0 - ¹⁵ N NR LP (² H).....	31
Table 1.3:	W1346253 - Filter HQ 300 ¹ H S.....	32
Table 1.4:	W1346369 - Filter HQ 300 ¹⁹ F S.....	33
Table 1.5:	W1346239 - Filter HQ 300 ¹ H BP (- ³¹ P, ¹⁹ F).....	34
Table 1.6:	W1346235 - Filter HQ 300 ¹ H BP (¹⁹ F).....	35
Table 1.7:	W1346249 - Filter HQ 300 ¹⁹ F BP (- ³¹ P, ¹ H).....	35
Table 1.8:	W1346236 - Filter HQ 300 ¹⁹ F BP (¹ H).....	36
Table 1.9:	W1346368 - Filter HQ 300 ¹⁹ F - ¹ H BP (- ³¹ P).....	36
Table 1.10:	W1346372 - Filter HQ 300 ⁷ Li - ³¹ P BP (- ¹³ C, ¹⁹ F-).....	37
Table 1.11:	W1346270 - Filter HQ 300 ⁷¹ Ga - ⁸⁷ Rb BP (- ¹³ C, ³¹ P-).....	37
Table 1.12:	W1346271 - Filter HQ 300 ⁵⁹ Co - ²³ Na BP (- ²⁹ Si, ⁸⁷ Rb-).....	38
Table 1.13:	W1346275 - Filter HQ 300 ⁵⁹ Co - ²³ Na BP (- ¹⁵ N).....	38
Table 1.14:	W1346272 - Filter HQ 300 ²⁹ Si BP (- ² H, ¹³ C-).....	39
Table 1.15:	W1346452 - Filter HQ 300 ¹⁷ O - ² H BP (¹¹ B-).....	39
Table 1.16:	W1346273 - Filter HQ 300 ¹³³ Cs - ² H BP (- ¹⁵ N, ²⁹ Si-).....	40
Table 1.17:	W1346274 - Filter HQ 300 ¹⁵ N BP (- ¹⁴ N, ² H-).....	40
Table 1.18:	W1346278 - Filter HQ 300 ¹⁵ N BP (¹³ C-).....	41
Table 1.19:	W1346237 - Diplexer HQ 300 ¹ H / ¹⁹ F.....	42
Table 1.20:	W1346636 - Filter HQ 400 0- ³¹ P LP (¹⁹ F- ¹ H).....	43
Table 1.21:	W1346628 - Filter HQ 400 0- ¹⁵ N NR LP (¹³ C).....	44
Table 1.22:	W1346370 - Filter HQ 400 ¹⁹ F S.....	45
Table 1.23:	W1346437 - Filter HQ 400 ² H P ⁷⁷ Se S.....	46
Table 1.24:	W1346240 - Filter HQ 400 ¹ H BP (- ³¹ P, ¹⁹ F).....	47
Table 1.25:	W1346204 - Filter HQ 400 ¹ H BP (¹⁹ F).....	48
Table 1.26:	W1346250 - Filter HQ 400 ¹⁹ F BP (- ³¹ P, ¹ H).....	48
Table 1.27:	W1346205 - Filter HQ 400 ¹⁹ F BP (¹ H).....	49
Table 1.28:	W1346373 - Filter HQ 400 ¹⁹ F - ¹ H BP (- ³¹ P).....	49
Table 1.29:	W1346256 - Filter HQ 400 ⁷ Li - ³¹ P BP (- ⁸⁷ Rb, ¹⁹ F-).....	50
Table 1.30:	W1346349 - Filter HQ 400 ⁷ Li - ³¹ P BP (- ¹³ C, ¹⁹ F-).....	50
Table 1.31:	W1346257 - Filter HQ 400 ⁷¹ Ga - ⁸⁷ Rb BP (- ¹³ C, ³¹ P-).....	51
Table 1.32:	W1346488 - Filter HQ 400 ¹³ C - ²³ Na BP (- ²⁹ Si, ⁷¹ Ga-).....	51
Table 1.33:	W1346258 - Filter HQ 400 ⁵⁹ Co - ²³ Na BP (- ²⁹ Si, ⁸⁷ Rb-).....	52
Table 1.34:	W1346440 - Filter HQ 400 ⁵⁹ Co - ²³ Na BP (- ² H, ¹⁹ F-).....	52
Table 1.35:	W1346226 - Filter HQ 400 ⁵⁹ Co - ²³ Na BP (- ¹⁵ N).....	53
Table 1.36:	W1346259 - Filter HQ 400 ²⁹ Si BP (- ² H, ¹³ C-).....	53
Table 1.37:	W1346442 - Filter HQ 400 ¹⁷ O - ² H BP (¹¹ B-).....	54
Table 1.38:	W1346260 - Filter HQ 400 ¹³³ Cs - ² H BP (- ¹⁵ N, ²⁹ Si-).....	54
Table 1.39:	W1346261 - Filter HQ 400 ¹⁵ N BP (- ¹⁴ N, ² H-).....	55

List of Tables

Table 1.40: W1346439 - Filter HQ 400 ^{23}N ^{65}Cu BP (^{-29}Si , ^{31}P -)	55
Table 1.41: W1346345 - Diplexer HQ 400 ^1H / ^{19}F	56
Table 1.42: W1346637 - Filter HQ 500 0 - ^{31}P LP (^{19}F - ^1H).....	57
Table 1.43: W1346629 - Filter HQ 500 0 - ^{15}N NR LP (^{13}C).....	58
Table 1.44: W1346232 - Filter HQ 500 ^1H S.....	59
Table 1.45: W1346358 - Filter HQ 500 ^{19}F S.....	60
Table 1.46: W1346430 - Filter HQ 500 ^7Li P ^{31}P S	60
Table 1.47: W1346431 - Filter HQ 500 ^{31}P P ^7Li S	61
Table 1.48: W1346434 - Filter HQ 500 ^6Li P ^2H S	61
Table 1.49: W1346435 - Filter HQ 500 ^2H P ^6Li S	62
Table 1.50: W130896- Filter HQ 500 ^{13}C -P ^{27}Al - S	62
Table 1.51: W130897- Filter HQ 500 ^{27}Al -P ^{13}C - S	63
Table 1.52: W1346223 - Filter HQ 500 ^1H BP (^{-31}P , ^{19}F)	64
Table 1.53: W1346247 - Filter HQ 500 ^1H BP (^{19}F)	65
Table 1.54: W1346251 - Filter HQ 500 ^{19}F BP (^{-31}P , ^1H)	65
Table 1.55: W1346248 - Filter HQ 500 ^{19}F BP (^1H)	66
Table 1.56: W1346374 - Filter HQ 500 ^{19}F - ^1H BP (^{-31}P)	66
Table 1.57: W1346262 - Filter HQ 500 ^7Li - ^{31}P BP (^{-87}Rb , ^{19}F -).....	67
Table 1.58: W1346351 - Filter HQ 500 ^7Li - ^{31}P BP (^{-13}C , ^{19}F -).....	67
Table 1.59: W1346263 - Filter HQ 500 ^{71}Ga - ^{87}Rb BP (^{-13}C , ^{31}P -)	68
Table 1.60: W1346489 - Filter HQ 500 ^{13}C - ^{23}Na BP (^{-29}Si , ^{71}Ga -).....	68
Table 1.61: W1346264 - Filter HQ 500 ^{59}Co - ^{23}Na BP (^{-29}Si , ^{87}Rb -).....	69
Table 1.62: W1346460 - Filter HQ 500 ^{59}Co - ^{23}Na BP (^{-2}H , ^{19}F -).....	69
Table 1.63: W1346224 - Filter HQ 500 ^{59}Co - ^{23}Na BP (^{-15}N)	70
Table 1.64: W1346265 - Filter HQ 500 ^{77}Se - ^{29}Si BP (^{-2}H , ^{13}C -).....	70
Table 1.65: W1346453 - Filter HQ 500 ^{17}O - ^2H BP (^{11}B -).....	71
Table 1.66: W1346266 - Filter HQ 500 ^{133}Cs - ^2H BP (^{-15}N , ^{29}Si -).....	71
Table 1.67: W1346267 - Filter HQ 500 ^{15}N BP (^{-14}N , ^2H -).....	72
Table 1.68: W1346692 - Filter HQ 500 ^1H BP (^3H)	72
Table 1.69: W1346693 - Filter HQ 500 ^1H BP (^{-31}P , ^3H)	73
Table 1.70: W1346694 - Filter HQ 500 ^3H BP (^1H)	73
Table 1.71: W1346695 - Filter HQ 500 ^3H BP (^{-31}P , ^1H)	74
Table 1.72: W1346730 - Filter HQ 500 ^{23}Na - ^{65}Cu BP (^{-29}Si , ^{31}P -).....	74
Table 1.73: W1346346 - Diplexer HQ 500 ^1H / ^{19}F	75
Table 1.74: W1346696 - Diplexer HQ 500 ^1H / ^3H	76
Table 1.75: W136807- Diplexer HQ 500 ^{13}C / ^{27}Al	76
Table 1.76: W1346638 - Filter HQ 600 0 - ^{31}P LP (^{19}F - ^1H).....	77
Table 1.77: W1346630 - Filter HQ 600 0 - ^{15}N NR LP (^{13}C).....	78
Table 1.78: W134543 - Filter HQ 600 0 - ^{14}N NR (^{109}Ag H^3) LP (^{13}C).....	78
Table 1.79: W1346255 - Filter HQ 600 ^1H S.....	79
Table 1.80: W1346371 - Filter HQ 600 ^{19}F S.....	80
Table 1.81: W1346679 - Filter HQ 600 ^6Li - P ^2H - S.....	80
Table 1.82: W1346680 - Filter HQ 600 ^2H - P ^6Li - S.....	81

Table 1.83: W1346241 - Filter HQ 600 ^1H BP (^{-31}P , ^{19}F)	82
Table 1.84: W1346228 - Filter HQ 600 ^1H BP (^{19}F)	83
Table 1.85: W1346252 - Filter HQ 600 ^{19}F BP (^{-31}P , ^1H)	83
Table 1.86: W1346229 - Filter HQ 600 ^{19}F BP (^1H)	84
Table 1.87: W1346375 - Filter HQ 600 ^{19}F - ^1H BP (^{31}P)	84
Table 1.88: W1346279 - Filter HQ 600 ^7Li - ^{31}P BP (^{-87}Rb , ^{19}F -).....	85
Table 1.89: W1346367 - Filter HQ 600 ^7Li - ^{31}P BP (^{-13}C , ^{19}F -).....	85
Table 1.90: W1346280 - Filter HQ 600 ^{71}Ga - ^{87}Rb BP (^{-13}C , ^{31}P -).....	86
Table 1.91: W1346490 - Filter HQ 600 ^{13}C - ^{23}Na BP (^{-29}Si , ^{71}Ga).....	86
Table 1.92: W1346276 - Filter HQ 600 ^{59}Co - ^{23}Na BP (^{-29}Si , ^{87}Rb -).....	87
Table 1.93: W1346461 - Filter HQ 600 ^{59}Co - ^{23}Na BP (^{-2}H , ^{19}F -).....	87
Table 1.94: W1346230 - Filter HQ 600 ^{59}Co - ^{23}Na BP (^{-15}N)	88
Table 1.95: W1346281 - Filter HQ 600 ^{29}Si BP (^{-2}H , ^{13}C -).....	88
Table 1.96: W1346454 - Filter HQ 600 ^{17}O - ^2H BP (^{11}B -).....	89
Table 1.97: W1346282 - Filter HQ 600 ^{133}Cs - ^2H BP (^{-15}N , ^{29}Si -).....	89
Table 1.98: W1346277 - Filter HQ 600 ^{15}N BP (^{-14}N , ^2H -).....	90
Table 1.99: W1346283 - Filter HQ 600 ^{14}N BP (^{15}N -).....	90
Table 1.100: W1346757 - Filter HQ 600 ^{29}Si - ^{207}Pb BP (^{-2}H , ^{13}C -).....	91
Table 1.101: W1346347 - Diplexer HQ 600 ^1H / ^{19}F	92
Table 1.102: W1346639 - Filter HQ 700 0- ^{31}P LP (^{19}F - ^1H).....	93
Table 1.103: W1346659 - Filter HQ 700 0- ^{15}N NR LP (^{13}C).....	94
Table 1.104: W1346284 - Filter HQ 700 ^1H S.....	95
Table 1.105: W1346387 - Filter HQ 700 ^{19}F S.....	96
Table 1.106: W1346394 - Filter HQ 700 ^2H S.....	96
Table 1.107: W1346296 - Filter HQ 700 ^1H BP (^{-31}P , ^{19}F)	97
Table 1.108: W1346288 - Filter HQ 700 ^1H BP (^{19}F)	98
Table 1.109: W1346300 - Filter HQ 700 ^{19}F BP (^{-31}P , ^1H)	98
Table 1.110: W1346292 - Filter HQ 700 ^{19}F BP (^1H)	99
Table 1.111: W1346304 - Filter HQ 700 ^7Li - ^{31}P BP (^{-87}Rb , ^{19}F -).....	99
Table 1.112: W1346308 - Filter HQ 700 ^{71}Ga - ^{87}Rb BP (^{-13}C , ^7Li -).....	100
Table 1.113: W1346491 - Filter HQ 700 ^{13}C - ^{23}Na BP (^{-29}Si , ^{71}Ga -).....	100
Table 1.114: W1346312 - Filter HQ 700 ^{59}Co - ^{23}Na BP (^{-29}Si , ^{87}Rb -).....	101
Table 1.115: W1346462 - Filter HQ 700 ^{59}Co - ^{23}Na BP (^{-2}H , ^{19}F -).....	101
Table 1.116: W1346316 - Filter HQ 700 ^{59}Co - ^{23}Na BP (^{-15}N)	102
Table 1.117: W1346320 - Filter HQ 700 ^{77}Se - ^{29}Si BP (^{-2}H , ^{13}C -).....	102
Table 1.118: W1346455 - Filter HQ 700 ^{17}O - ^2H BP (^{11}B -).....	103
Table 1.119: W1346492 - Filter HQ 700 ^{17}O - ^2H BP (^{-15}N , ^{29}Si -).....	103
Table 1.120: W1346324 - Filter HQ 700 ^{133}Cs - ^2H BP (^{-15}N , ^{29}Si -).....	104
Table 1.121: W1346328 - Filter HQ 700 ^{15}N BP (^{-14}N , ^2H).....	104
Table 1.122: W1346336 - Filter HQ 700 ^{14}N BP (^{15}N -).....	105
Table 1.123: W1346421 - Diplexer HQ 700 ^1H / ^{19}F	106
Table 1.124: W1346644 - Filter HQ 750 0- ^{31}P LP (^{19}F - ^1H).....	107
Table 1.125: W1346660 - Filter HQ 750 0- ^{15}N NR LP (^{13}C).....	108

List of Tables

Table 1.126: W1346285 Filter HQ 750 ^1H S.....	109
Table 1.127: W1346385 - Filter HQ 750 ^{19}F S.....	110
Table 1.128: W1346297 - Filter HQ 750 ^1H BP (^{-31}P , ^{19}F).....	111
Table 1.129: W1346289 - Filter HQ 750 ^1H BP (^{19}F).....	112
Table 1.130: W1346301 - Filter HQ 750 ^{19}F BP (^{-31}P , ^1H).....	112
Table 1.131: W1346293 - Filter HQ 750 ^{19}F BP (^1H).....	113
Table 1.132: W1346305 - Filter HQ 750 ^7Li - ^{31}P BP (^{-87}Rb , ^{19}F -).....	113
Table 1.133: W1346309 - Filter HQ 750 ^{71}Ga - ^{87}Rb BP (^{-13}C , ^7Li -).....	114
Table 1.134: W1346313 - Filter HQ 750 ^{13}C - ^{23}Na BP (^{-29}Si , ^{71}Ga -).....	114
Table 1.135: W1346317 - Filter HQ 750 ^{59}Co - ^{23}Na BP (^{-15}N).....	115
Table 1.136: W1346321 - Filter HQ 750 ^{77}Se - ^{29}Si BP (^{-2}H , ^{13}C -).....	115
Table 1.137: W1346493 - Filter HQ 750 ^{17}O - ^2H BP (^{-15}N , ^{29}Si -).....	116
Table 1.138: W1346325 - Filter HQ 750 ^{133}Cs - ^2H BP (^{-15}N , ^{29}Si -).....	116
Table 1.139: W1346329 - Filter HQ 750 ^{15}N BP (^{-14}N , ^2H -).....	117
Table 1.140: W1346378 - Diplexer HQ 750 ^1H / ^{19}F	118
Table 1.141: W1346640 - Filter HQ 800 0- ^{31}P LP (^{19}F - ^1H).....	119
Table 1.142: W1346661 - Filter HQ 800 0- ^{15}N NR LP (^{13}C).....	120
Table 1.143: W1346286 - Filter HQ 800 ^1H S.....	121
Table 1.144: W1346386 - Filter HQ 800 ^{19}F S.....	122
Table 1.145: W1346298 - Filter HQ 800 ^1H BP (^{-31}P , ^{19}F).....	123
Table 1.146: W1346290 - Filter HQ 800 ^1H BP (^{19}F).....	124
Table 1.147: W1346302 - Filter HQ 800 ^{19}F BP (^{-31}P , ^1H).....	124
Table 1.148: W1346294 - Filter HQ 800 ^{19}F BP (^1H).....	125
Table 1.149: W1346306 - Filter HQ 800 ^7Li - ^{31}P BP (^{-87}Rb , ^{19}F -).....	125
Table 1.150: W1346310 - Filter HQ 800 ^{71}Ga - ^{87}Rb BP (^{-13}C , ^7Li -).....	126
Table 1.151: W1346314 - Filter HQ 800 ^{13}C - ^{23}Na BP (^{-29}Si , ^{71}Ga -).....	126
Table 1.152: W1346318 - Filter HQ 800 ^{59}Co - ^{23}Na BP (^{-15}N).....	127
Table 1.153: W1346322 - Filter HQ 800 ^{77}Se - ^{29}Si BP (^{-2}H , ^{13}C -).....	127
Table 1.154: W1346326 - Filter HQ 800 ^{17}O - ^2H BP (^{-15}N , ^{29}Si -).....	128
Table 1.155: W1346330 - Filter HQ 800 ^{15}N BP (^{-14}N , ^2H -).....	128
Table 1.156: W1346677 - Filter HQ 850 0- ^{31}P LP (^{19}F - ^1H).....	129
Table 1.157: W1346656 - Filter HQ 850 0- ^{15}N NR LP (^{13}C).....	130
Table 1.158: W128330 - Filter HQ 850 0- ^{33}S NR LP (^{25}Mg , H3).....	130
Table 1.159: W128329 - Filter HQ 850 0- ^{89}Y NR LP (^{183}W , H3).....	131
Table 1.160: W128331 - Filter HQ 850 0- ^{15}N NR LP (^{91}Zr , H3).....	131
Table 1.161: W131538 - Filter HQ 850 0- ^{37}Cl NR LP (^{33}S , H3).....	132
Table 1.162: W119950 Filter HQ 850 ^{14}N BP.....	133
Table 1.163: W1346655 - Filter HQ 850 ^1H BP (^{-31}P , ^{19}F).....	134
Table 1.164: W1346674- Filter HQ 850 ^1H BP (^{19}F).....	134
Table 1.165: W1346675 - Filter HQ 850 ^{19}F BP (^{-31}P , ^1H).....	135
Table 1.166: W1346676 - Filter HQ 850 ^{19}F BP (^1H).....	135
Table 1.167: W1346715 - Filter HQ 850 ^7Li - ^{31}P BP (^{-87}Rb , ^{19}F -).....	136
Table 1.168: W1346714 - Filter HQ 850 ^{71}Ga - ^{87}Rb BP (^{-13}C , ^7Li -).....	136

Table 1.169:W1346710 - Filter HQ 850 ^{13}C - ^{23}Na BP (^{-29}Si , ^{71}Ga -).....	137
Table 1.170:W1346657 - Filter HQ 850 ^{59}Co - ^{23}Na BP (^{-15}N)	137
Table 1.171:W1346718 - Filter HQ 850 ^{77}Se - ^{29}Si BP (^{-2}H , ^{13}C -).....	138
Table 1.172:W1346697 - Filter HQ 850 ^{17}O - ^2H BP (^{-15}N , ^{29}Si -).....	138
Table 1.173:W1346716- Filter HQ 850 ^{15}N BP (^{-14}N , ^2H).....	139
Table 1.174:W116514 - Diplexer HQ 850 ^1H / ^{19}F	140
Table 1.175:W1346641 - Filter HQ 900 0 - ^{31}P LP (^{19}F - ^1H).....	141
Table 1.176:W1346662 - Filter HQ 900 0 - ^{15}N NR LP (^{13}C).....	142
Table 1.177:W1346287 - Filter HQ 900 ^1H S.....	143
Table 1.178:W1346299 - Filter HQ 900 ^1H BP (^{-31}P , ^{19}F)	144
Table 1.179:W1346291 - Filter HQ 900 ^1H BP (^{19}F)	145
Table 1.180:W1346303 - Filter HQ 900 ^{19}F BP (^{-31}P , ^1H)	145
Table 1.181:W1346295 - Filter HQ 900 ^{19}F BP (^1H)	146
Table 1.182:W1346307 - Filter HQ 900 ^7Li - ^{31}P BP (^{-87}Rb , ^{19}F -).....	146
Table 1.183:W1346311 - Filter HQ 900 ^{71}Ga - ^{87}Rb BP (^{-13}C , ^7Li -)	147
Table 1.184:W1346315 - Filter HQ 900 ^{13}C - ^{23}Na BP (^{-29}Si , ^{71}Ga -).....	147
Table 1.185:W1346319 - Filter HQ 900 ^{59}Co - ^{23}Na BP (^{-15}N)	148
Table 1.186:W1346323 - Filter HQ 900 ^{77}Se - ^{29}Si BP (^{-2}H , ^{13}C -).....	148
Table 1.187:W1346458 - Filter HQ 900 ^{17}O - ^2H BP (^{11}B -).....	149
Table 1.188:W1346327 - Filter HQ 900 ^{17}O - ^2H BP (^{-15}N , ^{29}Si -).....	149
Table 1.189:W1346331 - Filter HQ 900 ^{15}N BP (^{-14}N , ^2H -).....	150
Table 1.190:W1346335 - Filter HQ 900 ^{15}N BP (^{13}C -).....	150
Table 1.191:W1346423 - Diplexer HQ 900 ^1H / ^{19}F	151
Table 1.192:W143786 - Filter HQ 950 0 - ^{31}P LP (^{19}F - ^1H).....	152
Table 1.193:W143787 - Filter HQ 950 0 - ^{15}N NR LP (^{13}C).....	153
Table 1.194:W143773 - Filter HQ 950 ^1H S.....	154
Table 1.195:W1346707 - Filter HQ 950 ^1H BP (^{-31}P , ^{19}F)	155
Table 1.196:W143775 - Filter HQ 950 ^{19}F BP (^{-31}P , ^1H)	156
Table 1.197:W143774 - Filter HQ 950 ^{19}F BP (^1H)	156
Table 1.198:W143776 - Filter HQ 950 ^7Li - ^{31}P BP (^{-87}Rb , ^{19}F -).....	157
Table 1.199:W140205 - Filter HQ 950 ^{71}Ga - ^{87}Rb BP (^{-13}C , ^7Li -)	157
Table 1.200:W140204 - Filter HQ 950 ^{13}C - ^{23}Na BP (^{-29}Si , ^{71}Ga -).....	158
Table 1.201:W143783 - Filter HQ 950 ^{59}Co - ^{23}Na BP (^{-2}H , ^{31}P -).....	158
Table 1.202:W139785- Filter HQ 950 ^{17}O - ^2H BP (^{-15}N , ^{29}Si -).....	159
Table 1.203:W143778 - Filter HQ 950 ^{15}N BP (^{-14}N , ^2H -).....	159
Table 1.204:W143779 - Filter HQ 950 ^{15}N BP (^{13}C -).....	160
Table 1.205:W143780 - Filter HQ 950 ^{14}N BP (^{15}N -).....	160
Table 1.206:W1346869 - Filter HQ 950 ^1H BP (^{19}F)	161
Table 1.207:W140203- Filter HQ 950 ^{77}Se - ^{29}Si BP (^{-2}H , ^{13}C -).....	161
Table 1.208:W134781 - Diplexer HQ 950 ^1H / ^{19}F	162
Table 1.209:W1346678- Filter HQ 1000 0 - ^{31}P LP (^{19}F - ^1H).....	163
Table 1.210:W1346672- Filter HQ 1000 0 - ^{15}N NR LP (^{13}C).....	164
Table 1.211:W1346663 - Filter HQ 1000 ^1H S.....	165

List of Tables

Table 1.212: W1346664 - Filter HQ 1000 ^{19}F S	166
Table 1.213: W1346665 - Filter HQ 1000 ^1H BP (^{-31}P , ^{19}F)	167
Table 1.214: W1346666 - Filter HQ 1000 ^{19}F BP (^{-31}P , ^1H)	168
Table 1.215: W1346667 - Filter HQ 1000 ^7Li - ^{31}P BP (^{-87}Rb , ^{19}F -)	168
Table 1.216: W1346668- Filter HQ 1000 ^{71}Ga - ^{87}Rb BP (^{-13}C , ^7Li -)	169
Table 1.217: W1346669 - Filter HQ 1000 ^{13}C - ^{23}Na BP (^{-29}Si , ^{71}Ga -)	169
Table 1.218: W1346670 - Filter HQ 1000 ^{29}Si BP (^{-2}H , ^{13}C -)	170
Table 1.219: W1346671 - Filter HQ 1000 ^{17}O - ^2H BP (^{-15}N , ^{29}Si -)	170
Table 1.220: W1346717 - Filter HQ 1000 ^{15}N BP (^{-14}N , ^2H -)	171
Table 1.221: W1346856 - Filter HQ 1000 ^{59}Co - ^{23}Na BP (^{-2}H , ^{31}P -)	171
Table 1.222: W1346870 - Filter HQ 1000 ^1H BP (^{19}F)	172
Table 1.223: W1346871 - Filter HQ 1000 ^{19}F BP (^1H)	172
Table 1.224: W1346876 - Diplexer HQ 1000 ^1H / ^{19}F	173
Table 1.225: W1342624 - Filter HQ 1100 0 - ^{15}N NR LP (^{13}C)	174
Table 1.226: W1342373 - Filter HQ 1100 0 - ^{31}P LP (^{19}F - ^1H)	175
Table 1.227: W1342132- Filter HQ 1100 ^1H BP (^{-31}P , ^{19}F)	176
Table 1.228: W1346872- Filter HQ 1100 ^1H BP (^{19}F)	177
Table 1.229: W1346859 - Filter HQ 1100 ^7Li - ^{31}P BP (^{-87}Rb , ^{19}F -)	177
Table 1.230: W1346861- Filter HQ 1100 ^{71}Ga - ^{87}Rb BP (^{-13}C , ^7Li -)	178
Table 1.231: W1346863 - Filter HQ 1100 ^{77}Se - ^{29}Si BP (^{-2}H , ^{13}C -)	178
Table 1.232: W1346865 - Filter HQ 1100 ^{17}O - ^2H BP (^{-15}N , ^{29}Si -)	179
Table 1.233: W1346867- Filter HQ 1100 ^{15}N BP (^{-14}N , ^2H -)	179
Table 1.234: W1346873- Filter HQ 1100 ^{19}F BP (^1H)	180
Table 1.235: W1348165- Filter HQ 1100 ^{13}C - ^{23}Na BP (^{-29}Si , ^{71}Ga -)	180
Table 1.236: W1346473- Filter HQ 1100 ^{19}F BP (^{-31}P , ^1H)	181
Table 1.237: W1346857 - Filter HQ 1100 ^{59}Co - ^{23}Na BP (^{-2}H , ^{31}P -)	181
Table 1.238: W1346877- Diplexer HQ 1100 ^1H / ^{19}F	182
Table 1.239: W134926 - Filter HQ 1200 0 - ^{15}N NR LP (^{13}C)	183
Table 1.240: W134925 - Filter HQ 1200 0 - ^{31}P LP (^{19}F - ^1H)	184
Table 1.241: W1348174 Filter HQ 1200 ^1H BP (^{-31}P , ^{19}F)	185
Table 1.242: W1346874- Filter HQ 1200 ^1H BP (^{19}F)	186
Table 1.243: W1346860 - Filter HQ 1200 ^7Li - ^{31}P BP (^{-87}Rb , ^{19}F -)	186
Table 1.244: W1346862- Filter HQ 1200 ^{71}Ga - ^{87}Rb BP (^{-13}C , ^7Li -)	187
Table 1.245: W1346864 - Filter HQ 1200 ^{77}Se - ^{29}Si BP (^{-2}H , ^{13}C -)	187
Table 1.246: W1346866 - Filter HQ 1200 ^{17}O - ^2H BP (^{-15}N , ^{29}Si -)	188
Table 1.247: W1346868- Filter HQ 1200 ^{15}N BP (^{-14}N , ^2H -)	188
Table 1.248: W1346879- Filter HQ 1200 ^{19}F BP (^1H)	189
Table 1.249: W135028- Filter HQ 1200 ^{13}C - ^{23}Na BP (^{-29}Si , ^{71}Ga -)	189
Table 1.250: W135027 Filter HQ 1200 ^{19}F BP (^{-31}P , ^1H)	190
Table 1.251: W1346879 Filter HQ 1200 ^{19}F BP (^1H)	190
Table 1.252: W1346858 - Filter HQ 1200 ^{59}Co - ^{23}Na BP (^{-2}H , ^{31}P -)	191
Table 1.253: W1346878- Diplexer HQ 1200 ^1H / ^{19}F	192
Table 1.254: Filters HQ 0 - ^{31}P LP (^{19}F - ^1H)	194

Table 1.255: Filters HQ 0- ¹⁵ N NR LP (¹³ C).....	195
Table 1.256: Filters HQ ¹ H S	196
Table 1.257: Filters HQ ¹⁹ F S	197
Table 1.258: Filters HQ ⁷ Li P - ³¹ P S	198
Table 1.259: Filters HQ ³¹ P P - ⁷ Li S	198
Table 1.260: Filters HQ ² H P - ⁷⁷ Se S	198
Table 1.261: Filters HQ ² H S.....	198
Table 1.262: Filters HQ ⁶ Li P - ² H S	198
Table 1.263: Filters HQ ² H P - ⁶ Li S	199
Table 1.264: Filters HQ ¹ H BP (- ³¹ P, ¹⁹ F)	200
Table 1.265: Filters HQ ¹ H BP (¹⁹ F)	201
Table 1.266: Filters HQ ¹⁹ F BP (- ³¹ P, ¹ H)	202
Table 1.267: Filters HQ ¹⁹ F BP (¹ H)	203
Table 1.268: Filters HQ ¹⁹ F - ¹ H BP (- ³¹ P)	204
Table 1.269: Filters HQ ⁷ Li - ³¹ P BP (- ⁸⁷ Rb, ¹⁹ F-).....	204
Table 1.270: Filters HQ ⁷ Li - ³¹ P BP (- ¹³ C, ¹⁹ F-).....	205
Table 1.271: Filters HQ ⁷¹ Ga - ⁸⁷ Rb BP (- ¹³ C, ³¹ P-).....	205
Table 1.272: Filters HQ ⁷¹ Ga - ⁸⁷ Rb BP (- ¹³ C, ⁷ Li-).....	206
Table 1.273: Filters HQ ¹³ C - ²³ Na BP (- ²⁹ Si, ⁷¹ Ga-)	207
Table 1.274: Filters HQ ⁵⁹ Co - ²³ Na BP (- ²⁹ Si, ⁸⁷ Rb-).....	208
Table 1.275: Filters HQ ⁵⁹ Co - ²³ Na BP (- ² H, ¹⁹ F-).....	208
Table 1.276: Filters HQ ⁵⁹ Co - ²³ Na BP (- ¹⁵ N)	209
Table 1.277: Filters HQ ⁷⁷ Se- ²⁹ Si BP (- ² H, ¹³ C-).....	210
Table 1.278: Filters HQ ¹⁷ O - ² H BP (¹¹ B-)	211
Table 1.279: Filters HQ ¹⁷ O - ² H BP (- ¹⁵ N, ²⁹ Si-).....	211
Table 1.280: Filters HQ ¹³³ Cs - ² H BP (- ¹⁵ N, ²⁹ Si-).....	212
Table 1.281: Filters HQ ¹⁵ N BP (- ¹⁴ N, ² H-)	212
Table 1.282: Filters HQ ¹⁵ N BP (¹³ C).....	213
Table 1.283: Filters HQ ¹⁴ N BP (¹⁵ N-)	213
Table 1.284: Filters HQ ⁵⁹ Co - ²³ Na BP (- ² H, ³¹ P).....	213
Table 1.285: Diplexer HQ ¹ H / ¹⁹ F	214

Lastpage

1	NMR Filters	7
1.1	Quick Search - Typical Double Resonance Experiments	8
1.2	Typical Triple Resonance Experiments	12
1.2.1	Quick Search - Triple Resonance Experiments	12
1.2.1.1	300MHz Typical Experiments	14
1.2.1.2	400MHz Typical Experiments	15
1.2.1.3	500MHz Typical Experiments	16
1.2.1.4	600MHz Typical Experiments	17
1.2.1.5	700MHz Typical Experiments	18
1.2.1.6	750MHz Typical Experiments	19
1.2.1.7	800MHz Typical Experiments	20
1.2.1.8	850MHz Typical Experiments	21
1.2.1.9	900MHz Typical Experiments	22
1.2.1.10	950MHz Typical Experiments	23
1.2.1.11	1000MHz Typical Experiments	24
1.2.1.12	1100MHz Typical Experiments	25
1.2.1.13	1200MHz Typical Experiments	26
1.2.2	Triple Resonance 1H/X/Y	27
1.2.3	Triple resonance 1H/19F/X	28
1.3	AVANCE 300MHz	30
1.3.1	Lowpass Filter	30
1.3.2	Bandstop Filter	32
1.3.3	Bandpass Filter	34
1.3.4	Diplexer Filter	42
1.4	AVANCE 400MHz	43
1.4.1	Lowpass Filter	43
1.4.2	Bandstop Filter	45
1.4.3	Bandpass Filter	47
1.4.4	Diplexer Filter	56
1.5	AVANCE 500MHz	57
1.5.1	Lowpass Filter	57
1.5.2	Bandstop Filter	59
1.5.3	Bandpass Filter	64
1.5.4	Diplexer Filter	75
1.6	AVANCE 600MHz	77
1.6.1	Lowpass Filter	77
1.6.2	Bandstop Filter	79
1.6.3	Bandpass Filter	82
1.6.4	Diplexer Filter	92
1.7	AVANCE 700MHz	93
1.7.1	Lowpass Filter	93
1.7.2	Bandstop Filter	95
1.7.3	Bandpass Filter	97

1.7.4	Diplexer Filter	106
1.8	AVANCE 750MHz	107
1.8.1	Lowpass Filter	107
1.8.2	Bandstop Filter	109
1.8.3	Bandpass Filter	111
1.8.4	Diplexer Filter	118
1.9	AVANCE 800MHz	119
1.9.1	Lowpass Filter	119
1.9.2	Bandstop Filter	121
1.9.3	Bandpass Filter	123
1.10	AVANCE 850MHz	129
1.10.1	Lowpass Filter	129
1.10.2	Bandpass Filter	133
1.10.3	Diplexer Filter	140
1.11	AVANCE 900MHz	141
1.11.1	Lowpass Filter	141
1.11.2	Bandstop Filter	143
1.11.3	Bandpass Filter	144
1.11.4	Diplexer Filter	151
1.12	AVANCE 950MHz	152
1.12.1	Lowpass Filter	152
1.12.2	Bandstop Filter	154
1.12.3	Bandpass Filter	155
1.12.4	Diplexer Filter	162
1.13	AVANCE 1000MHz	163
1.13.1	Lowpass Filter	163
1.13.2	Bandstop Filter	165
1.13.3	Bandpass Filter	167
1.13.4	Diplexer Filter	173
1.14	AVANCE 1100MHz	174
1.14.1	Lowpass Filter	174
1.14.2	Bandpass Filter	176
1.14.3	Diplexer Filter	182
1.15	AVANCE 1200MHz	183
1.15.1	Lowpass Filter	183
1.15.2	Bandpass Filter	185
1.15.3	Diplexer Filter	192
1.16	Filters Category	193
1.16.1	Lowpass Filter	194
1.16.2	Bandstop Filter	196
1.16.3	Bandpass Filters	200
1.16.4	Diplexer Filter	214
2	Cases and Accessories	215
2.1	Notch Filters	215
2.1.1	Case 3S063	215
2.1.2	Case 4S063D	215
2.2	Bandpass Filters	216

2.2.1	Case 5R032	216
2.2.2	Case 5R045	216
2.2.3	Case 5R069	217
2.2.4	Case 5R100	217
2.2.5	Case 6R125	218
2.2.6	Case 7R155	218
2.2.7	Case 8R055P	219
2.2.8	Case 8R105P	219
2.2.9	Case 8R041P	220
2.3	Helical Filters	221
2.3.1	Case 4H073	221
2.3.2	Case 5H133	221
2.3.3	Case 5H069	222
2.3.4	Case 6H073	222
2.3.5	Case 8H073	223
2.3.6	Case 8H097	223
2.3.7	Case 8H139	224
2.3.8	Case 8H169	225
2.4	Low Pass Filters	226
2.4.1	Case 2D050	226
2.4.2	Case 2D095	226
2.5	Accessories	227
3	Typical Curves	229
3.1	Bandpass Filters Curves	229
4	Chemicool	237
4.1	Periodic Table	237
4.2	NMR Applications Frequencies	238
4.2.1	NMR Frequency Table	238
4.2.1.1	NMR Frequency Table (1)	239
4.2.1.2	NMR Frequency Table (2)	240
4.2.2	With increasing NMR Frequency	241
4.2.2.1	With increasing NMR Frequency (1)	242
4.2.2.2	With increasing NMR Frequency (2)	243
5	Contact	245
	List of Figures	247
	List of Tables	251

● **Bruker Corporation**

info@bruker.com
www.bruker.com



Order No: Z31481